

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 053 L
 Rev 1

CONTROL L

(Exp.)

2013-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.41	2.38							± 0.15
HB HGB	g/dl	6.0	6.0							± 0.4
	g/l	60	60							± 4
	mmol/l	3.73	3.73							± 0.25
HT HCT	%	16.9	16.9							± 2.0
	l/l	0.169	0.169							± 0.020
VGM MCV	µm ³ ; fl	70	71							± 5
TGMH MCH	pg	24.9	25.2							± 2.0
	fmol	1.55	1.57							± 0.12
CCMH MCHC	g/dl	35.6	35.5							± 3.0
	g/l	356	355							± 30
	mmol/l	22.09	22.05							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	60	65							± 20
VPM MPV	µm ³ ; fl	N/A	8.4							± 3.0
LYMPHO	#	N/A	0.95							± 0.40
	%	N/A	53.0							± 8.0
MONO	#	N/A	0.23							± 0.18
	%	N/A	12.5							± 6.0
GRANULO	#	N/A	0.62							± 0.40
	%	N/A	34.5							± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 053 N
 Rev 1

CONTROL N

(Exp.)

2013-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.0							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.70	4.65							± 0.18
HB HGB	g/dl	13.5	13.4							± 0.5
	g/l	135	134							± 5
	mmol/l	8.38	8.32							± 0.31
HT HCT	%	37.6	37.7							± 2.5
	l/l	0.376	0.377							± 0.025
VGM MCV	µm ³ ; fl	80	81							± 5
TGMH MCH	pg	28.7	28.8							± 2.0
	fmol	1.78	1.79							± 0.12
CCMH MCHC	g/dl	35.9	35.6							± 3.0
	g/l	359	356							± 30
	mmol/l	22.30	22.09							± 1.86
IDR RDW	%	N/A	13.2							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	242	240							± 40
VPM MPV	µm ³ ; fl	N/A	8.1							± 2.0
LYMPHO	#	N/A	1.75							± 0.50
	%	N/A	25.0							± 6.0
MONO	#	N/A	0.70							± 0.40
	%	N/A	10.0							± 5.0
GRANULO	#	N/A	4.55							± 0.60
	%	N/A	65.0							± 7.0

9930072-A
 Ref: TEMP-0630 Rev.31

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 053 H
 Rev 1

CONTROL H

(Exp.)

2013-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	18.9	19.0							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.70	5.64							± 0.20
HB HGB	g/dl	17.9	17.8							± 0.6
	g/l	179	178							± 6
	mmol/l	11.12	11.05							± 0.37
HT HCT	%	50.2	50.2							± 3.0
	l/l	0.502	0.502							± 0.030
VGM MCV	µm ³ ; fl	88	89							± 5
TGMH MCH	pg	31.4	31.6							± 2.0
	fmol	1.95	1.96							± 0.12
CCMH MCHC	g/dl	35.7	35.5							± 3.0
	g/l	357	355							± 30
	mmol/l	22.16	22.02							± 1.86
IDR RDW	%	N/A	12.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	525	495							± 45
VPM MPV	µm ³ ; fl	N/A	8.0							± 2.0
LYMPHO	#	N/A	2.09							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	0.95							± 0.76
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	15.96							± 1.40
	%	N/A	84.0							± 7.0