

ABX Minotrol 16



LOT MX 091 L

CONTROL L

(Exp.) 2011-12-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	2.0							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.40	2.47							± 0.15
HB HGB	g/dl	5.8	6.1							± 0.4
	g/l	58	61							± 4
	mmol/l	3.60	3.76							± 0.25
HT HCT	%	16.8	17.0							± 2.0
	l/l	0.168	0.170							± 0.020
VGM MCV	µm ³ ; fl	70	69							± 4
TGMH MCH	pg	24.2	24.5							± 2.0
	fmol	1.50	1.52							± 0.12
CCMH MCHC	g/dl	34.5	35.7							± 3.0
	g/l	345	357							± 30
	mmol/l	21.44	22.19							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	70	83							± 20
VPM MPV	µm ³ ; fl	N/A	7.0							± 2.0
LYMPHO	#	N/A	1.03							± 0.40
	%	N/A	51.8							± 7.0
MONO	#	N/A	0.27							± 0.20
	%	N/A	13.7							± 6.0
GRANULO	#	N/A	0.67							± 0.40
	%	N/A	33.7							± 7.0

Ref: FORM-0830 Rev.15 9930041-E



ABX Minotrol 16



LOT MX 091 N

CONTROL N

(Exp.) 2011-12-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.3	7.5							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.64	4.72							± 0.18
HB HGB	g/dl	13.5	13.3							± 0.5
	g/l	135	133							± 5
	mmol/l	8.38	8.28							± 0.31
HT HCT	%	38.0	37.9							± 2.5
	l/l	0.380	0.379							± 0.025
VGM MCV	µm ³ ; fl	82	80							± 4
TGMH MCH	pg	29.1	28.3							± 2.0
	f mol	1.81	1.75							± 0.12
CCMH MCHC	g/dl	35.5	35.2							± 3.0
	g/l	355	352							± 30
	mmol/l	22.03	21.84							± 1.86
IDR RDW	%	N/A	13.8							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	250	252							± 40
VPM MPV	µm ³ ; fl	N/A	7.2							± 2.0
LYMPHO	#	N/A	1.96							± 0.50
	%	N/A	26.0							± 6.0
MONO	#	N/A	0.71							± 0.40
	%	N/A	9.4							± 5.0
GRANULO	#	N/A	4.88							± 0.60
	%	N/A	64.8							± 7.0

Ref: FORM-0830 Rev.15 9930042-E



ABX Minotrol 16



LOT **MX 091 H**

CONTROL **H**

(Exp.) **2011-12-05**
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.4	20.5							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.69	5.76							± 0.20
HB HGB	g/dl	17.7	17.9							± 0.6
	g/l	177	179							± 6
	mmol/l	10.99	11.12							± 0.37
HT HCT	%	50.1	50.2							± 3.0
	l/l	0.501	0.502							± 0.030
VGM MCV	µm ³ ; fl	88	87							± 4
TGMH MCH	pg	31.1	31.1							± 2.0
	f mol	1.93	1.93							± 0.12
CCMH MCHC	g/dl	35.3	35.6							± 3.0
	g/l	353	356							± 30
	mmol/l	21.95	22.14							± 1.86
IDR RDW	%	N/A	12.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	520	510							± 55
VPM MPV	µm ³ ; fl	N/A	7.4							± 2.0
LYMPHO	#	N/A	2.21							± 1.00
	%	N/A	10.8							± 5.0
MONO	#	N/A	1.11							± 0.80
	%	N/A	5.4							± 4.0
GRANULO	#	N/A	17.22							± 1.40
	%	N/A	84.0							± 7.0

Ref: FORM-0830 Rev.15 9930043-E