



ABX Minotrol 16



LOT MX 071 L

CONTROL L



(Exp.) 2011-10-05

(YYYY - MM - DD)

Rev 1

		ABX VETPACK								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS								TOLERANCES TOLERANCE
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	2.0							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.36	2.45							± 0.15
HB HGB	g/dl	5.9	6.1							± 0.4
	g/l	59	61							± 4
	mmol/l	3.66	3.76							± 0.25
HT HCT	%	16.5	16.8							± 2.0
	l/l	0.165	0.168							± 0.020
VGM MCV	µm ³ ; fl	70	69							± 4
TGMH MCH	pg	25.0	24.7							± 2.0
	fmol	1.55	1.54							± 0.12
CCMH MCHC	g/dl	35.7	36.0							± 3.0
	g/l	357	360							± 30
	mmol/l	22.18	22.38							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	63	83							± 20
VPM MPV	µm ³ ; fl	N/A	7.0							± 2.0
LYMPHO	#	N/A	1.00							± 0.40
	%	N/A	50.6							± 7.0
MONO	#	N/A	0.26							± 0.20
	%	N/A	13.0							± 6.0
GRANULO	#	N/A	0.72							± 0.40
	%	N/A	36.3							± 7.0

Ref: FORM-0830 Rev.15 9930041-E



ABX Minotrol 16



LOT MX 071 N
 Rev 1

CONTROL N

(Exp.) 2011-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.3							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.65	4.71							± 0.18
HB HGB	g/dl	13.6	13.2							± 0.5
	g/l	136	132							± 5
	mmol/l	8.45	8.22							± 0.31
HT HCT	%	38.6	38.3							± 2.5
	l/l	0.386	0.383							± 0.025
VGM MCV	µm ³ ; fl	83	81							± 4
TGMH MCH	pg	29.2	28.1							± 2.0
	fmol	1.82	1.75							± 0.12
CCMH MCHC	g/dl	35.2	34.6							± 3.0
	g/l	352	346							± 30
	mmol/l	21.88	21.46							± 1.86
IDR RDW	%	N/A	13.8							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260	257							± 40
VPM MPV	µm ³ ; fl	N/A	7.2							± 2.0
LYMPHO	#	N/A	1.93							± 0.50
	%	N/A	26.4							± 6.0
MONO	#	N/A	0.64							± 0.40
	%	N/A	8.8							± 5.0
GRANULO	#	N/A	4.75							± 0.60
	%	N/A	64.8							± 7.0

Ref: FORM-0830 Rev.15 9930042-E



ABX Minotrol 16



LOT MX 071 H
 Rev 1

CONTROL H

(Exp.) 2011-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.0	20.6							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.70	5.80							± 0.20
HB HGB	g/dl	17.9	18.0							± 0.6
	g/l	179	180							± 6
	mmol/l	11.12	11.18							± 0.37
HT HCT	%	50.7	51.1							± 3.0
	l/l	0.507	0.511							± 0.030
VGM MCV	µm ³ ; fl	89	88							± 4
TGMH MCH	pg	31.4	31.0							± 2.0
	fmol	1.95	1.93							± 0.12
CCMH MCHC	g/dl	35.3	35.2							± 3.0
	g/l	353	352							± 30
	mmol/l	21.91	21.86							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	522	510							± 45
VPM MPV	µm ³ ; fl	N/A	7.6							± 2.0
LYMPHO	#	N/A	2.33							± 1.00
	%	N/A	11.3							± 5.0
MONO	#	N/A	1.07							± 0.80
	%	N/A	5.2							± 4.0
GRANULO	#	N/A	17.20							± 1.40
	%	N/A	83.5							± 7.0

Ref: FORM-0830 Rev.15 9930043-E