



ABX Minotrol 16



LOT MX 051 L

CONTROL L



(Exp.) 2011-08-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.9	2.1							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.41	2.48							± 0.15
HB HGB	g/dl	6.0	6.2							± 0.4
	g/l	60	62							± 4
	mmol/l	3.73	3.83							± 0.25
HT HCT	%	16.9	17.0							± 2.0
	l/l	0.169	0.170							± 0.020
VGM MCV	µm ³ ; fl	70	69							± 4
TGMH MCH	pg	24.9	24.8							± 2.0
	fmol	1.55	1.54							± 0.12
CCMH MCHC	g/dl	35.6	36.2							± 3.0
	g/l	356	362							± 30
	mmol/l	22.09	22.47							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	94							± 30
VPM MPV	µm ³ ; fl	N/A	7.4							± 3.0
LYMPHO	#	N/A	1.09							± 0.40
	%	N/A	52.2							± 7.0
MONO	#	N/A	0.30							± 0.20
	%	N/A	14.3							± 6.0
GRANULO	#	N/A	0.67							± 0.40
	%	N/A	32.3							± 7.0

Ref: FORM-0830 Rev.15 9930041-E



ABX Minotrol 16



LOT MX 051 N

CONTROL N

(Exp.) 2011-08-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.1	7.6							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.62	4.65							± 0.18
HB HGB	g/dl	13.6	13.2							± 0.5
	g/l	136	132							± 5
	mmol/l	8.45	8.22							± 0.31
HT HCT	%	38.3	37.8							± 2.5
	l/l	0.383	0.378							± 0.025
VGM MCV	µm ³ ; fl	83	81							± 4
TGMH MCH	pg	29.4	28.5							± 2.0
	fmol	1.83	1.77							± 0.12
CCMH MCHC	g/dl	35.5	35.0							± 3.0
	g/l	355	350							± 30
	mmol/l	22.02	21.74							± 1.86
IDR RDW	%	N/A	13.8							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	275	267							± 45
VPM MPV	µm ³ ; fl	N/A	7.4							± 3.0
LYMPHO	#	N/A	1.98							± 0.50
	%	N/A	26.0							± 6.0
MONO	#	N/A	0.71							± 0.40
	%	N/A	9.4							± 5.0
GRANULO	#	N/A	4.94							± 0.60
	%	N/A	64.8							± 7.0

Ref: FORM-0830 Rev.15 9930042-E



ABX Minotrol 16



LOT MX 051 H
 Rev 1

CONTROL H

(Exp.) 2011-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.3	20.8							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.63	5.69							± 0.20
HB HGB	g/dl	17.8	17.9							± 0.6
	g/l	178	179							± 6
	mmol/l	11.05	11.12							± 0.37
HT HCT	%	50.1	50.2							± 3.0
	l/l	0.501	0.502							± 0.030
VGM MCV	µm ³ ; fl	89	88							± 4
TGMH MCH	pg	31.6	31.5							± 2.0
	fmol	1.96	1.95							± 0.12
CCMH MCHC	g/dl	35.5	35.7							± 3.0
	g/l	355	357							± 30
	mmol/l	22.06	22.16							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	540	510							± 55
VPM MPV	µm ³ ; fl	N/A	7.5							± 3.0
LYMPHO	#	N/A	2.28							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	0.97							± 0.80
	%	N/A	4.7							± 4.0
GRANULO	#	N/A	17.59							± 1.40
	%	N/A	84.5							± 7.0

Ref: FORM-0830 Rev.15 9930043-E