



ABX Minotrol 16



LOT MX 031 L

CONTROL L

(Exp.) 2011-06-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8	2.0							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.46	2.50							± 0.15
HB HGB	g/dl	6.1	6.1							± 0.4
	g/l	61	61							± 4
	mmol/l	3.79	3.76							± 0.25
HT HCT	%	17.0	16.9							± 2.0
	l/l	0.170	0.169							± 0.020
VGM MCV	µm ³ ; fl	69	68							± 4
TGMH MCH	pg	24.8	24.2							± 2.0
	fmol	1.54	1.50							± 0.12
CCMH MCHC	g/dl	35.9	35.8							± 3.0
	g/l	359	358							± 30
	mmol/l	22.32	22.24							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	70	88							± 20
VPM MPV	µm ³ ; fl	N/A	7.4							± 2.0
LYMPHO	#	N/A	1.01							± 0.40
	%	N/A	51.0							± 6.0
MONO	#	N/A	0.28							± 0.20
	%	N/A	14.3							± 5.0
GRANULO	#	N/A	0.68							± 0.40
	%	N/A	34.3							± 6.0

Ref: FORM-0830 Rev.14 9930041-E



ABX Minotrol 16



LOT

MX 031 N

CONTROL N



(Exp.) 2011-06-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.5							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.80	4.85							± 0.18
HB HGB	g/dl	13.8	13.5							± 0.5
	g/l	138	135							± 5
	mmol/l	8.57	8.40							± 0.31
HT HCT	%	38.4	38.5							± 2.5
	l/l	0.384	0.385							± 0.025
VGM MCV	µm ³ ; fl	80	79							± 4
TGMH MCH	pg	28.8	27.9							± 2.0
	fmol	1.79	1.73							± 0.12
CCMH MCHC	g/dl	35.9	35.1							± 3.0
	g/l	359	351							± 30
	mmol/l	22.32	21.82							± 1.86
IDR RDW	%	N/A	13.9							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	252	257							± 40
VPM MPV	µm ³ ; fl	N/A	7.5							± 2.0
LYMPHO	#	N/A	1.99							± 0.40
	%	N/A	26.4							± 5.0
MONO	#	N/A	0.66							± 0.30
	%	N/A	8.8							± 4.0
GRANULO	#	N/A	4.88							± 0.60
	%	N/A	64.8							± 6.0

Ref: FORM-0830 Rev.14 9930042-E



ABX Minotrol 16



LOT MX 031 H
 Rev 1

CONTROL H

(Exp.) 2011-06-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.4	20.8							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.70	5.79							± 0.20
HB HGB	g/dl	18.0	18.1							± 0.6
	g/l	180	181							± 6
	mmol/l	11.18	11.24							± 0.37
HT HCT	%	50.7	50.5							± 3.0
	l/l	0.507	0.505							± 0.030
VGM MCV	µm ³ ; fl	89	87							± 4
TGMH MCH	pg	31.6	31.3							± 2.0
	fmol	1.96	1.94							± 0.12
CCMH MCHC	g/dl	35.5	35.9							± 3.0
	g/l	355	359							± 30
	mmol/l	22.03	22.27							± 1.86
IDR RDW	%	N/A	12.7							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	510	500							± 45
VPM MPV	µm ³ ; fl	N/A	7.4							± 2.0
LYMPHO	#	N/A	2.28							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	0.97							± 0.80
	%	N/A	4.7							± 4.0
GRANULO	#	N/A	17.59							± 1.40
	%	N/A	84.5							± 7.0

Ref: FORM-0630 Rev.14 9930043-E