



ABX Minotrol 16



LOT MX 011 L

CONTROL L

(Exp.) 2011-04-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8	2.0							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.46	2.51							± 0.15
HB HGB	g/dl	6.1	6.2							± 0.4
	g/l	61	62							± 4
	mmol/l	3.79	3.83							± 0.25
HT HCT	%	17.0	17.0							± 2.0
	l/l	0.170	0.170							± 0.020
VGM MCV	µm ³ ; fl	69	68							± 4
TGMH MCH	pg	24.8	24.5							± 2.0
	fmol	1.54	1.52							± 0.12
CCMH MCHC	g/dl	35.9	36.3							± 3.0
	g/l	359	363							± 30
	mmol/l	22.32	22.51							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	70	88							± 20
VPM MPV	µm ³ ; fl	N/A	7.3							± 3.0
LYMPHO	#	N/A	0.99							± 0.40
	%	N/A	50.2							± 7.0
MONO	#	N/A	0.27							± 0.20
	%	N/A	13.7							± 6.0
GRANULO	#	N/A	0.72							± 0.40
	%	N/A	36.3							± 7.0

Ref: FORM-0830 Rev.14 9930041-E

ABX Minotrol 16



LOT MX 011 N
 Rev 1

CONTROL N

(Exp.) 2011-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.2	7.4							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.65							± 0.18
HB HGB	g/dl	13.6	13.3							± 0.5
	g/l	136	133							± 5
	mmol/l	8.45	8.28							± 0.31
HT HCT	%	38.2	37.8							± 2.5
	l/l	0.382	0.378							± 0.025
VGM MCV	µm ³ ; fl	83	81							± 4
TGMH MCH	pg	29.6	28.7							± 2.0
	fmol	1.84	1.78							± 0.12
CCMH MCHC	g/dl	35.6	35.3							± 3.0
	g/l	356	353							± 30
	mmol/l	22.12	21.90							± 1.86
IDR RDW	%	N/A	13.8							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	266	262							± 40
VPM MPV	µm ³ ; fl	N/A	7.3							± 2.0
LYMPHO	#	N/A	2.00							± 0.50
	%	N/A	27.0							± 6.0
MONO	#	N/A	0.68							± 0.40
	%	N/A	9.1							± 5.0
GRANULO	#	N/A	4.73							± 0.60
	%	N/A	63.7							± 7.0

Ref: FORM-0630 Rev.14 9930042-E

ABX Minotrol 16



LOT MX 011 H
 Rev 1

CONTROL H

(Exp.) 2011-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.8	21.1							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.68	5.77							± 0.20
HB HGB	g/dl	18.0	18.2							± 0.6
	g/l	180	182							± 6
	mmol/l	11.18	11.30							± 0.37
HT HCT	%	50.6	50.9							± 3.0
	l/l	0.506	0.509							± 0.030
VGM MCV	µm ³ ; fl	89	88							± 4
TGMH MCH	pg	31.7	31.6							± 2.0
	fmol	1.97	1.96							± 0.12
CCMH MCHC	g/dl	35.6	35.8							± 3.0
	g/l	356	358							± 30
	mmol/l	22.11	22.22							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	540	519							± 45
VPM MPV	µm ³ ; fl	N/A	7.5							± 2.0
LYMPHO	#	N/A	2.31							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	0.99							± 0.80
	%	N/A	4.7							± 4.0
GRANULO	#	N/A	17.85							± 1.40
	%	N/A	84.5							± 7.0

Ref: FORM-0830 Rev.14 9930043-E