



ABX Minotrol 16



LOT MX 110 L

CONTROL L

(Exp.) 2011-02-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.7	1.9							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.34	2.36							± 0.15
HB HGB	g/dl	5.8	6.0							± 0.4
	g/l	58	60							± 4
	mmol/l	3.60	3.70							± 0.25
HT HCT	%	16.4	16.4							± 2.0
	l/l	0.164	0.164							± 0.020
VGM MCV	µm ³ ; fl	70	70							± 4
TGMH MCH	pg	24.8	25.3							± 2.0
	fmol	1.54	1.57							± 0.12
CCMH MCHC	g/dl	35.4	36.3							± 3.0
	g/l	354	363							± 30
	mmol/l	21.99	22.55							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	70	86							± 20
VPM MPV	µm ³ ; fl	N/A	7.4							± 3.0
LYMPHO	#	N/A	0.98							± 0.40
	%	N/A	51.8							± 7.0
MONO	#	N/A	0.24							± 0.20
	%	N/A	13.0							± 6.0
GRANULO	#	N/A	0.65							± 0.40
	%	N/A	34.3							± 7.0

Ref: FORM-0830 Rev.14 9930041-E



ABX Minotrol 16



LOT MX 110 N

CONTROL N

(Exp.) 2011-02-05
 (YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.3							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.65	4.69							± 0.18
HB HGB	g/dl	13.4	13.2							± 0.5
	g/l	134	132							± 5
	mmol/l	8.32	8.22							± 0.31
HT HCT	%	38.1	37.7							± 2.5
	l/l	0.381	0.377							± 0.025
VGM MCV	µm ³ ; fl	82	80							± 4
TGMH MCH	pg	28.8	28.2							± 2.0
	fmol	1.79	1.75							± 0.12
CCMH MCHC	g/dl	35.1	35.1							± 3.0
	g/l	351	351							± 30
	mmol/l	21.82	21.82							± 1.86
IDR RDW	%	N/A	13.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	257	259							± 40
VPM MPV	µm ³ ; fl	N/A	7.5							± 2.0
LYMPHO	#	N/A	1.93							± 0.50
	%	N/A	26.4							± 6.0
MONO	#	N/A	0.64							± 0.40
	%	N/A	8.8							± 5.0
GRANULO	#	N/A	4.75							± 0.60
	%	N/A	64.8							± 7.0

Ref: FORM-0630 Rev.14 9930042-E



ABX Minotrol 16



LOT MX 110 H
 Rev 1

CONTROL H

(Exp.) 2011-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	19.8	21.1							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.55	5.66							± 0.20
HB HGB	g/dl	17.6	17.8							± 0.6
	g/l	176	178							± 6
	mmol/l	10.93	11.05							± 0.37
HT HCT	%	49.4	49.9							± 3.0
	l/l	0.494	0.499							± 0.030
VGM MCV	µm ³ ; fl	89	88							± 4
TGMH MCH	pg	31.7	31.5							± 2.0
	fmol	1.97	1.95							± 0.12
CCMH MCHC	g/dl	35.6	35.7							± 3.0
	g/l	356	357							± 30
	mmol/l	22.13	22.16							± 1.86
IDR RDW	%	N/A	12.9							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	521	500							± 45
VPM MPV	µm ³ ; fl	N/A	7.6							± 2.0
LYMPHO	#	N/A	2.39							± 1.00
	%	N/A	11.3							± 5.0
MONO	#	N/A	0.99							± 0.80
	%	N/A	4.7							± 4.0
GRANULO	#	N/A	17.74							± 1.40
	%	N/A	84.0							± 7.0

Ref: FORM-0630 Rev.14 9930043-E