

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédécine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT MX 033 L
Rev 1

CONTROL L

(Exp.)

2013-06-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	MICROS		ABX VETPACK						TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.8								± 0.4
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.41	2.38								± 0.15
HB HGB	g/dl	6.0	6.0								± 0.4
	g/l	60	60								± 4
	mmol/l	3.73	3.73								± 0.25
HT HCT	%	16.6	16.7								± 2.0
	l/l	0.166	0.167								± 0.020
VGM MCV	µm ³ ; fl	69	70								± 5
TGMH MCH	pg	24.9	25.2								± 2.0
	fmol	1.55	1.57								± 0.12
CCMH MCHC	g/dl	36.1	36.0								± 3.0
	g/l	361	360								± 30
	mmol/l	22.41	22.36								± 1.86
IDR RDW	%	N/A	13.5								± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	70	74								± 20
VPM MPV	µm ³ ; fl	N/A	8.0								± 2.0
LYMPHO	#	N/A	0.99								± 0.40
	%	N/A	55.0								± 8.0
MONO	#	N/A	0.23								± 0.20
	%	N/A	13.0								± 6.0
GRANULO	#	N/A	0.58								± 0.40
	%	N/A	32.0								± 7.0

LOT MX 033 N
Rev 1

CONTROL N

(Exp.)

2013-06-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	MICROS		ABX VETPACK						TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.4	7.6								± 0.8
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.71	4.68								± 0.18
HB HGB	g/dl	13.6	13.4								± 0.5
	g/l	136	134								± 5
	mmol/l	8.45	8.32								± 0.31
HT HCT	%	37.7	37.4								± 2.5
	l/l	0.377	0.374								± 0.025
VGM MCV	µm ³ ; fl	80	80								± 5
TGMH MCH	pg	28.9	28.6								± 2.0
	fmol	1.79	1.78								± 0.12
CCMH MCHC	g/dl	36.1	35.8								± 3.0
	g/l	361	358								± 30
	mmol/l	22.41	22.23								± 1.86
IDR RDW	%	N/A	13.2								± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	277	264								± 40
VPM MPV	µm ³ ; fl	N/A	8.0								± 2.0
LYMPHO	#	N/A	1.98								± 0.40
	%	N/A	26.0								± 6.0
MONO	#	N/A	0.76								± 0.30
	%	N/A	10.0								± 5.0
GRANULO	#	N/A	4.86								± 0.60
	%	N/A	64.0								± 7.0

960072-A
Ref: TEMP-0330 Rev.30

LOT MX 033 H
Rev 1

CONTROL H

(Exp.)

2013-06-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	MICROS		ABX VETPACK						TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.3	20.5								± 1.6
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.52	5.52								± 0.20
HB HGB	g/dl	17.7	17.7								± 0.6
	g/l	177	177								± 6
	mmol/l	10.99	10.99								± 0.37
HT HCT	%	49.1	49.7								± 3.0
	l/l	0.491	0.497								± 0.030
VGM MCV	µm ³ ; fl	89	90								± 5
TGMH MCH	pg	32.1	32.1								± 2.0
	fmol	1.99	1.99								± 0.12
CCMH MCHC	g/dl	36.0	35.6								± 3.0
	g/l	360	356								± 30
	mmol/l	22.37	22.13								± 1.86
IDR RDW	%	N/A	12.8								± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	565	526								± 45
VPM MPV	µm ³ ; fl	N/A	7.8								± 2.0
LYMPHO	#	N/A	2.26								± 1.00
	%	N/A	11.0								± 5.0
MONO	#	N/A	1.03								± 0.82
	%	N/A	5.0								± 4.0
GRANULO	#	N/A	17.22								± 1.40
	%	N/A	84.0								± 7.0