




LOT		MX 055 L	CONTROL L							 (Exp.)	2015-08-05 (YYYY-MM-DD)
		Rev 1	ABX VETPACK								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS								TOLERANCES TOLERANCE	
		Sci1 Vet ABC	Sci1 Vet ABC Plus ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.45	2.38							± 0.15	
HB HGB	g/dl	6.0	6.1							± 0.4	
	g/l	60	61							± 4	
	mmol/l	3.73	3.79							± 0.25	
HT HCT	%	16.9	16.4							± 2.0	
	l/l	0.169	0.164							± 0.020	
VGM MCV	µm ³ ; fl	69	69							± 4	
TGMH MCH	pg	24.5	25.6							± 2.0	
	fmol	1.52	1.59							± 0.12	
CCMH MCHC	g/dl	35.5	37.1							± 3.0	
	g/l	355	371							± 30	
	mmol/l	22.04	23.07							± 1.86	
IDR RDW	%	N/A	13.5							± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	68	70							± 20	
VPM MPV	µm ³ ; fl	N/A	9.0							± 2.0	
LYMPHO	#	N/A	1.03							± 0.40	
	%	N/A	54.0							± 8.0	
MONO	#	N/A	0.26							± 0.20	
	%	N/A	13.5							± 6.0	
GRANULO	#	N/A	0.62							± 0.40	
	%	N/A	32.5							± 7.0	

LOT		MX 055 N	CONTROL N							 (Exp.)	2015-08-05 (YYYY-MM-DD)	
Rev 1		ABX VETPACK										
PARAMETRES PARAMETERS	UNITES UNITS	MICROS										TOLERANCES TOLERANCE
		Sci Vet ABC	Sci Vet ABC Plus ESV 60									
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.3									± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.63	4.59									± 0.18
HB HGB	g/dl	13.5	13.5									± 0.5
	g/l	135	135									± 5
	mmol/l	8.38	8.38									± 0.31
HT HCT	%	37.5	37.2									± 2.5
	l/l	0.375	0.372									± 0.025
VGM MCV	µm ³ ; fl	81	81									± 4
TGMH MCH	pg	29.2	29.4									± 2.0
	fmol	1.81	1.83									± 0.12
CCMH MCHC	g/dl	36.0	36.3									± 3.0
	g/l	360	363									± 30
	mmol/l	22.35	22.55									± 1.86
IDR RDW	%	N/A	13.0									± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	256	257									± 40
VPM MPV	µm ³ ; fl	N/A	8.2									± 2.0
LYMPHO	#	N/A	1.97									± 0.40
	%	N/A	27.0									± 5.0
MONO	#	N/A	0.66									± 0.30
	%	N/A	9.0									± 4.0
GRANULO	#	N/A	4.67									± 0.60
	%	N/A	64.0									± 6.0

9930072-A
Ref: TEMP-0830 Rev.36

LOT		MX 055 H	CONTROL H							 (Exp.)	2015-08-05 (YYYY-MM-DD)
		Rev 1	ABX VETPACK								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS									TOLERANCES TOLERANCE
		Scil Vet ABC	Scil Vet ABC Plus ESV 60								
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.4	19.8								± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.76	5.73								± 0.20
HB HGB	g/dl	18.4	18.4								± 0.6
	g/l	184	184								± 6
	mmol/l	11.43	11.43								± 0.37
HT HCT	%	50.7	51.0								± 3.0
	l/l	0.507	0.510								± 0.030
VGM MCV	µm ³ ; fl	88	89								± 4
TGMH MCH	pg	31.9	32.1								± 2.0
	fmol	1.98	1.99								± 0.12
CCMH MCHC	g/dl	36.3	36.1								± 3.0
	g/l	363	361								± 30
	mmol/l	22.54	22.41								± 1.86
IDR RDW	%	N/A	12.5								± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	515	503								± 45
VPM MPV	µm ³ ; fl	N/A	8.5								± 2.0
LYMPHO	#	N/A	2.48								± 1.00
	%	N/A	12.5								± 5.0
MONO	#	N/A	1.09								± 0.80
	%	N/A	5.5								± 4.0
GRANULO	#	N/A	16.24								± 1.40
	%	N/A	82.0								± 7.0