



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 015 L
 Rev 1

CONTROL L

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		SciI Vet ABC	SciI Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.7	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.36	2.34							± 0.15
HB HGB	g/dl	6.1	6.1							± 0.4
	g/l	61	61							± 4
	mmol/l	3.79	3.79							± 0.25
HT HCT	%	16.5	16.4							± 2.0
	l/l	0.165	0.164							± 0.020
VGM MCV	µm ³ ; fl	70	70							± 4
TGMH MCH	pg	25.8	26.1							± 2.0
	fmol	1.61	1.62							± 0.12
CCMH MCHC	g/dl	36.9	37.2							± 3.0
	g/l	369	372							± 30
	mmol/l	22.93	23.13							± 1.86
IDR RDW	%	N/A	13.5							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	69	73							± 20
VPM MPV	µm ³ ; fl	N/A	8.1							± 2.0
LYMPHO	#	N/A	0.97							± 0.40
	%	N/A	54.0							± 8.0
MONO	#	N/A	0.20							± 0.20
	%	N/A	11.0							± 6.0
GRANULO	#	N/A	0.63							± 0.40
	%	N/A	35.0							± 7.0



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 015 N
 Rev 1

CONTROL N

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	6.9	7.1							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.57							± 0.18
HB HGB	g/dl	13.6	13.5							± 0.5
	g/l	136	135							± 5
	mmol/l	8.45	8.38							± 0.31
HT HCT	%	37.3	37.0							± 2.5
	l/l	0.373	0.370							± 0.025
VGM MCV	µm ³ ; fl	81	81							± 4
TGMH MCH	pg	29.6	29.5							± 2.0
	fmol	1.84	1.83							± 0.12
CCMH MCHC	g/dl	36.5	36.5							± 3.0
	g/l	365	365							± 30
	mmol/l	22.67	22.65							± 1.86
IDR RDW	%	N/A	13.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	245	250							± 40
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	1.92							± 0.40
	%	N/A	27.0							± 5.0
MONO	#	N/A	0.57							± 0.30
	%	N/A	8.0							± 4.0
GRANULO	#	N/A	4.62							± 0.60
	%	N/A	65.0							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.36



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 015 H
 Rev 1

CONTROL H

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.5	19.6							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.67	5.65							± 0.20
HB HGB	g/dl	18.1	18.0							± 0.6
	g/l	181	180							± 6
	mmol/l	11.24	11.18							± 0.37
HT HCT	%	49.9	49.7							± 3.0
	l/l	0.499	0.497							± 0.030
VGM MCV	µm ³ ; fl	88	88							± 4
TGMH MCH	pg	31.9	31.9							± 2.0
	fmol	1.98	1.98							± 0.12
CCMH MCHC	g/dl	36.3	36.2							± 3.0
	g/l	363	362							± 30
	mmol/l	22.53	22.48							± 1.86
IDR RDW	%	N/A	12.5							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	505	505							± 45
VPM MPV	µm ³ ; fl	N/A	7.8							± 2.0
LYMPHO	#	N/A	2.16							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	0.98							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.46							± 1.40
	%	N/A	84.0							± 7.0