



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 093 L
 Rev 1

CONTROL L

(Exp.)

2013-12-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.44	2.42							± 0.15
HB HGB	g/dl	6.0	6.0							± 0.4
	g/l	60	60							± 4
	mmol/l	3.73	3.73							± 0.25
HT HCT	%	16.8	16.5							± 2.0
	l/l	0.168	0.165							± 0.020
VGM MCV	µm ³ ; fl	69	68							± 4
TGMH MCH	pg	24.6	24.8							± 2.0
	fmol	1.53	1.54							± 0.12
CCMH MCHC	g/dl	35.6	36.5							± 3.0
	g/l	356	365							± 30
	mmol/l	22.13	22.64							± 1.86
IDR RDW	%	N/A	14.2							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	65	75							± 20
VPM MPV	µm ³ ; fl	N/A	8.7							± 3.0
LYMPHO	#	N/A	0.99							± 0.40
	%	N/A	55.0							± 6.0
MONO	#	N/A	0.21							± 0.20
	%	N/A	11.5							± 5.0
GRANULO	#	N/A	0.60							± 0.40
	%	N/A	33.5							± 6.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 093 N
 Rev 1

CONTROL N

(Exp.)

2013-12-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.1							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.56	4.53							± 0.18
HB HGB	g/dl	13.5	13.4							± 0.5
	g/l	135	134							± 5
	mmol/l	8.38	8.32							± 0.31
HT HCT	%	37.4	37.1							± 2.5
	l/l	0.374	0.371							± 0.025
VGM MCV	µm ³ ; fl	82	82							± 4
TGMH MCH	pg	29.6	29.6							± 2.0
	fmol	1.84	1.84							± 0.12
CCMH MCHC	g/dl	36.1	36.1							± 3.0
	g/l	361	361							± 30
	mmol/l	22.42	22.40							± 1.86
IDR RDW	%	N/A	13.7							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260	260							± 40
VPM MPV	µm ³ ; fl	N/A	8.4							± 2.0
LYMPHO	#	N/A	1.88							± 0.40
	%	N/A	26.5							± 5.0
MONO	#	N/A	0.64							± 0.30
	%	N/A	9.0							± 4.0
GRANULO	#	N/A	4.58							± 0.60
	%	N/A	64.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.32



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 093 H
 Rev 1

CONTROL H

(Exp.)

2013-12-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.0	19.3							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.68	5.64							± 0.20
HB HGB	g/dl	17.8	17.7							± 0.6
	g/l	178	177							± 6
	mmol/l	11.05	10.99							± 0.37
HT HCT	%	50.0	50.2							± 3.0
	l/l	0.500	0.502							± 0.030
VGM MCV	µm ³ ; fl	88	89							± 4
TGMH MCH	pg	31.3	31.4							± 2.0
	fmol	1.95	1.95							± 0.12
CCMH MCHC	g/dl	35.6	35.3							± 3.0
	g/l	356	353							± 30
	mmol/l	22.11	21.90							± 1.86
IDR RDW	%	N/A	12.8							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	510	485							± 45
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	2.32							± 1.00
	%	N/A	12.0							± 5.0
MONO	#	N/A	0.97							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.02							± 1.40
	%	N/A	83.0							± 7.0