



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 114 L
 Rev 1

CONTROL L

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.30	2.29							± 0.15
HB HGB	g/dl	5.9	5.9							± 0.4
	g/l	59	59							± 4
	mmol/l	3.66	3.66							± 0.25
HT HCT	%	16.1	16.0							± 2.0
	l/l	0.161	0.160							± 0.020
VGM MCV	µm ³ ; fl	70	70							± 4
TGMH MCH	pg	25.7	25.8							± 2.0
	fmol	1.59	1.60							± 0.12
CCMH MCHC	g/dl	36.6	36.8							± 3.0
	g/l	366	368							± 30
	mmol/l	22.76	22.86							± 1.86
IDR RDW	%	N/A	13.6							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	72	80							± 20
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	1.02							± 0.40
	%	N/A	56.5							± 8.0
MONO	#	N/A	0.20							± 0.20
	%	N/A	9.5							± 6.0
GRANULO	#	N/A	0.61							± 0.40
	%	N/A	34.0							± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 114 N
 Rev 1

CONTROL N

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.2							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.62							± 0.18
HB HGB	g/dl	13.4	13.4							± 0.5
	g/l	134	134							± 5
	mmol/l	8.32	8.32							± 0.31
HT HCT	%	37.3	37.4							± 2.5
	l/l	0.373	0.374							± 0.025
VGM MCV	µm ³ ; fl	81	81							± 4
TGMH MCH	pg	29.1	29.0							± 2.0
	fmol	1.81	1.80							± 0.12
CCMH MCHC	g/dl	36.0	35.8							± 3.0
	g/l	360	358							± 30
	mmol/l	22.33	22.24							± 1.86
IDR RDW	%	N/A	13.5							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	246	255							± 40
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	2.02							± 0.40
	%	N/A	28.0							± 5.0
MONO	#	N/A	0.54							± 0.30
	%	N/A	7.5							± 4.0
GRANULO	#	N/A	4.64							± 0.60
	%	N/A	64.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.36



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 114 H
 Rev 1

CONTROL H

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.3	19.6							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.61	5.60							± 0.20
HB HGB	g/dl	17.9	17.9							± 0.6
	g/l	179	179							± 6
	mmol/l	11.12	11.12							± 0.37
HT HCT	%	49.9	49.8							± 3.0
	l/l	0.499	0.498							± 0.030
VGM MCV	µm ³ ; fl	89	89							± 4
TGMH MCH	pg	31.9	32.0							± 2.0
	fmol	1.98	1.98							± 0.12
CCMH MCHC	g/dl	35.9	35.9							± 3.0
	g/l	359	359							± 30
	mmol/l	22.26	22.30							± 1.86
IDR RDW	%	N/A	13.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	535	525							± 45
VPM MPV	µm ³ ; fl	N/A	7.9							± 2.0
LYMPHO	#	N/A	2.35							± 1.00
	%	N/A	12.0							± 5.0
MONO	#	N/A	0.80							± 0.80
	%	N/A	4.0							± 4.0
GRANULO	#	N/A	16.46							± 1.40
	%	N/A	84.0							± 7.0