

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 074 L
 Rev 1

CONTROL L

(Exp.)

2014-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.39	2.37							± 0.15
HB HGB	g/dl	5.9	5.9							± 0.4
	g/l	59	59							± 4
	mmol/l	3.66	3.66							± 0.25
HT HCT	%	16.5	16.4							± 2.0
	l/l	0.165	0.164							± 0.020
VGM MCV	µm ³ ; fl	69	69							± 5
TGMH MCH	pg	24.7	24.9							± 2.0
	fmol	1.53	1.55							± 0.12
CCMH MCHC	g/dl	35.8	36.1							± 3.0
	g/l	358	361							± 30
	mmol/l	22.22	22.41							± 1.86
IDR RDW	%	N/A	13.7							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	78	85							± 20
VPM MPV	µm ³ ; fl	N/A	8.3							± 2.0
LYMPHO	#	N/A	1.00							± 0.40
	%	N/A	53.0							± 7.0
MONO	#	N/A	0.23							± 0.20
	%	N/A	12.0							± 6.0
GRANULO	#	N/A	0.67							± 0.40
	%	N/A	35.0							± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 074 N
 Rev 1

CONTROL N

(Exp.)

2014-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.2	7.3							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.68	4.66							± 0.18
HB HGB	g/dl	13.4	13.4							± 0.5
	g/l	134	134							± 5
	mmol/l	8.32	8.32							± 0.31
HT HCT	%	37.0	37.3							± 2.5
	l/l	0.370	0.373							± 0.025
VGM MCV	µm ³ ; fl	79	80							± 5
TGMH MCH	pg	28.6	28.8							± 2.0
	fmol	1.78	1.79							± 0.12
CCMH MCHC	g/dl	36.2	35.9							± 3.0
	g/l	362	359							± 30
	mmol/l	22.51	22.32							± 1.86
IDR RDW	%	N/A	13.3							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	260	265							± 40
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	1.97							± 0.40
	%	N/A	27.0							± 5.0
MONO	#	N/A	0.69							± 0.30
	%	N/A	9.5							± 4.0
GRANULO	#	N/A	4.64							± 0.60
	%	N/A	63.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.35

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 074 H
 Rev 1

CONTROL H

(Exp.)

2014-10-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.4	19.6							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.60	5.55							± 0.20
HB HGB	g/dl	17.8	17.8							± 0.6
	g/l	178	178							± 6
	mmol/l	11.05	11.05							± 0.37
HT HCT	%	49.3	49.4							± 3.0
	l/l	0.493	0.494							± 0.030
VGM MCV	µm ³ ; fl	88	89							± 5
TGMH MCH	pg	31.8	32.1							± 2.0
	fmol	1.97	1.99							± 0.12
CCMH MCHC	g/dl	36.1	36.0							± 3.0
	g/l	361	360							± 30
	mmol/l	22.43	22.38							± 1.86
IDR RDW	%	N/A	12.8							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	525	515							± 45
VPM MPV	µm ³ ; fl	N/A	8.0							± 2.0
LYMPHO	#	N/A	2.25							± 1.00
	%	N/A	11.5							± 5.0
MONO	#	N/A	1.08							± 0.80
	%	N/A	5.5							± 4.0
GRANULO	#	N/A	16.27							± 1.40
	%	N/A	83.0							± 7.0