



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 054 L
 Rev 1

CONTROL L



2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.8							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.33	2.32							± 0.15
HB HGB	g/dl	5.9	5.9							± 0.4
	g/l	59	59							± 4
	mmol/l	3.66	3.66							± 0.25
HT HCT	%	16.3	16.5							± 2.0
	l/l	0.163	0.165							± 0.020
VGM MCV	µm ³ ; fl	70	71							± 5
TGMH MCH	pg	25.3	25.4							± 2.0
	fmol	1.57	1.58							± 0.12
CCMH MCHC	g/dl	36.2	35.8							± 3.0
	g/l	362	358							± 30
	mmol/l	22.46	22.24							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	68	77							± 20
VPM MPV	µm ³ ; fl	N/A	8.5							± 2.0
LYMPHO	#	N/A	0.97							± 0.40
	%	N/A	54.0							± 7.0
MONO	#	N/A	0.20							± 0.17
	%	N/A	11.0							± 6.0
GRANULO	#	N/A	0.63							± 0.40
	%	N/A	35.0							± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 054 N
 Rev 1

CONTROL N

(Exp.)

2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.3							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.52	4.51							± 0.18
HB HGB	g/dl	13.4	13.3							± 0.5
	g/l	134	133							± 5
	mmol/l	8.32	8.26							± 0.31
HT HCT	%	36.6	37.0							± 2.5
	l/l	0.366	0.370							± 0.025
VGM MCV	µm ³ ; fl	81	82							± 5
TGMH MCH	pg	29.6	29.5							± 2.0
	fmol	1.84	1.83							± 0.12
CCMH MCHC	g/dl	36.6	36.0							± 3.0
	g/l	366	360							± 30
	mmol/l	22.73	22.33							± 1.86
IDR RDW	%	N/A	13.5							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	255	260							± 40
VPM MPV	µm ³ ; fl	N/A	7.6							± 2.0
LYMPHO	#	N/A	1.90							± 0.40
	%	N/A	26.0							± 5.0
MONO	#	N/A	0.69							± 0.30
	%	N/A	9.5							± 4.0
GRANULO	#	N/A	4.71							± 0.60
	%	N/A	64.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.35



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT **MX 054** **H**
 Rev 1

CONTROL **H**

(Exp.)

2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.8	20.2							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.58	5.50							± 0.20
HB HGB	g/dl	17.8	17.6							± 0.6
	g/l	178	176							± 6
	mmol/l	11.05	10.93							± 0.37
HT HCT	%	49.7	49.5							± 3.0
	l/l	0.497	0.495							± 0.030
VGM MCV	µm ³ ; fl	89	90							± 5
TGMH MCH	pg	31.9	32.0							± 2.0
	fmol	1.98	1.99							± 0.12
CCMH MCHC	g/dl	35.8	35.6							± 3.0
	g/l	358	356							± 30
	mmol/l	22.26	22.08							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	515	495							± 45
VPM MPV	µm ³ ; fl	N/A	7.7							± 2.0
LYMPHO	#	N/A	2.42							± 1.00
	%	N/A	12.0							± 5.0
MONO	#	N/A	1.01							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.77							± 1.40
	%	N/A	83.0							± 7.0