



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 014 L
 Rev 1

CONTROL L

(Exp.)

2014-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.37	2.36							± 0.15
HB HGB	g/dl	5.9	6.0							± 0.4
	g/l	59	60							± 4
	mmol/l	3.66	3.73							± 0.25
HT HCT	%	16.4	16.5							± 2.0
	l/l	0.164	0.165							± 0.020
VGM MCV	µm ³ ; fl	69	70							± 4
TGMH MCH	pg	24.9	25.4							± 2.0
	fmol	1.55	1.58							± 0.12
CCMH MCHC	g/dl	36.1	36.3							± 3.0
	g/l	361	363							± 30
	mmol/l	22.41	22.55							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	68	70							± 20
VPM MPV	µm ³ ; fl	N/A	8.3							± 2.0
LYMPHO	#	N/A	1.03							± 0.40
	%	N/A	54.0							± 7.0
MONO	#	N/A	0.23							± 0.20
	%	N/A	12.0							± 6.0
GRANULO	#	N/A	0.65							± 0.40
	%	N/A	34.0							± 7.0



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 014 N
 Rev 1

CONTROL N



2014-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		SciI Vet ABC	SciI Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.2							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.64	4.57							± 0.18
HB HGB	g/dl	13.5	13.3							± 0.5
	g/l	135	133							± 5
	mmol/l	8.38	8.26							± 0.31
HT HCT	%	37.1	37.0							± 2.5
	l/l	0.371	0.370							± 0.025
VGM MCV	µm ³ ; fl	80	81							± 4
TGMH MCH	pg	29.1	29.1							± 2.0
	fmol	1.81	1.81							± 0.12
CCMH MCHC	g/dl	36.4	35.9							± 3.0
	g/l	364	359							± 30
	mmol/l	22.58	22.31							± 1.86
IDR RDW	%	N/A	13.2							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	241	235							± 40
VPM MPV	µm ³ ; fl	N/A	8.2							± 2.0
LYMPHO	#	N/A	1.94							± 0.40
	%	N/A	27.0							± 5.0
MONO	#	N/A	0.61							± 0.30
	%	N/A	8.5							± 4.0
GRANULO	#	N/A	4.64							± 0.60
	%	N/A	64.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.34



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 014 H
 Rev 1

CONTROL H

(Exp.)

2014-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.6	19.7							± 2.0
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.75	5.68							± 0.20
HB HGB	g/dl	18.0	17.7							± 0.6
	g/l	180	177							± 6
	mmol/l	11.18	10.99							± 0.37
HT HCT	%	50.0	50.0							± 3.0
	l/l	0.500	0.500							± 0.030
VGM MCV	µm ³ ; fl	87	88							± 4
TGMH MCH	pg	31.3	31.2							± 2.0
	fmol	1.94	1.94							± 0.12
CCMH MCHC	g/dl	36.0	35.4							± 3.0
	g/l	360	354							± 30
	mmol/l	22.34	21.99							± 1.86
IDR RDW	%	N/A	12.5							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	540	520							± 45
VPM MPV	µm ³ ; fl	N/A	8.0							± 2.0
LYMPHO	#	N/A	2.17							± 1.00
	%	N/A	11.0							± 5.0
MONO	#	N/A	1.08							± 0.80
	%	N/A	5.5							± 4.0
GRANULO	#	N/A	16.45							± 1.40
	%	N/A	83.5							± 7.0