



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 113 L
 Rev 1

CONTROL L



2014-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		SciI Vet ABC	SciI Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.40							± 0.15
HB HGB	g/dl	5.9	6.0							± 0.4
	g/l	59	60							± 4
	mmol/l	3.66	3.73							± 0.25
HT HCT	%	16.6	16.6							± 2.0
	l/l	0.166	0.166							± 0.020
VGM MCV	µm ³ ; fl	69	69							± 4
TGMH MCH	pg	24.6	25.0							± 2.0
	fmol	1.53	1.55							± 0.12
CCMH MCHC	g/dl	35.6	36.2							± 3.0
	g/l	356	362							± 30
	mmol/l	22.13	22.50							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	70	75							± 20
VPM MPV	µm ³ ; fl	N/A	8.0							± 2.0
LYMPHO	#	N/A	1.06							± 0.40
	%	N/A	56.0							± 7.0
MONO	#	N/A	0.21							± 0.20
	%	N/A	11.0							± 6.0
GRANULO	#	N/A	0.63							± 0.40
	%	N/A	33.0							± 7.0



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 113 N
 Rev 1

CONTROL N

(Exp.)

2014-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.2							± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.58	4.53							± 0.18
HB HGB	g/dl	13.4	13.3							± 0.5
	g/l	134	133							± 5
	mmol/l	8.32	8.26							± 0.31
HT HCT	%	37.1	36.7							± 2.5
	l/l	0.371	0.367							± 0.025
VGM MCV	µm ³ ; fl	81	81							± 4
TGMH MCH	pg	29.3	29.4							± 2.0
	fmol	1.82	1.82							± 0.12
CCMH MCHC	g/dl	36.1	36.2							± 3.0
	g/l	361	362							± 30
	mmol/l	22.43	22.51							± 1.86
IDR RDW	%	N/A	13.7							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	255	255							± 40
VPM MPV	µm ³ ; fl	N/A	7.8							± 2.0
LYMPHO	#	N/A	1.94							± 0.40
	%	N/A	27.0							± 5.0
MONO	#	N/A	0.61							± 0.30
	%	N/A	8.5							± 4.0
GRANULO	#	N/A	4.64							± 0.60
	%	N/A	64.5							± 6.0

9930072-A
 Ref: TEMP-0630 Rev.32

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 113 H
 Rev 1

CONTROL H

(Exp.)

2014-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Sci Vet ABC	Sci Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.0	19.3							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.67	5.60							± 0.20
HB HGB	g/dl	18.0	17.7							± 0.6
	g/l	180	177							± 6
	mmol/l	11.18	10.99							± 0.37
HT HCT	%	50.5	49.8							± 3.0
	l/l	0.505	0.498							± 0.030
VGM MCV	µm ³ ; fl	89	89							± 4
TGMH MCH	pg	31.7	31.6							± 2.0
	fmol	1.97	1.96							± 0.12
CCMH MCHC	g/dl	35.7	35.5							± 3.0
	g/l	357	355							± 30
	mmol/l	22.15	22.05							± 1.86
IDR RDW	%	N/A	12.8							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	500	480							± 45
VPM MPV	µm ³ ; fl	N/A	7.8							± 2.0
LYMPHO	#	N/A	2.32							± 1.00
	%	N/A	12.0							± 5.0
MONO	#	N/A	0.97							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.02							± 1.40
	%	N/A	83.0							± 7.0