

ABX Minotrol 16



LOT MX 092 L
 Rev 3

CONTROL L

(Exp.) 2012-12-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		PENTRA	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		DX NEXUS DF NEXUS	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60		60 60 C+ ES60	80 XL80	120 - 120 RETIC DF120 - DX120	
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1						2.0	2.1	2.1	± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40						2.40	2.35	2.40	± 0.15
HB HGB	g/dl	6.0						6.0	6.0	6.0	± 0.4
	g/l	60						60	60	60	± 4
	mmol/l	3.73						3.73	3.73	3.73	± 0.25
HT HCT	%	17.0						16.8	17.2	17.0	± 2.0
	l/l	0.170						0.168	0.172	0.170	± 0.020
VGM MCV	µm ³ ; fl	71						70	73	71	± 4
TGMH MCH	pg	25.0						25.0	25.5	25.0	± 2.0
	fmol	1.55						1.55	1.59	1.55	± 0.12
CCMH MCHC	g/dl	35.2						35.7	35.0	35.2	± 3.0
	g/l	352						357	350	352	± 30
	mmol/l	21.87						22.18	21.72	21.87	± 1.86
IDR RDW	%	17.5						15.0	15.0	17.5	± 4.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	77						* 80	* 75	77	± 20 *± 22
VPM MPV	µm ³ ; fl	9.5						9.7	9.7	9.5	± 2.0
LYMPHO	#	N/A						N/A	N/A	N/A	± 0.40
	%	N/A						N/A	N/A	N/A	± 9.0
MONO	#	N/A						N/A	N/A	N/A	± 0.20
	%	N/A						N/A	N/A	N/A	± 7.0
GRANULO	#	N/A						N/A	N/A	N/A	± 0.40
	%	N/A						N/A	N/A	N/A	± 8.0

Ref: TEMP-0830 Rev.17 BACK / VERSO 9930067-A

* For ABX Pentra 60 Pentra 80 instruments only

HORIBA ABX SAS
 Parc Euromédecine
 34184 Montpellier cedex 4 - France
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 092 N
 Rev 3

CONTROL N

(Exp.) 2012-12-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		PENTRA	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		DX NEXUS DF NEXUS	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60		60 60 C+ ES60	80 XL80	120 - 120 RETIC DF120 - DX 120	
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	8.1						7.6	7.4	8.1	± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.67						4.65	4.64	4.67	± 0.18
HB HGB	g/dl	13.5						13.4	13.5	13.5	± 0.5
	g/l	135						134	135	135	± 5
	mmol/l	8.38						8.32	8.38	8.38	± 0.31
HT HCT	%	37.4						37.2	37.1	37.4	± 2.5
	l/l	0.374						0.372	0.371	0.374	± 0.025
VGM MCV	µm ³ ; fl	80						80	80	80	± 4
TGMH MCH	pg	28.9						28.8	29.1	28.9	± 2.0
	fmol	1.80						1.79	1.81	1.80	± 0.12
CCMH MCHC	g/dl	36.1						36.0	36.4	36.1	± 3.0
	g/l	361						360	364	361	± 30
	mmol/l	22.44						22.37	22.58	22.44	± 1.86
IDR RDW	%	16.9						13.2	14.0	16.9	± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	263						* 260	* 265	263	± 40 *± 40
VPM MPV	µm ³ ; fl	8.2						8.3	8.3	8.2	± 2.0
LYMPHO	#	N/A						N/A	N/A	N/A	± 0.50
	%	N/A						N/A	N/A	N/A	± 6.0
MONO	#	N/A						N/A	N/A	N/A	± 0.40
	%	N/A						N/A	N/A	N/A	± 5.0
GRANULO	#	N/A						N/A	N/A	N/A	± 0.60
	%	N/A						N/A	N/A	N/A	± 7.0

* For ABX Pentra 60 Pentra 80 instruments only

Ref: TEMP-0830 Rev.17 BACK / VERSO 9930068-A



ABX Minotrol 16



LOT MX 092 H
 Rev 3

CONTROL H

(Exp.) 2012-12-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lysebio									TOLERANCES TOLERANCE
		PENTRA	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		DX NEXUS DF NEXUS	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60		60 60 C+ ES60	80 XL80	120 - 120 RETIC DF 120 - DX 120	
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	21.7						20.8	20.3	21.7	± 1.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.74						5.63	5.66	5.74	± 0.22
HB HGB	g/dl	17.9						18.1	18.1	17.9	± 0.7
	g/l	179						181	181	179	± 6
	mmol/l	11.12						11.24	11.24	11.12	± 0.37
HT HCT	%	49.4						49.5	49.2	49.4	± 3.0
	l/l	0.494						0.495	0.492	0.494	± 0.030
VGM MCV	µm ³ ; fl	86						88	87	86	± 4
TGMH MCH	pg	31.2						32.1	32.0	31.2	± 2.0
	f mol	1.94						2.00	1.99	1.94	± 0.12
CCMH MCHC	g/dl	36.3						36.5	36.8	36.3	± 3.0
	g/l	363						365	368	363	± 30
	mmol/l	22.52						22.69	22.83	22.52	± 1.86
IDR RDW	%	16.8						12.6	13.5	16.8	± 4.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	508						* 505	* 522	508	± 45 *± 55
VPM MPV	µm ³ ; fl	8.4						8.7	8.6	8.4	± 2.0
LYMPHO	#	N/A						N/A	N/A	N/A	± 1.00
	%	N/A						N/A	N/A	N/A	± 5.0
MONO	#	N/A						N/A	N/A	N/A	± 0.80
	%	N/A						N/A	N/A	N/A	± 4.0
GRANULO	#	N/A						N/A	N/A	N/A	± 1.40
	%	N/A						N/A	N/A	N/A	± 7.0

* For ABX Pentra 60 Pentra 80 instruments only

Ref: TEMP-0830 Rev.17 BACK / VERSO 9930069-A