



LOT		MX 420 L		Rev 1		CONTROL		L		 (Exp.)		2020-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	1.9	2.1	2.1	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.32	2.30	2.34	2.35	2.35						± 0.15	
HB HGB	g/dl	6.0	6.0	6.0	6.0	6.0						± 0.4	
	g/l	60	60	60	60	60						± 4	
	mmol/l	3.73	3.73	3.73	3.73	3.73						± 0.25	
HT HCT	%	16.2	16.1	16.4	16.7	16.5						± 2.0	
	l/l	0.162	0.161	0.164	0.167	0.165						± 0.020	
VGM MCV	µm ³ ; fl	70	70	70	71	70						± 4	
TGMH MCH	pg	25.9	26.1	25.6	25.5	25.5						± 2.0	
	fmol	1.61	1.62	1.59	1.59	1.59						± 0.12	
CCMH MCHC	g/dl	36.9	37.3	36.6	36.0	36.5						± 3.0	
	g/l	369	373	366	360	365						± 30	
	mmol/l	22.94	23.14	22.75	22.33	22.65						± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.5	14.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	77	80	71	73	76						± 20	
VPM MPV	µm ³ ; fl	8.5	8.4	N/A	9.0	8.0						± 2.0	
LYMPHO	#	1.28	1.26	N/A	1.11	1.27						± 0.40	
	%	64.0	66.5	N/A	53.0	63.5						± 8.0	
MONO	#	0.20	0.20	N/A	0.34	0.20						± 0.20	
	%	9.0	10.0	N/A	16.0	9.5						± 6.0	
GRANULO	#	0.54	0.45	N/A	0.65	0.54						± 0.40	
	%	27.0	23.5	N/A	31.0	27.0						± 7.0	

LOT		MX 420 N		CONTROL N		 (Exp.)		2020-02-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60											
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.7		7.1		7.6		7.5		7.6									± 0.8			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.65		4.61		4.68		4.64		4.65									± 0.18			
HB	HGB	g/dl	13.4		13.3		13.4		13.5		13.5									± 0.5			
		g/l	134		133		134		135		135									± 5			
		mmol/l	8.32		8.26		8.32		8.38		8.38									± 0.31			
HT	HCT	%	37.2		36.9		37.4		37.1		37.2									± 2.5			
		l/l	0.372		0.369		0.374		0.371		0.372									± 0.025			
VGM	MCV	µm ³ ; fl	80		80		80		80		80									± 4			
TGMH	MCH	pg	28.8		28.9		28.6		29.1		29.0									± 2.0			
		fmol	1.79		1.79		1.78		1.81		1.80									± 0.12			
CCMH	MCHC	g/dl	36.0		36.1		35.8		36.4		36.3									± 3.0			
		g/l	360		361		358		364		363									± 30			
		mmol/l	22.37		22.40		22.23		22.58		22.54									± 1.86			
IDR	RDW	%	14.0		13.7		N/A		14.5		13.5									± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	257		250		248		255		255									± 40			
VPM	MPV	µm ³ ; fl	7.4		7.5		N/A		8.2		7.2									± 2.0			
LYMPHO	#		2.66		2.45		N/A		2.03		2.58									± 0.40			
	%		34.5		34.5		N/A		27.0		34.0									± 5.0			
MONO	#		0.50		0.57		N/A		0.68		0.53									± 0.30			
	%		6.5		8.0		N/A		9.0		7.0									± 4.0			
GRANULO	#		4.54		4.08		N/A		4.80		4.48									± 0.60			
	%		59.0		57.5		N/A		64.0		59.0									± 6.0			

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 420 H		CONTROL		H		 (Exp.)		2020-02-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
				ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.3	19.0	20.0	20.1	20.4						± 1.6
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.66	5.60	5.64	5.66	5.68						± 0.20
HB	HGB	g/dl	18.0	17.7	17.9	18.2	18.1						± 0.6
		g/l	180	177	179	182	181						± 6
		mmol/l	11.18	10.99	11.12	11.30	11.24						± 0.37
HT	HCT	%	50.4	49.8	50.2	50.4	50.6						± 3.0
		l/l	0.504	0.498	0.502	0.504	0.506						± 0.030
VGM	MCV	µm ³ ; fl	89	89	89	89	89						± 4
TGMH	MCH	pg	31.8	31.6	31.7	32.2	31.9						± 2.0
		fmol	1.97	1.96	1.97	2.00	1.98						± 0.12
CCMH	MCHC	g/dl	35.7	35.5	35.7	36.1	35.8						± 3.0
		g/l	357	355	357	361	358						± 30
		mmol/l	22.19	22.05	22.14	22.44	22.23						± 1.86
IDR	RDW	%	13.5	13.2	N/A	13.5	12.8						± 3.0
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	510	490	510	510	506						± 45
VPM	MPV	µm ³ ; fl	7.6	7.7	N/A	8.5	7.4						± 2.0
LYMPHO	#	3.35	3.04	N/A	2.71	3.26							± 1.00
	%	16.5	16.0	N/A	13.5	16.0							± 5.0
MONO	#	1.22	1.24	N/A	0.80	1.02							± 0.80
	%	6.0	6.5	N/A	4.0	5.0							± 4.0
GRANULO	#	15.73	14.73	N/A	16.58	16.12							± 1.40
	%	77.5	77.5	N/A	82.5	79.0							± 7.0