



LOT		MX 421 L		Rev 1		CONTROL		L		 (Exp.)		2020-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8	1.8	1.9	1.9	1.9						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.35	2.40	2.44	2.40						± 0.15	
HB HGB	g/dl	6.0	6.0	6.0	6.1	6.1						± 0.4	
	g/l	60	60	60	61	61						± 4	
	mmol/l	3.73	3.73	3.73	3.79	3.79						± 0.25	
HT HCT	%	16.6	16.2	16.6	16.8	16.6						± 2.0	
	l/l	0.166	0.162	0.166	0.168	0.166						± 0.020	
VGM MCV	µm ³ ; fl	69	69	69	69	69						± 4	
TGMH MCH	pg	25.0	25.5	25.0	25.0	25.4						± 2.0	
	fmol	1.55	1.59	1.55	1.55	1.58						± 0.12	
CCMH MCHC	g/dl	36.2	37.0	36.2	36.2	36.8						± 3.0	
	g/l	362	370	362	362	368						± 30	
	mmol/l	22.50	22.98	22.50	22.50	22.88						± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.5	13.5						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	75	73	65	72	74						± 20	
VPM MPV	µm ³ ; fl	8.7	8.5	N/A	8.7	7.8						± 2.0	
LYMPHO	#	1.14	1.17	N/A	0.91	1.21						± 0.40	
	%	63.5	65.0	N/A	48.0	63.5						± 8.0	
MONO	#	0.20	0.20	N/A	0.36	0.20						± 0.20	
	%	10.0	10.0	N/A	19.0	10.0						± 6.0	
GRANULO	#	0.48	0.45	N/A	0.63	0.50						± 0.40	
	%	26.5	25.0	N/A	33.0	26.5						± 7.0	



LOT		MX 421 N		Rev 1		CONTROL N		 (Exp.)		2020-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7	7.3	7.6	7.6	7.7					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.55	4.55	4.59	4.58					± 0.18
HB HGB	g/dl	13.4	13.4	13.4	13.5	13.5					± 0.5
	g/l	134	134	134	135	135					± 5
	mmol/l	8.32	8.32	8.32	8.38	8.38					± 0.31
HT HCT	%	37.3	36.9	36.9	37.2	37.1					± 2.5
	l/l	0.373	0.369	0.369	0.372	0.371					± 0.025
VGM MCV	µm ³ ; fl	81	81	81	81	81					± 4
TGMH MCH	pg	29.1	29.5	29.5	29.4	29.5					± 2.0
	fmol	1.81	1.83	1.83	1.83	1.83					± 0.12
CCMH MCHC	g/dl	36.0	36.4	36.4	36.3	36.4					± 3.0
	g/l	360	364	364	363	364					± 30
	mmol/l	22.33	22.58	22.58	22.55	22.60					± 1.86
IDR RDW	%	14.0	13.5	N/A	14.0	13.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	268	260	255	263	265					± 40
VPM MPV	µm ³ ; fl	8.0	8.0	N/A	8.5	7.5					± 2.0
LYMPHO	#	2.70	2.56	N/A	2.05	2.66					± 0.40
	%	35.0	35.0	N/A	27.0	34.5					± 5.0
MONO	#	0.50	0.51	N/A	0.72	0.54					± 0.30
	%	6.5	7.0	N/A	9.5	7.0					± 4.0
GRANULO	#	4.50	4.23	N/A	4.83	4.50					± 0.60
	%	58.5	58.0	N/A	63.5	58.5					± 6.0



LOT		MX 421 H		CONTROL		H		 (Exp.)		2020-04-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE								
		MICROS		MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60														
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.7	19.6	20.2	20.6	20.7													± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.71	5.62	5.68	5.73	5.70													± 0.20	
HB HGB	g/dl	17.8	17.8	17.8	18.1	18.1													± 0.6	
	g/l	178	178	178	181	181													± 6	
	mmol/l	11.05	11.05	11.05	11.24	11.24													± 0.37	
HT HCT	%	50.2	49.5	50.6	51.0	50.7													± 3.0	
	l/l	0.502	0.495	0.506	0.510	0.507													± 0.030	
VGM MCV	µm ³ ; fl	88	88	89	89	89													± 4	
TGMH MCH	pg	31.2	31.7	31.3	31.6	31.8													± 2.0	
	fmol	1.94	1.97	1.95	1.96	1.97													± 0.12	
CCMH MCHC	g/dl	35.4	36.0	35.2	35.5	35.7													± 3.0	
	g/l	354	360	352	355	357													± 30	
	mmol/l	22.00	22.35	21.87	22.04	22.16													± 1.86	
IDR RDW	%	13.5	13.5	N/A	13.5	12.5													± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	572	545	560	580	565													± 45	
VPM MPV	µm ³ ; fl	7.9	8.0	N/A	8.4	7.5													± 2.0	
LYMPHO	#	3.62	3.43	N/A	2.88	3.52													± 1.00	
	%	17.5	17.5	N/A	14.0	17.0													± 5.0	
MONO	#	0.93	1.18	N/A	0.82	0.93													± 0.80	
	%	4.5	6.0	N/A	4.0	4.5													± 4.0	
GRANULO	#	16.15	14.99	N/A	16.89	16.25													± 1.40	
	%	78.0	76.5	N/A	82.0	78.5													± 7.0	