



LOT		MX 431 L		CONTROL		L		 (Exp.)		2021-12-05 (YYYY - MM - DD)											
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse							TOLERANCES TOLERANCE										
				MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		
				ES60 Care ST		CRP 200							60 8P - 60 18P ADVIA 60								
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.2	2.1					2.2										± 0.4		
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27	2.22					2.27										± 0.15		
HB	HGB	g/dl	5.9	5.9					6.0										± 0.4		
		g/l	59	59					60									± 4			
		mmol/l	3.66	3.66					3.73									± 0.25			
HT	HCT	%	16.1	15.8					16.1									± 2.0			
		l/l	0.161	0.158					0.161									± 0.020			
VGM	MCV	µm ³ ; fl	71	71					71									± 4			
TGMH	MCH	pg	26.0	26.6					26.4									± 2.0			
		fmol	1.61	1.65					1.64									± 0.12			
CCMH	MCHC	g/dl	36.6	37.4					37.2									± 3.0			
		g/l	366	374					372									± 30			
		mmol/l	22.73	23.25					23.12									± 1.86			
IDR	RDW	%	14.0	13.5					13.0								± 3.0				
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	79	78					77									± 20			
VPM	MPV	µm ³ ; fl	8.8	9.0					8.5									± 2.0			
LYMPHO	#	1.42	1.39						1.43									± 0.40			
	%	64.5	66.0						65.0									± 8.0			
MONO	#	0.20	0.21						0.20									± 0.20			
	%	9.0	10.0						8.5									± 6.0			
GRANULO	#	0.58	0.50						0.58									± 0.40			
	%	26.5	24.0						26.5									± 7.0			

LOT		MX 431 N		CONTROL		N		 (Exp.)		2021-12-05 (YYYY - MM - DD)				
		Rev 1		ABX Lyse										
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE	
				ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60						
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.8	7.7			7.9						± 0.8	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.55	4.50			4.58						± 0.18	
HB	HGB	g/dl	13.5	13.4			13.6						± 0.5	
		g/l	135	134			136					± 5		
		mmol/l	8.38	8.32			8.45					± 0.31		
HT	HCT	%	37.3	36.9			37.6						± 2.5	
		l/l	0.373	0.369			0.376					± 0.025		
VGM	MCV	µm ³ ; fl	82	82			82						± 4	
TGMH	MCH	pg	29.7	29.8			29.7						± 2.0	
		fmol	1.84	1.85			1.84					± 0.12		
CCMH	MCHC	g/dl	36.2	36.3			36.2						± 3.0	
		g/l	362	363			362					± 30		
		mmol/l	22.47	22.55			22.49					± 1.86		
IDR	RDW	%	13.5	13.0			12.5						± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	266	260			266						± 40	
VPM	MPV	µm ³ ; fl	8.2	8.4			8.0						± 2.0	
LYMPHO	#	2.61	2.66			2.65							± 0.40	
	%	33.5	34.5			33.5							± 5.0	
MONO	#	0.62	0.69			0.59							± 0.30	
	%	8.0	9.0			7.5							± 4.0	
GRANULO	#	4.56	4.35			4.66							± 0.60	
	%	58.5	56.5			59.0							± 6.0	



LOT		MX 431 H		CONTROL		H		 (Exp.)		2021-12-05 (YYYY - MM - DD)									
		Rev 1		ABX Lyse															
PARAMETRES PARAMETERS		UNITES UNITS										TOLERANCES TOLERANCE							
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA	
		ES60 Care ST		CRP 200						60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.0	19.7						20.3									± 1.6
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.80	5.84						5.84									± 0.20
HB	HGB	g/dl	18.0	18.1						18.2									± 0.6
		g/l	180	181						182									± 6
		mmol/l	11.18	11.24						11.30									± 0.37
HT	HCT	%	51.0	51.4						51.4									± 3.0
		l/l	0.510	0.514						0.514									± 0.030
VGM	MCV	µm ³ ; fl	88	88						88									± 4
TGMH	MCH	pg	31.0	31.0						31.2									± 2.0
		fmol	1.93	1.92						1.94									± 0.12
CCMH	MCHC	g/dl	35.3	35.2						35.4									± 3.0
		g/l	353	352						354									± 30
		mmol/l	21.90	21.87						21.99									± 1.86
IDR	RDW	%	13.5	13.5						12.5									± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	508	495						497									± 45
VPM	MPV	µm ³ ; fl	7.6	7.8						7.8									± 2.0
LYMPHO	#		3.10	3.05						3.05									± 1.00
	%		15.5	15.5						15.0									± 5.0
MONO	#		0.90	0.99						0.91									± 0.80
	%		4.5	5.0						4.5									± 4.0
GRANULO	#		16.00	15.66						16.34									± 1.40
	%		80.0	79.5						80.5									± 7.0