



LOT		MX 412 L		Rev 1		CONTROL		L		 (Exp.)		2018-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.1	2.2	2.1	2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27	2.28	2.29	2.30	2.27						± 0.15	
HB HGB	g/dl	6.0	6.0	6.0	6.0	6.0						± 0.4	
	g/l	60	60	60	60	60						± 4	
	mmol/l	3.73	3.73	3.73	3.73	3.73						± 0.25	
HT HCT	%	16.1	16.2	16.3	16.3	16.1						± 2.0	
	l/l	0.161	0.162	0.163	0.163	0.161						± 0.020	
VGM MCV	µm ³ ; fl	71	71	71	71	71						± 4	
TGMH MCH	pg	26.4	26.3	26.2	26.1	26.4						± 2.0	
	fmol	1.64	1.63	1.63	1.62	1.64						± 0.12	
CCMH MCHC	g/dl	37.2	37.1	36.9	36.7	37.2						± 3.0	
	g/l	372	371	369	367	372						± 30	
	mmol/l	23.12	23.02	22.92	22.82	23.12						± 1.86	
IDR RDW	%	14.0	13.5	N/A	14.0	13.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	80	71	80	79						± 20	
VPM MPV	µm ³ ; fl	8.7	8.3	N/A	8.7	7.8						± 2.0	
LYMPHO	#	1.29	1.37	N/A	1.05	1.32						± 0.40	
	%	64.5	65.0	N/A	50.0	63.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.34	0.21						± 0.20	
	%	8.5	9.0	N/A	16.0	10.0						± 6.0	
GRANULO	#	0.54	0.55	N/A	0.71	0.57						± 0.40	
	%	27.0	26.0	N/A	34.0	27.0						± 7.0	

LOT		MX 412 N		CONTROL N		 (Exp.)		2018-10-05 (YYYY - MM - DD)			
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.4	7.4	7.4	7.5	7.6				± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.55	4.57	4.62	4.59	4.59				± 0.18	
HB HGB	g/dl	13.3	13.3	13.4	13.4	13.3				± 0.5	
	g/l	133	133	134	134	133				± 5	
	mmol/l	8.26	8.26	8.32	8.32	8.26				± 0.31	
HT HCT	%	36.4	36.6	37.0	36.7	36.7				± 2.5	
	l/l	0.364	0.366	0.370	0.367	0.367				± 0.025	
VGM MCV	µm ³ ; fl	80	80	80	80	80				± 4	
TGMH MCH	pg	29.2	29.1	29.0	29.2	29.0				± 2.0	
	fmol	1.82	1.81	1.80	1.81	1.80				± 0.12	
CCMH MCHC	g/dl	36.5	36.4	36.3	36.5	36.2				± 3.0	
	g/l	365	364	363	365	362				± 30	
	mmol/l	22.69	22.59	22.51	22.66	22.49				± 1.86	
IDR RDW	%	14.0	13.5	N/A	14.0	13.0				± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	260	265	260	270	263				± 40	
VPM MPV	µm ³ ; fl	7.8	7.8	N/A	8.2	7.4				± 2.0	
LYMPHO	#	2.52	2.44	N/A	1.95	2.43				± 0.40	
	%	34.0	33.0	N/A	26.0	32.0				± 5.0	
MONO	#	0.52	0.59	N/A	0.75	0.61				± 0.30	
	%	7.0	8.0	N/A	10.0	8.0				± 4.0	
GRANULO	#	4.37	4.37	N/A	4.80	4.56				± 0.60	
	%	59.0	59.0	N/A	64.0	60.0				± 6.0	

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B



LOT		MX 412 H		Rev 1		CONTROL H		 (Exp.)		2018-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.5	20.4	19.9	20.4	20.5					± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.60	5.60	5.60	5.64	5.62					± 0.20
HB HGB	g/dl	18.0	17.8	17.9	18.1	17.9					± 0.6
	g/l	180	178	179	181	179					± 6
	mmol/l	11.18	11.05	11.12	11.24	11.12					± 0.37
HT HCT	%	49.8	49.8	49.8	50.2	50.0					± 3.0
	l/l	0.498	0.498	0.498	0.502	0.500					± 0.030
VGM MCV	µm ³ ; fl	89	89	89	89	89					± 4
TGMH MCH	pg	32.1	31.8	32.0	32.1	31.9					± 2.0
	fmol	2.00	1.97	1.98	1.99	1.98					± 0.12
CCMH MCHC	g/dl	36.1	35.7	35.9	36.1	35.8					± 3.0
	g/l	361	357	359	361	358					± 30
	mmol/l	22.43	22.18	22.30	22.39	22.22					± 1.86
IDR RDW	%	13.0	13.0	N/A	13.5	12.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	525	525	520	540	520					± 45
VPM MPV	µm ³ ; fl	8.0	7.9	N/A	8.5	7.5					± 2.0
LYMPHO	#	3.08	2.96	N/A	2.35	2.87					± 1.00
	%	15.0	14.5	N/A	11.5	14.0					± 5.0
MONO	#	1.13	1.22	N/A	0.92	1.13					± 0.80
	%	5.5	6.0	N/A	4.5	5.5					± 4.0
GRANULO	#	16.30	16.22	N/A	17.14	16.50					± 1.40
	%	79.5	79.5	N/A	84.0	80.5					± 7.0