



| LOT                      |      | MX 422 L                                               |       | CONTROL         |       | L       |       |  (Exp.) |  | 2020-06-05<br>( YYYY - MM - DD ) |  |                            |  |                         |  |                |  |        |         |        |  |        |  |
|--------------------------|------|--------------------------------------------------------|-------|-----------------|-------|---------|-------|--------------------------------------------------------------------------------------------|--|----------------------------------|--|----------------------------|--|-------------------------|--|----------------|--|--------|---------|--------|--|--------|--|
|                          |      | Rev 1                                                  |       |                 |       |         |       |                                                                                            |  |                                  |  |                            |  |                         |  |                |  |        |         |        |  |        |  |
| PARAMETRES<br>PARAMETERS |      | UNITES<br>UNITS                                        |       | ABX Lyse        |       |         |       |                                                                                            |  |                                  |  |                            |  | TOLERANCES<br>TOLERANCE |  |                |  |        |         |        |  |        |  |
|                          |      |                                                        |       | MICROS          |       | MICROS  |       | MICROS                                                                                     |  | MICROS                           |  | MICROS                     |  |                         |  | ARGOS - HELIOS |  | PENTRA |         | PENTRA |  | PENTRA |  |
|                          |      |                                                        |       | ES60<br>Care ST |       | CRP 200 |       | 8P                                                                                         |  | 18P                              |  | 60 8P - 60 18P<br>ADVIA 60 |  |                         |  |                |  |        |         |        |  |        |  |
| GB                       | WBC  | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 2.0   | 2.0             | 2.0   | 2.0     | 2.0   | 2.0                                                                                        |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.4  |  |        |  |
| GR                       | RBC  | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 2.35  | 2.35            | 2.35  | 2.34    | 2.34  | 2.34                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.15 |  |        |  |
| HB                       | HGB  | g/dl                                                   | 6.1   | 5.9             | 6.0   | 6.1     | 6.1   | 6.1                                                                                        |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.4  |  |        |  |
|                          |      | g/l                                                    | 61    | 59              | 60    | 61      | 61    | 61                                                                                         |  |                                  |  |                            |  |                         |  |                |  |        | ± 4     |        |  |        |  |
|                          |      | mmol/l                                                 | 3.79  | 3.66            | 3.73  | 3.79    | 3.79  | 3.79                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        | ± 0.25  |        |  |        |  |
| HT                       | HCT  | %                                                      | 16.5  | 16.5            | 16.5  | 16.6    | 16.4  | 16.4                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 2.0  |  |        |  |
|                          |      | l/l                                                    | 0.165 | 0.165           | 0.165 | 0.166   | 0.164 | 0.164                                                                                      |  |                                  |  |                            |  |                         |  |                |  |        | ± 0.020 |        |  |        |  |
| VGM                      | MCV  | µm <sup>3</sup> ; fl                                   | 70    | 70              | 70    | 71      | 70    | 70                                                                                         |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 4    |  |        |  |
| TGMH                     | MCH  | pg                                                     | 26.0  | 25.1            | 25.5  | 26.1    | 26.1  | 26.1                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 2.0  |  |        |  |
|                          |      | fmol                                                   | 1.61  | 1.56            | 1.59  | 1.62    | 1.62  | 1.62                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.12 |  |        |  |
| CCMH                     | MCHC | g/dl                                                   | 37.1  | 35.9            | 36.5  | 36.7    | 37.2  | 37.2                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 3.0  |  |        |  |
|                          |      | g/l                                                    | 371   | 359             | 365   | 367     | 372   | 372                                                                                        |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 30   |  |        |  |
|                          |      | mmol/l                                                 | 23.03 | 22.27           | 22.65 | 22.80   | 23.13 | 23.13                                                                                      |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 1.86 |  |        |  |
| IDR                      | RDW  | %                                                      | 14.0  | 14.0            | N/A   | 14.5    | 13.5  | 13.5                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 3.0  |  |        |  |
| PLAQ.                    | PLTS | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 76    | 73              | 69    | 75      | 76    | 76                                                                                         |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 20   |  |        |  |
| VPM                      | MPV  | µm <sup>3</sup> ; fl                                   | 8.4   | 8.7             | N/A   | 8.4     | 8.0   | 8.0                                                                                        |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 2.0  |  |        |  |
| LYMPHO                   | #    |                                                        | 1.29  | 1.30            | N/A   | 0.97    | 1.28  | 1.28                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.40 |  |        |  |
|                          | %    |                                                        | 64.5  | 65.0            | N/A   | 48.5    | 64.0  | 64.0                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 8.0  |  |        |  |
| MONO                     | #    |                                                        | 0.20  | 0.22            | N/A   | 0.37    | 0.21  | 0.21                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.20 |  |        |  |
|                          | %    |                                                        | 10.0  | 11.0            | N/A   | 18.5    | 10.5  | 10.5                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 6.0  |  |        |  |
| GRANULO                  | #    |                                                        | 0.51  | 0.48            | N/A   | 0.66    | 0.51  | 0.51                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 0.40 |  |        |  |
|                          | %    |                                                        | 25.5  | 24.0            | N/A   | 33.0    | 25.5  | 25.5                                                                                       |  |                                  |  |                            |  |                         |  |                |  |        |         | ± 7.0  |  |        |  |

| LOT                      |      | MX 422 N                                               |  | CONTROL N       |  |  (Exp.) |  | 2020-06-05<br>(YYYY - MM - DD) |  |        |  |                            |  |                         |  |                |  |        |  |         |  |
|--------------------------|------|--------------------------------------------------------|--|-----------------|--|--------------------------------------------------------------------------------------------|--|--------------------------------|--|--------|--|----------------------------|--|-------------------------|--|----------------|--|--------|--|---------|--|
|                          |      | Rev 1                                                  |  | ABX Lyse        |  |                                                                                            |  |                                |  |        |  |                            |  |                         |  |                |  |        |  |         |  |
| PARAMETRES<br>PARAMETERS |      | UNITES<br>UNITS                                        |  | ABX Lyse        |  |                                                                                            |  |                                |  |        |  |                            |  | TOLERANCES<br>TOLERANCE |  |                |  |        |  |         |  |
|                          |      |                                                        |  | MICROS          |  | MICROS                                                                                     |  | MICROS                         |  | MICROS |  | MICROS                     |  |                         |  | ARGOS - HELIOS |  | PENTRA |  | PENTRA  |  |
|                          |      |                                                        |  | ES60<br>Care ST |  | CRP 200                                                                                    |  | 8P                             |  | 18P    |  | 60 8P - 60 18P<br>ADVIA 60 |  |                         |  |                |  |        |  |         |  |
| GB                       | WBC  | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  |  | 7.6             |  | 7.5                                                                                        |  | 7.5                            |  | 7.7    |  | 7.8                        |  |                         |  |                |  |        |  | ± 0.8   |  |
| GR                       | RBC  | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l |  | 4.65            |  | 4.68                                                                                       |  | 4.67                           |  | 4.65   |  | 4.68                       |  |                         |  |                |  |        |  | ± 0.18  |  |
| HB                       | HGB  | g/dl                                                   |  | 13.3            |  | 13.2                                                                                       |  | 13.2                           |  | 13.4   |  | 13.4                       |  |                         |  |                |  |        |  | ± 0.5   |  |
|                          |      | g/l                                                    |  | 133             |  | 132                                                                                        |  | 132                            |  | 134    |  | 134                        |  |                         |  |                |  |        |  | ± 5     |  |
|                          |      | mmol/l                                                 |  | 8.26            |  | 8.20                                                                                       |  | 8.20                           |  | 8.32   |  | 8.32                       |  |                         |  |                |  |        |  | ± 0.31  |  |
| HT                       | HCT  | %                                                      |  | 36.7            |  | 37.0                                                                                       |  | 36.9                           |  | 36.7   |  | 37.0                       |  |                         |  |                |  |        |  | ± 2.5   |  |
|                          |      | l/l                                                    |  | 0.367           |  | 0.370                                                                                      |  | 0.369                          |  | 0.367  |  | 0.370                      |  |                         |  |                |  |        |  | ± 0.025 |  |
| VGM                      | MCV  | µm <sup>3</sup> ; fl                                   |  | 79              |  | 79                                                                                         |  | 79                             |  | 79     |  | 79                         |  |                         |  |                |  |        |  | ± 4     |  |
| TGMH                     | MCH  | pg                                                     |  | 28.6            |  | 28.2                                                                                       |  | 28.3                           |  | 28.8   |  | 28.6                       |  |                         |  |                |  |        |  | ± 2.0   |  |
|                          |      | fmol                                                   |  | 1.78            |  | 1.75                                                                                       |  | 1.76                           |  | 1.79   |  | 1.78                       |  |                         |  |                |  |        |  | ± 0.12  |  |
| CCMH                     | MCHC | g/dl                                                   |  | 36.2            |  | 35.7                                                                                       |  | 35.8                           |  | 36.5   |  | 36.2                       |  |                         |  |                |  |        |  | ± 3.0   |  |
|                          |      | g/l                                                    |  | 362             |  | 357                                                                                        |  | 358                            |  | 365    |  | 362                        |  |                         |  |                |  |        |  | ± 30    |  |
|                          |      | mmol/l                                                 |  | 22.48           |  | 22.17                                                                                      |  | 22.22                          |  | 22.65  |  | 22.51                      |  |                         |  |                |  |        |  | ± 1.86  |  |
| IDR                      | RDW  | %                                                      |  | 14.0            |  | 13.5                                                                                       |  | N/A                            |  | 14.0   |  | 13.5                       |  |                         |  |                |  |        |  | ± 3.0   |  |
| PLAQ.                    | PLTS | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  |  | 265             |  | 260                                                                                        |  | 260                            |  | 260    |  | 265                        |  |                         |  |                |  |        |  | ± 40    |  |
| VPM                      | MPV  | µm <sup>3</sup> ; fl                                   |  | 7.7             |  | 7.8                                                                                        |  | N/A                            |  | 8.2    |  | 7.4                        |  |                         |  |                |  |        |  | ± 2.0   |  |
| LYMPHO                   |      | #                                                      |  | 2.62            |  | 2.66                                                                                       |  | N/A                            |  | 2.08   |  | 2.65                       |  |                         |  |                |  |        |  | ± 0.40  |  |
|                          |      | %                                                      |  | 34.5            |  | 35.5                                                                                       |  | N/A                            |  | 27.0   |  | 34.0                       |  |                         |  |                |  |        |  | ± 5.0   |  |
| MONO                     |      | #                                                      |  | 0.46            |  | 0.49                                                                                       |  | N/A                            |  | 0.73   |  | 0.51                       |  |                         |  |                |  |        |  | ± 0.30  |  |
|                          |      | %                                                      |  | 6.0             |  | 6.5                                                                                        |  | N/A                            |  | 9.5    |  | 6.5                        |  |                         |  |                |  |        |  | ± 4.0   |  |
| GRANULO                  |      | #                                                      |  | 4.52            |  | 4.35                                                                                       |  | N/A                            |  | 4.89   |  | 4.64                       |  |                         |  |                |  |        |  | ± 0.60  |  |
|                          |      | %                                                      |  | 59.5            |  | 58.0                                                                                       |  | N/A                            |  | 63.5   |  | 59.5                       |  |                         |  |                |  |        |  | ± 6.0   |  |

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



| LOT                      |      | MX 422 H                                               |  | CONTROL         |  | H       |  |  (Exp.) |  | 2020-06-05<br>( YYYY - MM - DD ) |  |                         |  |                            |  |                |  |        |  |         |  |        |  |
|--------------------------|------|--------------------------------------------------------|--|-----------------|--|---------|--|--------------------------------------------------------------------------------------------|--|----------------------------------|--|-------------------------|--|----------------------------|--|----------------|--|--------|--|---------|--|--------|--|
|                          |      | Rev 1                                                  |  |                 |  |         |  |                                                                                            |  |                                  |  |                         |  |                            |  |                |  |        |  |         |  |        |  |
| PARAMETRES<br>PARAMETERS |      | UNITES<br>UNITS                                        |  | ABX Lyse        |  |         |  |                                                                                            |  |                                  |  | TOLERANCES<br>TOLERANCE |  |                            |  |                |  |        |  |         |  |        |  |
|                          |      |                                                        |  | MICROS          |  | MICROS  |  | MICROS                                                                                     |  | MICROS                           |  |                         |  | MICROS                     |  | ARGOS - HELIOS |  | PENTRA |  | PENTRA  |  | PENTRA |  |
|                          |      |                                                        |  | ES60<br>Care ST |  | CRP 200 |  | 8P                                                                                         |  | 18P                              |  |                         |  | 60 8P - 60 18P<br>ADVIA 60 |  |                |  |        |  |         |  |        |  |
| GB                       | WBC  | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  |  | 20.8            |  | 20.3    |  | 20.1                                                                                       |  | 20.5                             |  | 21.0                    |  |                            |  |                |  |        |  | ± 1.6   |  |        |  |
| GR                       | RBC  | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l |  | 5.72            |  | 5.67    |  | 5.67                                                                                       |  | 5.70                             |  | 5.69                    |  |                            |  |                |  |        |  | ± 0.20  |  |        |  |
| HB                       | HGB  | g/dl                                                   |  | 18.1            |  | 17.9    |  | 17.9                                                                                       |  | 18.2                             |  | 18.2                    |  |                            |  |                |  |        |  | ± 0.6   |  |        |  |
|                          |      | g/l                                                    |  | 181             |  | 179     |  | 179                                                                                        |  | 182                              |  | 182                     |  |                            |  |                |  |        |  | ± 6     |  |        |  |
|                          |      | mmol/l                                                 |  | 11.24           |  | 11.12   |  | 11.12                                                                                      |  | 11.30                            |  | 11.30                   |  |                            |  |                |  |        |  | ± 0.37  |  |        |  |
| HT                       | HCT  | %                                                      |  | 50.3            |  | 49.9    |  | 49.9                                                                                       |  | 50.2                             |  | 50.6                    |  |                            |  |                |  |        |  | ± 3.0   |  |        |  |
|                          |      | l/l                                                    |  | 0.503           |  | 0.499   |  | 0.499                                                                                      |  | 0.502                            |  | 0.506                   |  |                            |  |                |  |        |  | ± 0.030 |  |        |  |
| VGM                      | MCV  | µm <sup>3</sup> ; fl                                   |  | 88              |  | 88      |  | 88                                                                                         |  | 88                               |  | 89                      |  |                            |  |                |  |        |  | ± 4     |  |        |  |
| TGMH                     | MCH  | pg                                                     |  | 31.6            |  | 31.6    |  | 31.6                                                                                       |  | 31.9                             |  | 32.0                    |  |                            |  |                |  |        |  | ± 2.0   |  |        |  |
|                          |      | fmol                                                   |  | 1.97            |  | 1.96    |  | 1.96                                                                                       |  | 1.98                             |  | 1.99                    |  |                            |  |                |  |        |  | ± 0.12  |  |        |  |
| CCMH                     | MCHC | g/dl                                                   |  | 36.0            |  | 35.9    |  | 35.9                                                                                       |  | 36.3                             |  | 35.9                    |  |                            |  |                |  |        |  | ± 3.0   |  |        |  |
|                          |      | g/l                                                    |  | 360             |  | 359     |  | 359                                                                                        |  | 363                              |  | 359                     |  |                            |  |                |  |        |  | ± 30    |  |        |  |
|                          |      | mmol/l                                                 |  | 22.33           |  | 22.28   |  | 22.28                                                                                      |  | 22.53                            |  | 22.32                   |  |                            |  |                |  |        |  | ± 1.86  |  |        |  |
| IDR                      | RDW  | %                                                      |  | 13.0            |  | 12.5    |  | N/A                                                                                        |  | 12.5                             |  | 12.0                    |  |                            |  |                |  |        |  | ± 3.0   |  |        |  |
| PLAQ.                    | PLTS | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  |  | 530             |  | 510     |  | 530                                                                                        |  | 530                              |  | 525                     |  |                            |  |                |  |        |  | ± 45    |  |        |  |
| VPM                      | MPV  | µm <sup>3</sup> ; fl                                   |  | 8.0             |  | 8.2     |  | N/A                                                                                        |  | 8.7                              |  | 7.6                     |  |                            |  |                |  |        |  | ± 2.0   |  |        |  |
| LYMPHO                   | #    |                                                        |  | 3.33            |  | 3.35    |  | N/A                                                                                        |  | 2.67                             |  | 3.26                    |  |                            |  |                |  |        |  | ± 1.00  |  |        |  |
|                          | %    |                                                        |  | 16.0            |  | 16.5    |  | N/A                                                                                        |  | 13.0                             |  | 15.5                    |  |                            |  |                |  |        |  | ± 5.0   |  |        |  |
| MONO                     | #    |                                                        |  | 0.83            |  | 0.91    |  | N/A                                                                                        |  | 0.82                             |  | 0.84                    |  |                            |  |                |  |        |  | ± 0.80  |  |        |  |
|                          | %    |                                                        |  | 4.0             |  | 4.5     |  | N/A                                                                                        |  | 4.0                              |  | 4.0                     |  |                            |  |                |  |        |  | ± 4.0   |  |        |  |
| GRANULO                  | #    |                                                        |  | 16.64           |  | 16.04   |  | N/A                                                                                        |  | 17.02                            |  | 16.91                   |  |                            |  |                |  |        |  | ± 1.40  |  |        |  |
|                          | %    |                                                        |  | 80.0            |  | 79.0    |  | N/A                                                                                        |  | 83.0                             |  | 80.5                    |  |                            |  |                |  |        |  | ± 7.0   |  |        |  |