



LOT		MX 432 L		CONTROL		L		 (Exp.)		2022-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0			2.0					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.32	2.31			2.31					± 0.17
HB HGB	g/dl	6.1	6.0			6.1					± 0.5
	g/l	61	60			61					± 5
	mmol/l	3.79	3.73			3.79					± 0.31
HT HCT	%	16.5	16.2			16.4					± 2.0
	l/l	0.165	0.162			0.164					± 0.020
VGM MCV	µm ³ ; fl	71	70			71					± 4
TGMH MCH	pg	26.3	26.0			26.4					± 2.0
	fmol	1.63	1.61			1.64					± 0.12
CCMH MCHC	g/dl	37.0	37.1			37.2					± 3.0
	g/l	370	371			372					± 30
	mmol/l	23.00	23.04			23.10					± 1.86
IDR RDW	%	15.0	14.0			14.0					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	80			78					± 22
VPM MPV	µm ³ ; fl	8.3	8.3			7.6					± 2.0
LYMPHO	#	1.25	1.28			1.27					± 0.40
	%	62.5	64.0			63.5					± 8.0
MONO	#	0.22	0.26			0.22					± 0.20
	%	11.0	13.0			11.0					± 6.0
GRANULO	#	0.53	0.46			0.51					± 0.40
	%	26.5	23.0			25.5					± 7.0

LOT		MX 432 N		CONTROL		N		 (Exp.)		2022-02-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200								60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l		7.5		7.3						7.5								± 0.8			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.54		4.53						4.54								± 0.26			
HB	HGB	g/dl		13.4		13.4						13.5								± 0.7			
		g/l		134		134						135								± 7			
		mmol/l		8.32		8.32						8.38								± 0.43			
HT	HCT	%		37.2		37.1						37.2								± 2.7			
		l/l		0.372		0.371						0.372								± 0.027			
VGM	MCV	µm ³ ; fl		82		82						82								± 4			
TGMH	MCH	pg		29.5		29.6						29.7								± 2.0			
		fmol		1.83		1.84						1.85								± 0.12			
CCMH	MCHC	g/dl		36.0		36.1						36.3								± 3.0			
		g/l		360		361						363								± 30			
		mmol/l		22.35		22.40						22.52								± 1.86			
IDR	RDW	%		13.5		13.5						13.0								± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		255		250						250								± 40			
VPM	MPV	µm ³ ; fl		7.6		7.7						7.4								± 2.0			
LYMPHO		#		2.66		2.63						2.63								± 0.40			
		%		35.5		36.0						35.0								± 5.0			
MONO		#		0.41		0.51						0.38								± 0.30			
		%		5.5		7.0						5.0								± 4.0			
GRANULO		#		4.43		4.16						4.50								± 0.60			
		%		59.0		57.0						60.0								± 6.0			

Ref: TEMP-0830 Rev.41 FRONT / RECTO 130090307



LOT		MX 432 H		Rev 1		CONTROL H		 (Exp.)		2022-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.5	19.7			20.2					± 1.7
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.66	5.63			5.65					± 0.32
HB HGB	g/dl	18.1	18.1			18.2					± 0.9
	g/l	181	181			182					± 9
	mmol/l	11.24	11.24			11.30					± 0.56
HT HCT	%	50.9	50.1			50.9					± 3.4
	l/l	0.509	0.501			0.509					± 0.034
VGM MCV	µm ³ ; fl	90	89			90					± 4
TGMH MCH	pg	32.0	32.1			32.2					± 2.0
	fmol	1.99	2.00			2.00					± 0.12
CCMH MCHC	g/dl	35.5	36.1			35.8					± 3.0
	g/l	355	361			358					± 30
	mmol/l	22.07	22.43			22.23					± 1.86
IDR RDW	%	15.0	14.0			13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	502	485			495					± 65
VPM MPV	µm ³ ; fl	7.7	7.7			7.3					± 2.0
LYMPHO	#	3.08	2.76			2.83					± 1.00
	%	15.0	14.0			14.0					± 5.0
MONO	#	1.23	1.38			1.21					± 0.80
	%	6.0	7.0			6.0					± 4.0
GRANULO	#	16.20	15.56			16.16					± 1.40
	%	79.0	79.0			80.0					± 7.0