



LOT		MX 403 L		Rev 1		CONTROL		L		 (Exp.)		2017-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.9	2.0	2.1	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.40	2.39	2.42	2.39						± 0.15	
HB HGB	g/dl	6.2	6.2	6.2	6.2	6.2						± 0.4	
	g/l	62	62	62	62	62						± 4	
	mmol/l	3.85	3.85	3.85	3.85	3.85						± 0.25	
HT HCT	%	16.6	16.6	17.0	17.2	16.7						± 2.0	
	l/l	0.166	0.166	0.170	0.172	0.167						± 0.020	
VGM MCV	µm ³ ; fl	69	69	71	71	70						± 4	
TGMH MCH	pg	25.8	25.8	25.9	25.6	25.9						± 2.0	
	fmol	1.60	1.60	1.61	1.59	1.61						± 0.12	
CCMH MCHC	g/dl	37.4	37.4	36.5	36.1	37.1						± 3.0	
	g/l	374	374	365	361	371						± 30	
	mmol/l	23.25	23.25	22.69	22.41	23.01						± 1.86	
IDR RDW	%	14.0	13.5	N/A	14.5	13.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	73	70	63	70	73						± 20	
VPM MPV	µm ³ ; fl	9.0	9.0	N/A	9.0	8.2						± 2.0	
LYMPHO	#	1.22	1.29	N/A	1.08	1.28						± 0.40	
	%	64.0	64.5	N/A	54.0	64.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.26	0.20						± 0.20	
	%	9.0	9.5	N/A	13.0	10.0						± 6.0	
GRANULO	#	0.51	0.52	N/A	0.66	0.52						± 0.40	
	%	27.0	26.0	N/A	33.0	26.0						± 7.0	

LOT		MX 403 N		CONTROL N				 (Exp.)		2017-04-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		7.5		7.6		7.5		7.5		7.6									± 0.8		
GR RBC		10 ⁹ /mm ³ ; 10 ¹² /l		4.70		4.68		4.64		4.68		4.66									± 0.18		
HB HGB		g/dl		13.6		13.6		13.6		13.7		13.6									± 0.5		
		g/l		136		136		136		137		136									± 5		
		mmol/l		8.45		8.45		8.45		8.51		8.45									± 0.31		
HT HCT		%		37.6		37.4		37.6		37.4		37.3									± 2.5		
		l/l		0.376		0.374		0.376		0.374		0.373									± 0.025		
VGM MCV		µm ³ ; fl		80		80		81		80		80									± 4		
TGMH MCH		pg		28.9		29.1		29.3		29.3		29.2									± 2.0		
		fmol		1.80		1.80		1.82		1.82		1.81									± 0.12		
CCMH MCHC		g/dl		36.2		36.3		36.2		36.6		36.5									± 3.0		
		g/l		362		363		362		366		365									± 30		
		mmol/l		22.46		22.56		22.47		22.72		22.65									± 1.86		
IDR RDW		%		13.8		13.0		N/A		14.0		12.8									± 3.0		
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		247		252		246		249		249									± 40		
VPM MPV		µm ³ ; fl		8.6		8.8		N/A		9.0		8.0									± 2.0		
LYMPHO		#		2.63		2.70		N/A		2.36		2.66									± 0.40		
		%		35.0		35.5		N/A		31.5		35.0									± 5.0		
MONO		#		0.45		0.49		N/A		0.45		0.49									± 0.30		
		%		6.0		6.5		N/A		6.0		6.5									± 4.0		
GRANULO		#		4.43		4.41		N/A		4.69		4.45									± 0.60		
		%		59.0		58.0		N/A		62.5		58.5									± 6.0		

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B

LOT		MX 403 H		Rev 1		CONTROL		H		 (Exp.)		2017-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.0	20.0	19.6	20.0	20.0						± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.68	5.62	5.53	5.61	5.62						± 0.20	
HB HGB	g/dl	17.8	17.8	17.5	18.0	17.8						± 0.6	
	g/l	178	178	175	180	178						± 6	
	mmol/l	11.05	11.05	10.87	11.18	11.05						± 0.37	
HT HCT	%	50.0	49.5	48.7	49.4	49.5						± 3.0	
	l/l	0.500	0.495	0.487	0.494	0.495						± 0.030	
VGM MCV	µm ³ ; fl	88	88	88	88	88						± 4	
TGMH MCH	pg	31.3	31.7	31.6	32.1	31.7						± 2.0	
	fmol	1.95	1.97	1.97	1.99	1.97						± 0.12	
CCMH MCHC	g/dl	35.6	36.0	36.0	36.5	36.0						± 3.0	
	g/l	356	360	360	365	360						± 30	
	mmol/l	22.11	22.35	22.33	22.64	22.35						± 1.86	
IDR RDW	%	13.5	12.5	N/A	13.5	12.3						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	510	515	520	520	510						± 45	
VPM MPV	µm ³ ; fl	8.6	8.7	N/A	9.0	8.0						± 2.0	
LYMPHO	#	3.20	3.00	N/A	2.60	3.10						± 1.00	
	%	16.0	15.0	N/A	13.0	15.5						± 5.0	
MONO	#	1.10	1.00	N/A	0.80	0.90						± 0.80	
	%	5.5	5.0	N/A	4.0	4.5						± 4.0	
GRANULO	#	15.70	16.00	N/A	16.60	16.00						± 1.40	
	%	78.5	80.0	N/A	83.0	80.0						± 7.0	