

LOT		MX 413 L		Rev 1		CONTROL		L		 (Exp.)		2018-12-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.1	2.1	2.1	2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27	2.29	2.31	2.30	2.27						± 0.15	
HB HGB	g/dl	5.9	5.9	5.8	5.9	5.9						± 0.4	
	g/l	59	59	58	59	59						± 4	
	mmol/l	3.66	3.66	3.60	3.66	3.66						± 0.25	
HT HCT	%	15.7	15.8	16.2	16.1	15.9						± 2.0	
	l/l	0.157	0.158	0.162	0.161	0.159						± 0.020	
VGM MCV	µm ³ ; fl	69	69	70	70	70						± 4	
TGMH MCH	pg	26.0	25.8	25.1	25.7	26.0						± 2.0	
	fmol	1.61	1.60	1.56	1.59	1.61						± 0.12	
CCMH MCHC	g/dl	37.7	37.3	35.9	36.6	37.1						± 3.0	
	g/l	377	373	359	366	371						± 30	
	mmol/l	23.39	23.19	22.27	22.76	23.06						± 1.86	
IDR RDW	%	14.5	13.5	N/A	14.5	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	77	80	69	76	77						± 20	
VPM MPV	µm ³ ; fl	9.0	8.8	N/A	9.1	8.2						± 2.0	
LYMPHO	#	1.31	1.41	N/A	1.07	1.35						± 0.40	
	%	65.5	67.0	N/A	51.0	64.5						± 8.0	
MONO	#	0.20	0.20	N/A	0.36	0.21						± 0.20	
	%	8.5	9.5	N/A	17.0	10.0						± 6.0	
GRANULO	#	0.52	0.49	N/A	0.67	0.54						± 0.40	
	%	26.0	23.5	N/A	32.0	25.5						± 7.0	

LOT		MX 413 N		CONTROL N				 (Exp.)		2018-12-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		7.5		7.5		7.5		7.5		7.6									± 0.8		
GR RBC		10 ⁹ /mm ³ ; 10 ¹² /l		4.53		4.54		4.57		4.58		4.55									± 0.18		
HB HGB		g/dl		13.3		13.4		13.4		13.5		13.5									± 0.5		
		g/l		133		134		134		135		135									± 5		
		mmol/l		8.26		8.32		8.32		8.38		8.38									± 0.31		
HT HCT		%		36.2		36.8		36.6		36.6		36.9									± 2.5		
		l/l		0.362		0.368		0.366		0.366		0.369									± 0.025		
VGM MCV		µm ³ ; fl		80		81		80		80		81									± 4		
TGMH MCH		pg		29.4		29.5		29.3		29.5		29.7									± 2.0		
		fmol		1.82		1.83		1.82		1.83		1.84									± 0.12		
CCMH MCHC		g/dl		36.7		36.4		36.7		36.8		36.6									± 3.0		
		g/l		367		364		367		368		366									± 30		
		mmol/l		22.79		22.63		22.76		22.88		22.75									± 1.86		
IDR RDW		%		14.5		13.5		N/A		14.5		13.0									± 3.0		
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		262		262		253		260		260									± 40		
VPM MPV		µm ³ ; fl		8.6		8.4		N/A		9.0		7.9									± 2.0		
LYMPHO		#		2.63		2.70		N/A		2.10		2.62									± 0.40		
		%		35.0		36.0		N/A		28.0		34.5									± 5.0		
MONO		#		0.45		0.45		N/A		0.71		0.53									± 0.30		
		%		6.0		6.0		N/A		9.5		7.0									± 4.0		
GRANULO		#		4.43		4.35		N/A		4.69		4.45									± 0.60		
		%		59.0		58.0		N/A		62.5		58.5									± 6.0		

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B



LOT		MX 413 H		CONTROL		H		 (Exp.)		2018-12-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
				ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.3	20.3	20.0	20.1	20.6						± 1.6
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.60	5.56	5.59	5.63	5.55						± 0.20
HB	HGB	g/dl	17.7	17.5	17.7	17.9	17.8						± 0.6
		g/l	177	175	177	179	178						± 6
		mmol/l	10.99	10.87	10.99	11.12	11.05						± 0.37
HT	HCT	%	49.3	49.5	49.8	49.5	49.4						± 3.0
		l/l	0.493	0.495	0.498	0.495	0.494						± 0.030
VGM	MCV	µm ³ ; fl	88	89	89	88	89						± 4
TGMH	MCH	pg	31.6	31.5	31.7	31.8	32.1						± 2.0
		fmol	1.96	1.95	1.97	1.97	1.99						± 0.12
CCMH	MCHC	g/dl	35.9	35.4	35.6	36.1	36.0						± 3.0
		g/l	359	354	356	361	360						± 30
		mmol/l	22.30	21.96	22.09	22.44	22.38						± 1.86
IDR	RDW	%	13.5	13.0	N/A	13.5	12.5						± 3.0
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	507	513	500	515	510						± 55
VPM	MPV	µm ³ ; fl	8.3	8.2	N/A	8.8	7.7						± 2.0
LYMPHO	#	3.45	3.35	N/A	2.61	3.30							± 1.00
	%	17.0	16.5	N/A	13.0	16.0							± 5.0
MONO	#	1.22	1.32	N/A	1.01	1.24							± 0.80
	%	6.0	6.5	N/A	5.0	6.0							± 4.0
GRANULO	#	15.63	15.63	N/A	16.48	16.07							± 1.40
	%	77.0	77.0	N/A	82.0	78.0							± 7.0