



LOT		MX 423 L		CONTROL		L		 (Exp.)		2020-08-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE								
		MICROS		MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
		ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1		2.0		2.1		2.1		2.0								± 0.4		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.34		2.32		2.36		2.37		2.32								± 0.15		
HB HGB	g/dl	6.0		6.0		5.9		6.0		6.0								± 0.4		
	g/l	60		60		59		60		60								± 4		
	mmol/l	3.73		3.73		3.66		3.73		3.73								± 0.25		
HT HCT	%	16.1		16.0		16.3		16.6		16.0								± 2.0		
	l/l	0.161		0.160		0.163		0.166		0.160								± 0.020		
VGM MCV	µm ³ ; fl	69		69		69		70		69								± 4		
TGMH MCH	pg	25.6		25.9		25.0		25.3		25.9								± 2.0		
	fmol	1.59		1.61		1.55		1.57		1.61								± 0.12		
CCMH MCHC	g/dl	37.2		37.5		36.2		36.2		37.5								± 3.0		
	g/l	372		375		362		362		375								± 30		
	mmol/l	23.08		23.28		22.50		22.46		23.28								± 1.86		
IDR RDW	%	14.0		14.0		N/A		14.0		13.5								± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	74		74		70		73		74								± 20		
VPM MPV	µm ³ ; fl	8.4		8.3		N/A		8.6		7.6								± 2.0		
LYMPHO	#	1.34		1.31		N/A		1.03		1.28								± 0.40		
	%	64.0		65.5		N/A		49.0		64.0								± 8.0		
MONO	#	0.22		0.25		N/A		0.44		0.22								± 0.20		
	%	10.5		12.5		N/A		21.0		11.0								± 6.0		
GRANULO	#	0.54		0.44		N/A		0.63		0.50								± 0.40		
	%	25.5		22.0		N/A		30.0		25.0								± 7.0		

LOT		MX 423 N		Rev 1		CONTROL N		 (Exp.)		2020-08-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.8	7.7	7.9	7.8	7.9					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.68	4.65	4.66	4.65	4.66					± 0.18
HB HGB	g/dl	13.6	13.5	13.5	13.5	13.6					± 0.5
	g/l	136	135	135	135	136					± 5
	mmol/l	8.45	8.38	8.38	8.38	8.45					± 0.31
HT HCT	%	37.9	37.7	37.7	37.7	37.7					± 2.5
	l/l	0.379	0.377	0.377	0.377	0.377					± 0.025
VGM MCV	µm ³ ; fl	81	81	81	81	81					± 4
TGMH MCH	pg	29.1	29.0	29.0	29.0	29.2					± 2.0
	fmol	1.80	1.80	1.80	1.80	1.81					± 0.12
CCMH MCHC	g/dl	35.9	35.8	35.8	35.8	36.0					± 3.0
	g/l	359	358	358	358	360					± 30
	mmol/l	22.28	22.26	22.21	22.26	22.37					± 1.86
IDR RDW	%	14.0	14.0	N/A	14.0	13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	263	260	262	260	264					± 40
VPM MPV	µm ³ ; fl	8.1	8.1	N/A	8.5	7.7					± 2.0
LYMPHO	#	2.65	2.62	N/A	2.11	2.69					± 0.40
	%	34.0	34.0	N/A	27.0	34.0					± 5.0
MONO	#	0.59	0.69	N/A	0.78	0.55					± 0.30
	%	7.5	9.0	N/A	10.0	7.0					± 4.0
GRANULO	#	4.56	4.39	N/A	4.91	4.66					± 0.60
	%	58.5	57.0	N/A	63.0	59.0					± 6.0

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 423 H		CONTROL		H		 (Exp.)		2020-08-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE									
		MICROS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		
		ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60										
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	21.0		20.3		20.0		20.5		20.9										± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.70		5.63		5.62		5.67		5.65										± 0.20
HB HGB	g/dl	18.2		17.9		17.7		18.2		18.1										± 0.6
	g/l	182		179		177		182		181										± 6
	mmol/l	11.30		11.12		10.99		11.30		11.24										± 0.37
HT HCT	%	50.7		50.1		50.0		50.5		50.9										± 3.0
	l/l	0.507		0.501		0.500		0.505		0.509										± 0.030
VGM MCV	µm ³ ; fl	89		89		89		89		90										± 4
TGMH MCH	pg	31.9		31.8		31.5		32.1		32.0										± 2.0
	fmol	1.98		1.97		1.96		1.99		1.99										± 0.12
CCMH MCHC	g/dl	35.9		35.7		35.4		36.1		35.6										± 3.0
	g/l	359		357		354		361		356										± 30
	mmol/l	22.28		22.18		21.98		22.40		22.10										± 1.86
IDR RDW	%	13.5		13.5		N/A		13.0		12.5										± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	507		496		505		500		500										± 45
VPM MPV	µm ³ ; fl	7.9		8.0		N/A		8.6		7.5										± 2.0
LYMPHO	#	3.47		3.25		N/A		2.56		3.24										± 1.00
	%	16.5		16.0		N/A		12.5		15.5										± 5.0
MONO	#	1.05		1.22		N/A		0.92		0.94										± 0.80
	%	5.0		6.0		N/A		4.5		4.5										± 4.0
GRANULO	#	16.49		15.83		N/A		17.02		16.72										± 1.40
	%	78.5		78.0		N/A		83.0		80.0										± 7.0