

LOT		MX 433 L		Rev 1		CONTROL		L		 (Exp.)		2022-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1	2.1			2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.30	2.30			2.28						± 0.17	
HB HGB	g/dl	6.2	6.2			6.2						± 0.5	
	g/l	62	62			62						± 5	
	mmol/l	3.85	3.85			3.85						± 0.31	
HT HCT	%	16.6	16.6			16.4						± 2.0	
	l/l	0.166	0.166			0.164						± 0.020	
VGM MCV	µm ³ ; fl	72	72			72						± 4	
TGMH MCH	pg	27.0	27.0			27.2						± 2.0	
	fmol	1.67	1.67			1.69						± 0.12	
CCMH MCHC	g/dl	37.4	37.4			37.8						± 3.0	
	g/l	374	374			378						± 30	
	mmol/l	23.25	23.25			23.45						± 1.86	
IDR RDW	%	15.0	14.0			14.0						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	79	78			80						± 22	
VPM MPV	µm ³ ; fl	9.0	8.5			8.5						± 2.0	
LYMPHO	#	1.37	1.39			1.38						± 0.40	
	%	65.0	66.0			65.5						± 8.0	
MONO	#	0.20	0.22			0.19						± 0.19	
	%	9.5	10.5			9.0						± 6.0	
GRANULO	#	0.54	0.49			0.54						± 0.40	
	%	25.5	23.5			25.5						± 7.0	

LOT		MX 433 N		CONTROL		N		 (Exp.)		2022-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.8	7.6			7.8					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.47	4.47			4.49					± 0.26
HB HGB	g/dl	13.4	13.4			13.5					± 0.7
	g/l	134	134			135					± 7
	mmol/l	8.32	8.32			8.38					± 0.43
HT HCT	%	37.1	37.1			37.3					± 2.7
	l/l	0.371	0.371			0.373					± 0.027
VGM MCV	µm ³ ; fl	83	83			83					± 4
TGMH MCH	pg	30.0	30.0			30.1					± 2.0
	fmol	1.86	1.86			1.87					± 0.12
CCMH MCHC	g/dl	36.1	36.1			36.2					± 3.0
	g/l	361	361			362					± 30
	mmol/l	22.43	22.43			22.50					± 1.86
IDR RDW	%	14.5	14.0			13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	267	264			267					± 40
VPM MPV	µm ³ ; fl	8.1	8.0			7.8					± 2.0
LYMPHO	#	2.61	2.62			2.61					± 0.40
	%	33.5	34.5			33.5					± 5.0
MONO	#	0.55	0.57			0.51					± 0.30
	%	7.0	7.5			6.5					± 4.0
GRANULO	#	4.64	4.41			4.68					± 0.60
	%	59.5	58.0			60.0					± 6.0

Ref: TEMP-0830 Rev.41 FRONT / RECTO 130090307

LOT		MX 433 H		CONTROL		H		 (Exp.)		2022-04-05 (YYYY - MM - DD)									
		Rev 1		ABX Lyse															
PARAMETRES PARAMETERS		UNITES UNITS										TOLERANCES TOLERANCE							
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA	
		ES60 Care ST		CRP 200						60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.5	19.8						20.1									± 1.7
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.60	5.54						5.63									± 0.32
HB	HGB		17.9	17.8						17.9									± 0.9
		g/dl	179	178						179									± 9
		mmol/l	11.12	11.05						11.12									± 0.56
HT	HCT	%	49.8	49.3						50.1									± 3.4
		l/l	0.498	0.493						0.501									± 0.034
VGM	MCV	µm ³ ; fl	89	89						89									± 4
TGMH	MCH	pg	32.0	32.1						31.8									± 2.0
		fmol	1.98	2.00						1.97									± 0.12
CCMH	MCHC	g/dl	35.9	36.1						35.7									± 3.0
		g/l	359	361						357									± 30
		mmol/l	22.30	22.42						22.18									± 1.86
IDR	RDW	%	13.0	13.0						12.0									± 3.0
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	486	475						485									± 65
VPM	MPV	µm ³ ; fl	8.3	8.2						8.0									± 2.0
LYMPHO	#		3.59	3.37						3.32									± 1.00
	%		17.5	17.0						16.5									± 5.0
MONO	#		0.92	0.99						0.90									± 0.80
	%		4.5	5.0						4.5									± 4.0
GRANULO	#		15.99	15.44						15.88									± 1.40
	%		78.0	78.0						79.0									± 7.0