



LOT		MX 414 L		CONTROL		L		 (Exp.)				2019-02-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse											
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE				
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA					
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60									
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.1	2.1	2.1	2.1					± 0.4				
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.34	2.35	2.36	2.38	2.35					± 0.15				
HB HGB	g/dl	6.0	6.1	5.9	6.1	6.1					± 0.4				
	g/l	60	61	59	61	61					± 4				
	mmol/l	3.73	3.79	3.66	3.79	3.79					± 0.25				
HT HCT	%	16.1	16.5	16.5	16.7	16.5					± 2.0				
	l/l	0.161	0.165	0.165	0.167	0.165					± 0.020				
VGM MCV	µm ³ ; fl	69	70	70	70	70					± 4				
TGMH MCH	pg	25.6	26.0	25.0	25.6	26.0					± 2.0				
	fmol	1.59	1.61	1.55	1.59	1.61					± 0.12				
CCMH MCHC	g/dl	37.2	37.1	35.7	36.6	37.1					± 3.0				
	g/l	372	371	357	366	371					± 30				
	mmol/l	23.08	23.03	22.18	22.74	23.03					± 1.86				
IDR RDW	%	14.0	14.0	N/A	14.0	13.5					± 3.0				
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	78	77	70	78	78					± 20				
VPM MPV	µm ³ ; fl	9.1	8.8	N/A	9.1	8.0					± 2.0				
LYMPHO	#	1.30	1.35	N/A	1.07	1.32					± 0.40				
	%	65.0	64.5	N/A	51.0	63.0					± 8.0				
MONO	#	0.20	0.21	N/A	0.36	0.21					± 0.20				
	%	8.5	10.0	N/A	17.0	10.0					± 6.0				
GRANULO	#	0.53	0.54	N/A	0.67	0.57					± 0.40				
	%	26.5	25.5	N/A	32.0	27.0					± 7.0				

LOT		MX 414 N		CONTROL N				 (Exp.)		2019-02-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
				ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.6	7.5	7.5	7.5	7.8						± 0.8
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.62	4.60	4.62	4.61	4.60						± 0.18
HB	HGB	g/dl	13.3	13.3	13.3	13.5	13.4						± 0.5
		g/l	133	133	133	135	134						± 5
		mmol/l	8.26	8.26	8.26	8.38	8.32						± 0.31
HT	HCT	%	37.0	36.8	37.0	37.3	36.8						± 2.5
		l/l	0.370	0.368	0.370	0.373	0.368						± 0.025
VGM	MCV	µm ³ ; fl	80	80	80	81	80						± 4
TGMH	MCH	pg	28.8	28.9	28.8	29.3	29.1						± 2.0
		fmol	1.79	1.80	1.79	1.82	1.81						± 0.12
CCMH	MCHC	g/dl	36.0	36.1	36.0	36.2	36.4						± 3.0
		g/l	360	361	360	362	364						± 30
		mmol/l	22.35	22.44	22.35	22.45	22.61						± 1.86
IDR	RDW	%	14.0	13.5	N/A	14.0	13.0						± 3.0
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260	257	248	256	255						± 40
VPM	MPV	µm ³ ; fl	8.3	8.1	N/A	8.8	7.5						± 2.0
LYMPHO	#	2.66	2.55	N/A	1.99	2.61							± 0.40
	%	35.0	34.0	N/A	26.5	33.5							± 5.0
MONO	#	0.57	0.56	N/A	0.75	0.62							± 0.30
	%	7.5	7.5	N/A	10.0	8.0							± 4.0
GRANULO	#	4.37	4.39	N/A	4.76	4.56							± 0.60
	%	57.5	58.5	N/A	63.5	58.5							± 6.0



LOT		MX 414 H		Rev 1		CONTROL		H		 (Exp.)		2019-02-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.3	20.3	20.0	20.2	20.6						± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.69	5.65	5.65	5.68	5.68						± 0.20	
HB HGB	g/dl	17.9	17.8	17.9	18.3	18.1						± 0.6	
	g/l	179	178	179	183	181						± 6	
	mmol/l	11.12	11.05	11.12	11.36	11.24						± 0.37	
HT HCT	%	50.6	50.3	50.3	50.6	50.6						± 3.0	
	l/l	0.506	0.503	0.503	0.506	0.506						± 0.030	
VGM MCV	µm ³ ; fl	89	89	89	89	89						± 4	
TGMH MCH	pg	31.5	31.5	31.7	32.2	31.9						± 2.0	
	fmol	1.95	1.96	1.97	2.00	1.98						± 0.12	
CCMH MCHC	g/dl	35.3	35.4	35.6	36.2	35.8						± 3.0	
	g/l	353	354	356	362	358						± 30	
	mmol/l	21.95	21.98	22.11	22.48	22.23						± 1.86	
IDR RDW	%	13.5	13.0	N/A	13.5	12.5						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	488	475	475	480	476						± 45	
VPM MPV	µm ³ ; fl	8.2	8.0	N/A	8.8	7.5						± 2.0	
LYMPHO	#	3.35	3.15	N/A	2.32	3.09						± 1.00	
	%	16.5	15.5	N/A	11.5	15.0						± 5.0	
MONO	#	1.32	1.22	N/A	1.01	1.03						± 0.80	
	%	6.5	6.0	N/A	5.0	5.0						± 4.0	
GRANULO	#	15.63	15.94	N/A	16.87	16.48						± 1.40	
	%	77.0	78.5	N/A	83.5	80.0						± 7.0	