



LOT		MX 424 L		Rev 1		CONTROL		L		 (Exp.)		2020-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.0	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.36	2.34	2.36	2.37	2.36						± 0.15	
HB HGB	g/dl	6.0	5.9	6.1	6.0	6.1						± 0.4	
	g/l	60	59	61	60	61						± 4	
	mmol/l	3.73	3.66	3.79	3.73	3.79						± 0.25	
HT HCT	%	16.3	16.1	16.5	16.6	16.3						± 2.0	
	l/l	0.163	0.161	0.165	0.166	0.163						± 0.020	
VGM MCV	µm ³ ; fl	69	69	70	70	69						± 4	
TGMH MCH	pg	25.4	25.2	25.8	25.3	25.8						± 2.0	
	fmol	1.58	1.57	1.61	1.57	1.61						± 0.12	
CCMH MCHC	g/dl	36.8	36.5	36.9	36.2	37.5						± 3.0	
	g/l	368	365	369	362	375						± 30	
	mmol/l	22.88	22.69	22.93	22.46	23.26						± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.5	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	75	74	66	74	74						± 20	
VPM MPV	µm ³ ; fl	8.4	8.0	N/A	8.0	7.6						± 2.0	
LYMPHO	#	1.26	1.26	N/A	0.92	1.25						± 0.40	
	%	63.0	63.0	N/A	46.0	62.5						± 8.0	
MONO	#	0.20	0.22	N/A	0.38	0.20						± 0.20	
	%	9.5	11.0	N/A	19.0	10.0						± 6.0	
GRANULO	#	0.55	0.52	N/A	0.70	0.55						± 0.40	
	%	27.5	26.0	N/A	35.0	27.5						± 7.0	

LOT		MX 424 N		CONTROL N		 (Exp.)		2020-10-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60											
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.8	7.7	7.9	7.7	7.9													± 0.8			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.54	4.56	4.56	4.50	4.56													± 0.18			
HB	HGB	g/dl	13.4	13.3	13.5	13.4	13.5													± 0.5			
		g/l	134	133	135	134	135												± 5				
		mmol/l	8.32	8.26	8.38	8.32	8.38												± 0.31				
HT	HCT	%	36.8	36.9	36.9	36.5	36.9													± 2.5			
		l/l	0.368	0.369	0.369	0.365	0.369													± 0.025			
VGM	MCV	µm ³ ; fl	81	81	81	81	81													± 4			
TGMH	MCH	pg	29.5	29.2	29.6	29.8	29.6													± 2.0			
		fmol	1.83	1.81	1.84	1.85	1.84													± 0.12			
CCMH	MCHC	g/dl	36.4	36.0	36.5	36.8	36.5													± 3.0			
		g/l	364	360	365	368	365													± 30			
		mmol/l	22.63	22.36	22.70	22.83	22.70													± 1.86			
IDR	RDW	%	14.2	13.8	N/A	14.3	13.4													± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	267	261	260	255	262													± 40			
VPM	MPV	µm ³ ; fl	7.7	7.6	N/A	8.0	7.1													± 2.0			
LYMPHO	#		2.65	2.62	N/A	1.96	2.65													± 0.40			
	%		34.0	34.0	N/A	25.5	33.5													± 5.0			
MONO	#		0.55	0.58	N/A	0.77	0.55													± 0.30			
	%		7.0	7.5	N/A	10.0	7.0													± 4.0			
GRANULO	#		4.60	4.50	N/A	4.97	4.70													± 0.60			
	%		59.0	58.5	N/A	64.5	59.5													± 6.0			

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 424 H		CONTROL		H		(Exp.)				2020-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.8	20.1	20.5	20.1	20.7					± 1.6		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.63	5.62	5.64	5.58	5.65					± 0.20		
HB HGB	g/dl	17.9	17.8	17.9	18.1	18.2					± 0.6		
	g/l	179	178	179	181	182					± 6		
	mmol/l	11.12	11.05	11.12	11.24	11.30					± 0.37		
HT HCT	%	50.1	50.0	50.2	49.7	50.3					± 3.0		
	l/l	0.501	0.500	0.502	0.497	0.503					± 0.030		
VGM MCV	µm ³ ; fl	89	89	89	89	89					± 4		
TGMH MCH	pg	31.8	31.7	31.7	32.4	32.2					± 2.0		
	fmol	1.97	1.97	1.97	2.01	2.00					± 0.12		
CCMH MCHC	g/dl	35.7	35.6	35.7	36.4	36.2					± 3.0		
	g/l	357	356	357	364	362					± 30		
	mmol/l	22.18	22.10	22.14	22.63	22.48					± 1.86		
IDR RDW	%	13.2	13.2	N/A	13.0	12.2					± 3.0		
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	523	510	530	510	518					± 45		
VPM MPV	µm ³ ; fl	7.4	7.4	N/A	7.9	7.0					± 2.0		
LYMPHO	#	3.54	3.32	N/A	2.61	3.21					± 1.00		
	%	17.0	16.5	N/A	13.0	15.5					± 5.0		
MONO	#	0.94	1.01	N/A	0.90	0.93					± 0.80		
	%	4.5	5.0	N/A	4.5	4.5					± 4.0		
GRANULO	#	16.33	15.78	N/A	16.58	16.56					± 1.40		
	%	78.5	78.5	N/A	82.5	80.0					± 7.0		