

LOT		MX 434 L		Rev 1		CONTROL		L		 (Exp.)		2022-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0			2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.23	2.26			2.23						± 0.17	
HB HGB	g/dl	6.0	6.1			6.0						± 0.5	
	g/l	60	61			60						± 5	
	mmol/l	3.73	3.79			3.73						± 0.31	
HT HCT	%	16.1	16.3			16.1						± 2.0	
	l/l	0.161	0.163			0.161						± 0.020	
VGM MCV	µm ³ ; fl	72	72			72						± 4	
TGMH MCH	pg	26.9	27.0			26.9						± 2.0	
	fmol	1.67	1.68			1.67						± 0.12	
CCMH MCHC	g/dl	37.4	37.5			37.4						± 3.0	
	g/l	374	375			374						± 30	
	mmol/l	23.21	23.28			23.21						± 1.86	
IDR RDW	%	15.0	13.5			13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	79			78						± 22	
VPM MPV	µm ³ ; fl	8.7	8.5			8.3						± 2.0	
LYMPHO	#	1.23	1.26			1.26						± 0.40	
	%	61.5	63.0			63.0						± 8.0	
MONO	#	0.23	0.25			0.22						± 0.20	
	%	11.5	12.5			11.0						± 6.0	
GRANULO	#	0.54	0.49			0.52						± 0.40	
	%	27.0	24.5			26.0						± 7.0	

LOT		MX 434 N		CONTROL		N		 (Exp.)		2022-06-05 (YYYY - MM - DD)			
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.8	7.8			7.8					± 0.8		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.55	4.61			4.59					± 0.26		
HB HGB	g/dl	13.4	13.6			13.6					± 0.7		
	g/l	134	136			136					± 7		
	mmol/l	8.32	8.45			8.45					± 0.43		
HT HCT	%	37.3	37.3			37.6					± 2.7		
	l/l	0.373	0.373			0.376					± 0.027		
VGM MCV	µm ³ ; fl	82	81			82					± 4		
TGMH MCH	pg	29.5	29.5			29.6					± 2.0		
	fmol	1.83	1.83			1.84					± 0.12		
CCMH MCHC	g/dl	35.9	36.4			36.1					± 3.0		
	g/l	359	364			361					± 30		
	mmol/l	22.30	22.62			22.44					± 1.86		
IDR RDW	%	14.0	13.0			13.0					± 3.0		
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	259	260			258					± 40		
VPM MPV	µm ³ ; fl	8.4	7.9			7.8					± 2.0		
LYMPHO	#	2.57	2.65			2.57					± 0.40		
	%	33.0	34.0			33.0					± 5.0		
MONO	#	0.55	0.62			0.51					± 0.30		
	%	7.0	8.0			6.5					± 4.0		
GRANULO	#	4.68	4.52			4.72					± 0.60		
	%	60.0	58.0			60.5					± 6.0		

Rel: TEMP-0830 Rev.41 FRONT / RECTO 1300090307

LOT		MX 434 H		CONTROL		H		 (Exp.)		2022-06-05 (YYYY - MM - DD)								
		Rev 1		ABX Lyse														
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE						
		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS			PENTRA		PENTRA		PENTRA	
		ES60 Care ST		CRP 200				60 8P - 60 18P ADVIA 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.3		20.2				20.5								± 1.7		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.64		5.66				5.70								± 0.32		
HB HGB	g/dl	18.2		18.3				18.3								± 0.9		
	g/l	182		183				183								± 9		
	mmol/l	11.30		11.36				11.36								± 0.56		
HT HCT	%	50.2		50.9				50.9								± 3.4		
	l/l	0.502		0.509				0.509								± 0.034		
VGM MCV	µm ³ ; fl	89		90				89								± 4		
TGMH MCH	pg	32.3		32.3				32.1								± 2.0		
	fmol	2.00		2.01				1.99								± 0.12		
CCMH MCHC	g/dl	36.3		35.9				36.1								± 3.0		
	g/l	363		359				361								± 30		
	mmol/l	22.52		22.31				22.40								± 1.86		
IDR RDW	%	13.0		13.0				12.0								± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	518		515				515								± 65		
VPM MPV	µm ³ ; fl	8.3		7.9				7.8								± 2.0		
LYMPHO	#	3.15		3.23				3.08								± 1.00		
	%	15.5		16.0				15.0								± 5.0		
MONO	#	1.02		1.11				1.03								± 0.80		
	%	5.0		5.5				5.0								± 4.0		
GRANULO	#	16.14		15.86				16.40								± 1.40		
	%	79.5		78.5				80.0								± 7.0		