



LOT		MX 415 L		Rev 1		CONTROL		L		 (Exp.)		2019-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.9	2.0	2.0	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.35	2.36	2.38	2.38	2.37						± 0.15	
HB HGB	g/dl	6.0	6.1	6.1	6.1	6.1						± 0.4	
	g/l	60	61	61	61	61						± 4	
	mmol/l	3.73	3.79	3.79	3.79	3.79						± 0.25	
HT HCT	%	16.2	16.5	16.7	16.7	16.6						± 2.0	
	l/l	0.162	0.165	0.167	0.167	0.166						± 0.020	
VGM MCV	µm ³ ; fl	69	70	70	70	70						± 4	
TGMH MCH	pg	25.5	25.8	25.6	25.6	25.7						± 2.0	
	fmol	1.59	1.61	1.59	1.59	1.60						± 0.12	
CCMH MCHC	g/dl	37.0	36.9	36.6	36.6	36.8						± 3.0	
	g/l	370	369	366	366	368						± 30	
	mmol/l	22.98	22.93	22.74	22.74	22.83						± 1.86	
IDR RDW	%	14.0	14.0	N/A	14.5	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	74	75	66	75	71						± 20	
VPM MPV	µm ³ ; fl	8.5	8.5	N/A	8.7	7.8						± 2.0	
LYMPHO	#	1.20	1.30	N/A	1.00	1.26						± 0.40	
	%	63.0	65.0	N/A	50.0	63.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.36	0.20						± 0.20	
	%	10.0	9.5	N/A	18.0	10.0						± 6.0	
GRANULO	#	0.51	0.51	N/A	0.64	0.54						± 0.40	
	%	27.0	25.5	N/A	32.0	27.0						± 7.0	

LOT		MX 415 N		CONTROL N				 (Exp.)		2019-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.6	7.6	7.6	7.6	7.8					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.71	4.73	4.75	4.72	4.73					± 0.18
HB HGB	g/dl	13.4	13.5	13.5	13.5	13.5					± 0.5
	g/l	134	135	135	135	135					± 5
	mmol/l	8.32	8.38	8.38	8.38	8.38					± 0.31
HT HCT	%	37.2	37.4	37.5	37.3	37.4					± 2.5
	l/l	0.372	0.374	0.375	0.373	0.374					± 0.025
VGM MCV	µm ³ ; fl	79	79	79	79	79					± 4
TGMH MCH	pg	28.5	28.5	28.4	28.6	28.5					± 2.0
	fmol	1.77	1.77	1.76	1.78	1.77					± 0.12
CCMH MCHC	g/dl	36.0	36.1	36.0	36.2	36.1					± 3.0
	g/l	360	361	360	362	361					± 30
	mmol/l	22.36	22.44	22.34	22.48	22.44					± 1.86
IDR RDW	%	14.0	13.5	N/A	14.0	13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	262	265	260	265	260					± 40
VPM MPV	µm ³ ; fl	8.4	8.3	N/A	8.8	7.9					± 2.0
LYMPHO	#	2.58	2.58	N/A	2.01	2.61					± 0.40
	%	34.0	34.0	N/A	26.5	33.5					± 5.0
MONO	#	0.57	0.57	N/A	0.72	0.59					± 0.30
	%	7.5	7.5	N/A	9.5	7.5					± 4.0
GRANULO	#	4.45	4.45	N/A	4.86	4.60					± 0.60
	%	58.5	58.5	N/A	64.0	59.0					± 6.0

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 415 H		CONTROL		H		 (Exp.)		2019-04-05 (YYYY - MM - DD)											
		Rev 1		ABX Lyse																	
PARAMETRES PARAMETERS		UNITES UNITS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE	
				ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.3		20.2		20.3		20.5		20.6								± 1.6	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.70		5.70		5.72		5.75		5.74								± 0.20	
HB	HGB	g/dl		17.9		17.9		18.0		18.3		18.2								± 0.6	
		g/l		179		179		180		183		182								± 6	
		mmol/l		11.12		11.12		11.18		11.36		11.30								± 0.37	
HT	HCT	%		50.2		50.7		50.9		51.2		51.1								± 3.0	
		l/l		0.502		0.507		0.509		0.512		0.511								± 0.030	
VGM	MCV	µm ³ ; fl		88		89		89		89		89								± 4	
TGMH	MCH	pg		31.4		31.4		31.5		31.8		31.7								± 2.0	
		fmol		1.95		1.95		1.95		1.98		1.97								± 0.12	
CCMH	MCHC	g/dl		35.7		35.3		35.4		35.8		35.6								± 3.0	
		g/l		357		353		354		358		356								± 30	
		mmol/l		22.16		21.91		21.96		22.21		22.12								± 1.86	
IDR	RDW	%		13.5		13.0		N/A		13.5		12.5								± 3.0	
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		513		510		525		530		515								± 45	
VPM	MPV	µm ³ ; fl		8.0		7.9		N/A		8.6		7.5								± 2.0	
LYMPHO		#		3.15		3.03		N/A		2.36		3.09								± 1.00	
		%		15.5		15.0		N/A		11.5		15.0								± 5.0	
MONO		#		1.02		1.11		N/A		0.92		1.03								± 0.80	
		%		5.0		5.5		N/A		4.5		5.0								± 4.0	
GRANULO		#		16.14		16.06		N/A		17.22		16.48								± 1.40	
		%		79.5		79.5		N/A		84.0		80.0								± 7.0	