



LOT		MX 425 L		CONTROL		L		(Exp.)				2020-12-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse											
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE		
				ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1	2.1	2.1	2.1	2.1						± 0.4		
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.22	2.23	2.23	2.21	2.22						± 0.15		
HB	HGB	g/dl	5.9	5.9	6.0	5.8	5.9						± 0.4		
		g/l	59	59	60	58	59						± 4		
		mmol/l	3.66	3.66	3.73	3.60	3.66						± 0.25		
HT	HCT	%	16.0	15.8	16.3	15.9	15.8						± 2.0		
		l/l	0.160	0.158	0.163	0.159	0.158						± 0.020		
VGM	MCV	µm ³ ; fl	72	71	73	72	71						± 4		
TGMH	MCH	pg	26.6	26.5	26.9	26.2	26.6						± 2.0		
		fmol	1.65	1.64	1.67	1.63	1.65						± 0.12		
CCMH	MCHC	g/dl	36.9	37.3	36.9	36.5	37.4						± 3.0		
		g/l	369	373	369	365	374						± 30		
		mmol/l	22.92	23.14	22.89	22.64	23.25						± 1.86		
IDR	RDW	%	14.0	13.5		14.1	13.5						± 3.0		
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	78	77	75	75	77						± 20		
VPM	MPV	µm ³ ; fl	8.0	8.0		8.6	7.5						± 2.0		
LYMPHO	#	1.34	1.39		1.00	1.34							± 0.40		
	%	63.7	66.0		47.5	64.0							± 8.0		
MONO	#	0.21	0.23		0.45	0.22							± 0.20		
	%	9.8	11.0		21.5	10.5							± 6.0		
GRANULO	#	0.56	0.48		0.65	0.54							± 0.40		
	%	26.5	23.0		31.0	25.5							± 7.0		

LOT		MX 425 N		CONTROL N		 (Exp.)		2020-12-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60											
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	8.0		7.9		7.7		7.7		8.0									± 0.8			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.66		4.66		4.64		4.62		4.65									± 0.18			
HB	HGB	g/dl	13.5		13.3		13.3		13.4		13.4									± 0.5			
		g/l	135		133		133		134		134									± 5			
		mmol/l	8.38		8.26		8.26		8.32		8.32									± 0.31			
HT	HCT	%	36.8		36.8		36.7		36.5		36.7									± 2.5			
		l/l	0.368		0.368		0.367		0.365		0.367									± 0.025			
VGM	MCV	µm ³ ; fl	79		79		79		79		79									± 4			
TGMH	MCH	pg	29.0		28.5		28.7		29.0		28.8									± 2.0			
		fmol	1.80		1.77		1.78		1.80		1.79									± 0.12			
CCMH	MCHC	g/dl	36.7		36.1		36.3		36.7		36.5									± 3.0			
		g/l	367		361		363		367		365									± 30			
		mmol/l	22.77		22.44		22.53		22.80		22.65									± 1.86			
IDR	RDW	%	13.5		13.5		N/A		13.5		13.0									± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	255		246		250		250		246									± 40			
VPM	MPV	µm ³ ; fl	7.8		8.0		N/A		8.5		7.4									± 2.0			
LYMPHO	#		2.71		2.73		N/A		1.96		2.76									± 0.40			
	%		33.9		34.5		N/A		25.5		34.5									± 5.0			
MONO	#		0.53		0.59		N/A		0.81		0.56									± 0.30			
	%		6.6		7.5		N/A		10.5		7.0									± 4.0			
GRANULO	#		4.76		4.58		N/A		4.93		4.68									± 0.60			
	%		59.5		58.0		N/A		64.0		58.5									± 6.0			

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 425 H		CONTROL		H		 (Exp.)		2020-12-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	21.2	20.2	20.0	20.0	20.9					± 1.6		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.58	5.56	5.54	5.50	5.59					± 0.20		
HB HGB	g/dl	17.8	17.6	17.6	17.7	17.9					± 0.6		
	g/l	178	176	176	177	179					± 6		
	mmol/l	11.05	10.93	10.93	10.99	11.12					± 0.37		
HT HCT	%	49.7	49.5	49.3	49.0	49.8					± 3.0		
	l/l	0.497	0.495	0.493	0.490	0.498					± 0.030		
VGM MCV	µm ³ ; fl	89	89	89	89	89					± 4		
TGMH MCH	pg	31.9	31.7	31.8	32.2	32.0					± 2.0		
	fmol	1.98	1.97	1.97	2.00	1.99					± 0.12		
CCMH MCHC	g/dl	35.8	35.6	35.7	36.2	36.0					± 3.0		
	g/l	358	356	357	362	360					± 30		
	mmol/l	22.26	22.09	22.17	22.45	22.34					± 1.86		
IDR RDW	%	13.0	12.6	N/A	13.0	12.0					± 3.0		
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	515	497	510	510	500					± 45		
VPM MPV	µm ³ ; fl	7.5	7.7	N/A	8.3	7.2					± 2.0		
LYMPHO	#	3.43	3.23	N/A	2.50	3.45					± 1.00		
	%	16.2	16.0	N/A	12.5	16.5					± 5.0		
MONO	#	1.00	1.11	N/A	0.90	0.94					± 0.80		
	%	4.7	5.5	N/A	4.5	4.5					± 4.0		
GRANULO	#	16.77	15.86	N/A	16.60	16.51					± 1.40		
	%	79.1	78.5	N/A	83.0	79.0					± 7.0		