

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 435 L		CONTROL		L		(Exp.)			2022-08-05 (YYYY-MM-DD)	
Rev 1		ABX Lyse										TOLERANCES TOLERANCE
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		
		ES60 Care ST	CRP 200			60 BP - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.9	1.9			1.9						± 0.4
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.22	2.23			2.22						± 0.17
HB HGB	g/dl	5.9	5.9			5.9						± 0.5
	g/l	59	59			59						± 5
	mmol/l	3.66	3.66			3.66						± 0.25
HT HCT	%	16.0	16.1			16.0						± 2.0
	l/l	0.160	0.161			0.160						± 0.020
VGM MCV	µm ³ ; fl	72	72			72						± 4
TGMH MCH	pg	26.6	26.5			26.6						± 2.0
	fmol	1.65	1.64			1.65						± 0.12
CCMH MCHC	g/dl	36.9	36.7			36.9						± 3.0
	g/l	369	367			369						± 30
	mmol/l	22.92	22.82			22.92						± 1.86
IDR RDW	%	15.0	14.0			13.5						± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	82	81			80						± 22
VPM MPV	µm ³ ; fl	9.0	8.7			8.2						± 2.0
LYMPHO	#	1.23	1.24			1.22						± 0.40
	%	64.5	65.5			64.0						± 8.0
MONO	#	0.19	0.20			0.20						± 0.20
	%	10.0	10.5			10.5						± 6.0
GRANULO	#	0.48	0.46			0.48						± 0.40
	%	25.5	24.0			25.5						± 7.0

LOT		MX 435 N		CONTROL		N		(Exp.)			2022-08-05 (YYYY-MM-DD)	
Rev 1		ABX Lyse										TOLERANCES TOLERANCE
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		
		ES60 Care ST	CRP 200			60 BP - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.6	7.4			7.6						± 0.8
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.67	4.67			4.67						± 0.26
HB HGB	g/dl	13.6	13.7			13.7						± 0.7
	g/l	136	137			137						± 7
	mmol/l	8.45	8.51			8.51						± 0.31
HT HCT	%	37.8	37.8			37.8						± 2.7
	l/l	0.378	0.378			0.378						± 0.027
VGM MCV	µm ³ ; fl	81	81			81						± 4
TGMH MCH	pg	29.1	29.3			29.3						± 2.0
	fmol	1.81	1.82			1.82						± 0.12
CCMH MCHC	g/dl	36.0	36.2			36.2						± 3.0
	g/l	360	362			362						± 30
	mmol/l	22.33	22.49			22.49						± 1.86
IDR RDW	%	14.0	13.5			13.0						± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	262	257			259						± 40
VPM MPV	µm ³ ; fl	8.6	8.1			7.8						± 2.0
LYMPHO	#	2.62	2.70			2.66						± 0.40
	%	34.5	36.5			35.0						± 5.0
MONO	#	0.57	0.59			0.53						± 0.30
	%	7.5	8.0			7.0						± 4.0
GRANULO	#	4.41	4.11			4.41						± 0.60
	%	58.0	55.5			58.0						± 6.0

LOT		MX 435 H		CONTROL		H		(Exp.)			2022-08-05 (YYYY-MM-DD)	
Rev 1		ABX Lyse										TOLERANCES TOLERANCE
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		
		ES60 Care ST	CRP 200			60 BP - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.1	19.4			20.1						± 1.7
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.68	5.68			5.70						± 0.32
HB HGB	g/dl	18.1	18.2			18.2						± 0.9
	g/l	181	182			182						± 9
	mmol/l	11.24	11.30			11.30						± 0.37
HT HCT	%	50.6	50.6			51.3						± 3.4
	l/l	0.506	0.506			0.513						± 0.034
VGM MCV	µm ³ ; fl	89	89			90						± 4
TGMH MCH	pg	31.9	32.0			31.9						± 2.0
	fmol	1.98	1.99			1.98						± 0.12
CCMH MCHC	g/dl	35.8	36.0			35.5						± 3.0
	g/l	358	360			355						± 30
	mmol/l	22.23	22.36			22.03						± 1.86
IDR RDW	%	13.5	13.0			12.0						± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	525	505			510						± 65
VPM MPV	µm ³ ; fl	8.0	7.7			7.5						± 2.0
LYMPHO	#	3.12	3.20			3.12						± 1.00
	%	15.5	16.5			15.5						± 5.0
MONO	#	1.11	1.07			1.01						± 0.80
	%	5.5	5.5			5.0						± 4.0
GRANULO	#	15.88	15.13			15.98						± 1.40
	%	79.0	78.0			79.5						± 7.0