



LOT		MX 406 L		Rev 1		CONTROL		L		 (Exp.)		2017-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.1	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.41	2.45	2.43	2.44	2.43						± 0.15	
HB HGB	g/dl	6.3	6.3	6.2	6.2	6.3						± 0.4	
	g/l	63	63	62	62	63						± 4	
	mmol/l	3.91	3.91	3.85	3.85	3.91						± 0.25	
HT HCT	%	16.6	17.2	17.0	17.1	17.0						± 2.0	
	l/l	0.166	0.172	0.170	0.171	0.170						± 0.020	
VGM MCV	µm ³ ; fl	69	70	70	70	70						± 4	
TGMH MCH	pg	26.1	25.7	25.5	25.4	25.9						± 2.0	
	fmol	1.62	1.60	1.58	1.58	1.61						± 0.12	
CCMH MCHC	g/dl	37.9	36.7	36.4	36.3	37.0						± 3.0	
	g/l	379	367	364	363	370						± 30	
	mmol/l	23.53	22.81	22.63	22.54	23.00						± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.5	12.8						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	77	75	70	80	75						± 20	
VPM MPV	µm ³ ; fl	8.8	8.8	N/A	8.5	8.0						± 2.0	
LYMPHO	#	1.31	1.32	N/A	1.13	1.28						± 0.40	
	%	65.5	66.0	N/A	56.5	64.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.22	0.20						± 0.20	
	%	8.0	8.5	N/A	11.0	9.0						± 6.0	
GRANULO	#	0.53	0.51	N/A	0.65	0.54						± 0.40	
	%	26.5	25.5	N/A	32.5	27.0						± 7.0	

LOT		MX 406 N		CONTROL N				 (Exp.)		2017-10-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB WBC		10 ³ /mm ³ ; 10 ⁹ /l		7.4		7.4		7.5		7.3		7.5									± 0.8		
GR RBC		10 ⁶ /mm ³ ; 10 ¹² /l		4.74		4.77		4.75		4.75		4.79									± 0.18		
HB HGB		g/dl		13.4		13.4		13.4		13.4		13.5									± 0.5		
		g/l		134		134		134		134		135									± 5		
		mmol/l		8.32		8.32		8.32		8.32		8.38									± 0.31		
HT HCT		%		37.9		38.2		38.0		38.0		38.3									± 2.5		
		l/l		0.379		0.382		0.380		0.380		0.383									± 0.025		
VGM MCV		µm ³ ; fl		80		80		80		80		80									± 4		
TGMH MCH		pg		28.3		28.1		28.2		28.2		28.2									± 2.0		
		fmol		1.76		1.74		1.75		1.75		1.75									± 0.12		
CCMH MCHC		g/dl		35.3		35.1		35.3		35.3		35.2									± 3.0		
		g/l		353		351		353		353		352									± 30		
		mmol/l		21.94		21.81		21.90		21.90		21.88									± 1.86		
IDR RDW		%		14.0		13.8		N/A		14.2		12.5									± 3.0		
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		251		250		250		253		250									± 40		
VPM MPV		µm ³ ; fl		8.3		8.3		N/A		8.5		7.7									± 2.0		
LYMPHO		#		2.52		2.44		N/A		2.08		2.48									± 0.40		
		%		34.0		33.0		N/A		28.5		33.0									± 5.0		
MONO		#		0.63		0.67		N/A		0.58		0.60									± 0.30		
		%		8.5		9.0		N/A		8.0		8.0									± 4.0		
GRANULO		#		4.26		4.29		N/A		4.64		4.43									± 0.60		
		%		57.5		58.0		N/A		63.5		59.0									± 6.0		

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B

LOT		MX 406 H		CONTROL		H		 (Exp.)		2017-10-05 (YYYY - MM - DD)													
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.1	20.1	20.1	19.9	20.3													± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.78	5.76	5.76	5.78	5.80													± 0.20			
HB	HGB	g/dl	18.0	18.0	17.8	18.0	18.1													± 0.6			
		g/l	180	180	178	180	181												± 6				
		mmol/l	11.18	11.18	11.05	11.18	11.24												± 0.37				
HT	HCT	%	50.9	50.7	50.7	50.9	51.0													± 3.0			
		l/l	0.509	0.507	0.507	0.509	0.510													± 0.030			
VGM	MCV	µm ³ ; fl	88	88	88	88	88													± 4			
TGMH	MCH	pg	31.1	31.3	30.9	31.1	31.2													± 2.0			
		fmol	1.93	1.94	1.92	1.93	1.94													± 0.12			
CCMH	MCHC	g/dl	35.4	35.5	35.1	35.4	35.5													± 3.0			
		g/l	354	355	351	354	355													± 30			
		mmol/l	21.98	22.05	21.81	21.98	22.02													± 1.86			
IDR	RDW	%	13.5	13.5	N/A	13.8	12.0													± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	500	500	500	505	495													± 45			
VPM	MPV	µm ³ ; fl	7.5	7.4	N/A	7.6	7.0													± 2.0			
LYMPHO	#		3.12	3.12	N/A	2.59	3.05													± 1.00			
	%		15.5	15.5	N/A	13.0	15.0													± 5.0			
MONO	#		1.31	1.31	N/A	0.80	1.12													± 0.80			
	%		6.5	6.5	N/A	4.0	5.5													± 4.0			
GRANULO	#		15.68	15.68	N/A	16.52	16.14													± 1.40			
	%		78.0	78.0	N/A	83.0	79.5													± 7.0			