



| LOT                      |                                                        | MX 426 L        |         | Rev 1  |        | CONTROL                    |                | L      |        |  (Exp.) |                         | 2021-02-05<br>( YYYY - MM - DD ) |  |
|--------------------------|--------------------------------------------------------|-----------------|---------|--------|--------|----------------------------|----------------|--------|--------|--------------------------------------------------------------------------------------------|-------------------------|----------------------------------|--|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        | ABX Lyse        |         |        |        |                            |                |        |        |                                                                                            | TOLERANCES<br>TOLERANCE |                                  |  |
|                          |                                                        | MICROS          | MICROS  | MICROS | MICROS | MICROS                     | ARGOS - HELIOS | PENTRA | PENTRA | PENTRA                                                                                     |                         |                                  |  |
|                          |                                                        | ES60<br>Care ST | CRP 200 | 8P     | 18P    | 60 8P - 60 18P<br>ADVIA 60 |                |        |        |                                                                                            |                         |                                  |  |
| GB WBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 2.0             | 2.0     | 2.0    | 2.0    | 2.0                        |                |        |        |                                                                                            |                         | ± 0.4                            |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 2.39            | 2.38    | 2.42   | 2.41   | 2.40                       |                |        |        |                                                                                            |                         | ± 0.15                           |  |
| HB HGB                   | g/dl                                                   | 5.9             | 5.9     | 5.9    | 5.9    | 6.0                        |                |        |        |                                                                                            |                         | ± 0.4                            |  |
|                          | g/l                                                    | 59              | 59      | 59     | 59     | 60                         |                |        |        |                                                                                            |                         | ± 4                              |  |
|                          | mmol/l                                                 | 3.66            | 3.66    | 3.66   | 3.66   | 3.73                       |                |        |        |                                                                                            |                         | ± 0.25                           |  |
| HT HCT                   | %                                                      | 16.7            | 16.7    | 16.9   | 16.6   | 16.8                       |                |        |        |                                                                                            |                         | ± 2.0                            |  |
|                          | l/l                                                    | 0.167           | 0.167   | 0.169  | 0.166  | 0.168                      |                |        |        |                                                                                            |                         | ± 0.020                          |  |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 70              | 70      | 70     | 69     | 70                         |                |        |        |                                                                                            |                         | ± 4                              |  |
| TGMH MCH                 | pg                                                     | 24.7            | 24.8    | 24.4   | 24.5   | 25.0                       |                |        |        |                                                                                            |                         | ± 2.0                            |  |
|                          | fmol                                                   | 1.53            | 1.54    | 1.51   | 1.52   | 1.55                       |                |        |        |                                                                                            |                         | ± 0.12                           |  |
| CCMH MCHC                | g/dl                                                   | 35.3            | 35.4    | 34.8   | 35.5   | 35.7                       |                |        |        |                                                                                            |                         | ± 3.0                            |  |
|                          | g/l                                                    | 353             | 354     | 348    | 355    | 357                        |                |        |        |                                                                                            |                         | ± 30                             |  |
|                          | mmol/l                                                 | 21.90           | 21.99   | 21.63  | 22.03  | 22.18                      |                |        |        |                                                                                            |                         | ± 1.86                           |  |
| IDR RDW                  | %                                                      | 14.5            | 13.8    | N/A    | 14.5   | 13.5                       |                |        |        |                                                                                            |                         | ± 3.0                            |  |
| PLAQ. PLTS               | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 73              | 74      | 70     | 71     | 73                         |                |        |        |                                                                                            |                         | ± 20                             |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 8.1             | 8.2     | N/A    | 8.6    | 7.5                        |                |        |        |                                                                                            |                         | ± 2.0                            |  |
| LYMPHO                   | #                                                      | 1.27            | 1.32    | N/A    | 1.00   | 1.28                       |                |        |        |                                                                                            |                         | ± 0.40                           |  |
|                          | %                                                      | 63.7            | 66.1    | N/A    | 50.0   | 64.0                       |                |        |        |                                                                                            |                         | ± 8.0                            |  |
| MONO                     | #                                                      | 0.20            | 0.20    | N/A    | 0.36   | 0.20                       |                |        |        |                                                                                            |                         | ± 0.20                           |  |
|                          | %                                                      | 8.9             | 9.1     | N/A    | 18.0   | 8.5                        |                |        |        |                                                                                            |                         | ± 6.0                            |  |
| GRANULO                  | #                                                      | 0.55            | 0.50    | N/A    | 0.64   | 0.55                       |                |        |        |                                                                                            |                         | ± 0.40                           |  |
|                          | %                                                      | 27.4            | 24.8    | N/A    | 32.0   | 27.5                       |                |        |        |                                                                                            |                         | ± 7.0                            |  |

| LOT                      |                                                        | MX 426 N        |         | CONTROL |        | N                          |                |  (Exp.) |        | 2021-02-05<br>(YYYY - MM - DD) |                         |
|--------------------------|--------------------------------------------------------|-----------------|---------|---------|--------|----------------------------|----------------|--------------------------------------------------------------------------------------------|--------|--------------------------------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        | ABX Lyse        |         |         |        |                            |                |                                                                                            |        |                                | TOLERANCES<br>TOLERANCE |
|                          |                                                        | MICROS          | MICROS  | MICROS  | MICROS | MICROS                     | ARGOS - HELIOS | PENTRA                                                                                     | PENTRA | PENTRA                         |                         |
|                          |                                                        | ES60<br>Care ST | CRP 200 | 8P      | 18P    | 60 8P - 60 18P<br>ADVIA 60 |                |                                                                                            |        |                                |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 7.5             | 7.3     | 7.4     | 7.3    | 7.6                        |                |                                                                                            |        |                                | ± 0.8                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 4.58            | 4.56    | 4.62    | 4.60   | 4.62                       |                |                                                                                            |        |                                | ± 0.18                  |
| HB HGB                   | g/dl                                                   | 13.5            | 13.5    | 13.5    | 13.6   | 13.6                       |                |                                                                                            |        |                                | ± 0.5                   |
|                          | g/l                                                    | 135             | 135     | 135     | 136    | 136                        |                |                                                                                            |        |                                | ± 5                     |
|                          | mmol/l                                                 | 8.38            | 8.38    | 8.38    | 8.45   | 8.45                       |                |                                                                                            |        |                                | ± 0.31                  |
| HT HCT                   | %                                                      | 38.0            | 37.8    | 37.9    | 37.7   | 38.4                       |                |                                                                                            |        |                                | ± 2.5                   |
|                          | l/l                                                    | 0.380           | 0.378   | 0.379   | 0.377  | 0.384                      |                |                                                                                            |        |                                | ± 0.025                 |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 83              | 83      | 82      | 82     | 83                         |                |                                                                                            |        |                                | ± 4                     |
| TGMH MCH                 | pg                                                     | 29.5            | 29.6    | 29.2    | 29.6   | 29.4                       |                |                                                                                            |        |                                | ± 2.0                   |
|                          | fmol                                                   | 1.83            | 1.84    | 1.81    | 1.84   | 1.83                       |                |                                                                                            |        |                                | ± 0.12                  |
| CCMH MCHC                | g/dl                                                   | 35.5            | 35.7    | 35.6    | 36.1   | 35.5                       |                |                                                                                            |        |                                | ± 3.0                   |
|                          | g/l                                                    | 355             | 357     | 356     | 361    | 355                        |                |                                                                                            |        |                                | ± 30                    |
|                          | mmol/l                                                 | 22.05           | 22.15   | 22.13   | 22.39  | 22.02                      |                |                                                                                            |        |                                | ± 1.86                  |
| IDR RDW                  | %                                                      | 13.8            | 13.5    | N/A     | 13.5   | 13.0                       |                |                                                                                            |        |                                | ± 3.0                   |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 252             | 248     | 245     | 254    | 250                        |                |                                                                                            |        |                                | ± 40                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 7.3             | 7.4     | N/A     | 8.1    | 6.9                        |                |                                                                                            |        |                                | ± 2.0                   |
| LYMPHO                   | #                                                      | 2.70            | 2.65    | N/A     | 2.08   | 2.70                       |                |                                                                                            |        |                                | ± 0.40                  |
|                          | %                                                      | 36.0            | 36.3    | N/A     | 28.5   | 35.5                       |                |                                                                                            |        |                                | ± 5.0                   |
| MONO                     | #                                                      | 0.54            | 0.58    | N/A     | 0.73   | 0.49                       |                |                                                                                            |        |                                | ± 0.30                  |
|                          | %                                                      | 7.2             | 7.9     | N/A     | 10.0   | 6.5                        |                |                                                                                            |        |                                | ± 4.0                   |
| GRANULO                  | #                                                      | 4.26            | 4.07    | N/A     | 4.49   | 4.41                       |                |                                                                                            |        |                                | ± 0.60                  |
|                          | %                                                      | 56.8            | 55.8    | N/A     | 61.5   | 58.0                       |                |                                                                                            |        |                                | ± 6.0                   |

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513

| LOT                      |                                                        | MX 426 H        |         | CONTROL  |       | H                          |  |  (Exp.) |  | 2021-02-05<br>( YYYY - MM - DD ) |                         |                |  |        |  |        |  |        |         |  |
|--------------------------|--------------------------------------------------------|-----------------|---------|----------|-------|----------------------------|--|--------------------------------------------------------------------------------------------|--|----------------------------------|-------------------------|----------------|--|--------|--|--------|--|--------|---------|--|
|                          |                                                        | Rev 1           |         | ABX Lyse |       |                            |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        |         |  |
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        |                 |         |          |       |                            |  |                                                                                            |  |                                  | TOLERANCES<br>TOLERANCE |                |  |        |  |        |  |        |         |  |
|                          |                                                        | MICROS          |         | MICROS   |       | MICROS                     |  | MICROS                                                                                     |  | MICROS                           |                         | ARGOS - HELIOS |  | PENTRA |  | PENTRA |  | PENTRA |         |  |
|                          |                                                        | ES60<br>Care ST | CRP 200 | 8P       | 18P   | 60 8P - 60 18P<br>ADVIA 60 |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        |         |  |
| GB WBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 20.2            | 19.5    | 19.4     | 19.3  | 20.1                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 1.6   |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 5.70            | 5.67    | 5.72     | 5.70  | 5.75                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.20  |  |
| HB HGB                   | g/dl                                                   | 18.0            | 18.0    | 17.9     | 18.2  | 18.2                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.6   |  |
|                          | g/l                                                    | 180             | 180     | 179      | 182   | 182                        |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 6     |  |
|                          | mmol/l                                                 | 11.18           | 11.18   | 11.12    | 11.30 | 11.30                      |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.37  |  |
| HT HCT                   | %                                                      | 51.9            | 51.6    | 51.5     | 50.7  | 52.3                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 3.0   |  |
|                          | l/l                                                    | 0.519           | 0.516   | 0.515    | 0.507 | 0.523                      |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.030 |  |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 91              | 91      | 90       | 89    | 91                         |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 4     |  |
| TGMH MCH                 | pg                                                     | 31.6            | 31.7    | 31.3     | 31.9  | 31.7                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 2.0   |  |
|                          | fmol                                                   | 1.96            | 1.97    | 1.94     | 1.98  | 1.97                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.12  |  |
| CCMH MCHC                | g/dl                                                   | 34.7            | 34.9    | 34.8     | 35.9  | 34.8                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 3.0   |  |
|                          | g/l                                                    | 347             | 349     | 348      | 359   | 348                        |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 30    |  |
|                          | mmol/l                                                 | 21.55           | 21.66   | 21.59    | 22.28 | 21.60                      |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 1.86  |  |
| IDR RDW                  | %                                                      | 13.0            | 12.8    | N/A      | 13.0  | 12.0                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 3.0   |  |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 515             | 505     | 520      | 532   | 520                        |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 45    |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 7.2             | 7.2     | N/A      | 8.0   | 6.9                        |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 2.0   |  |
| LYMPHO                   | #                                                      | 3.37            | 3.22    | N/A      | 2.51  | 3.22                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 1.00  |  |
|                          | %                                                      | 16.7            | 16.5    | N/A      | 13.0  | 16.0                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 5.0   |  |
| MONO                     | #                                                      | 0.93            | 0.99    | N/A      | 0.87  | 0.90                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 0.80  |  |
|                          | %                                                      | 4.6             | 5.1     | N/A      | 4.5   | 4.5                        |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 4.0   |  |
| GRANULO                  | #                                                      | 15.90           | 15.29   | N/A      | 15.92 | 15.98                      |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 1.40  |  |
|                          | %                                                      | 78.7            | 78.4    | N/A      | 82.5  | 79.5                       |  |                                                                                            |  |                                  |                         |                |  |        |  |        |  |        | ± 7.0   |  |