



LOT		MX 427 L		Rev 1		CONTROL		L		 (Exp.)		2021-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.1	2.0	2.1						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.35	2.32	2.41	2.34	2.35						± 0.15	
HB HGB	g/dl	6.0	5.9	6.1	6.0	6.1						± 0.4	
	g/l	60	59	61	60	61						± 4	
	mmol/l	3.73	3.66	3.79	3.73	3.79						± 0.25	
HT HCT	%	16.2	16.0	16.9	16.4	16.5						± 2.0	
	l/l	0.162	0.160	0.169	0.164	0.165						± 0.020	
VGM MCV	µm ³ ; fl	69	69	70	70	70						± 4	
TGMH MCH	pg	25.5	25.4	25.3	25.6	26.0						± 2.0	
	fmol	1.59	1.58	1.57	1.59	1.61						± 0.12	
CCMH MCHC	g/dl	37.0	36.9	36.2	36.6	37.1						± 3.0	
	g/l	370	369	362	366	371						± 30	
	mmol/l	22.98	22.89	22.45	22.75	23.03						± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.0	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	83	83	73	79	81						± 20	
VPM MPV	µm ³ ; fl	8.3	8.1	N/A	8.7	7.6						± 2.0	
LYMPHO	#	1.27	1.34	N/A	0.99	1.34						± 0.40	
	%	63.7	67.0	N/A	49.3	63.9						± 8.0	
MONO	#	0.20	0.19	N/A	0.37	0.20						± 0.20	
	%	9.8	9.5	N/A	18.7	9.5						± 6.0	
GRANULO	#	0.53	0.47	N/A	0.64	0.56						± 0.40	
	%	26.5	23.5	N/A	32.0	26.6						± 7.0	

LOT		MX 427 N		CONTROL		N		 (Exp.)		2021-04-05 (YYYY-MM-DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.7	7.5	7.7	7.6	7.7					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.69	4.67	4.75	4.70	4.71					± 0.18
HB HGB	g/dl	13.7	13.5	13.8	13.8	13.8					± 0.5
	g/l	137	135	138	138	138					± 5
	mmol/l	8.51	8.38	8.57	8.57	8.57					± 0.31
HT HCT	%	38.0	37.8	38.5	38.1	38.2					± 2.5
	l/l	0.380	0.378	0.385	0.381	0.382					± 0.025
VGM MCV	µm ³ ; fl	81	81	81	81	81					± 4
TGMH MCH	pg	29.2	28.9	29.1	29.4	29.3					± 2.0
	fmol	1.81	1.80	1.80	1.82	1.82					± 0.12
CCMH MCHC	g/dl	36.1	35.7	35.9	36.2	36.2					± 3.0
	g/l	361	357	359	362	362					± 30
	mmol/l	22.40	22.16	22.27	22.51	22.46					± 1.86
IDR RDW	%	13.8	13.5	N/A	13.5	13.2					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	280	274	268	273	275					± 40
VPM MPV	µm ³ ; fl	7.9	8.1	N/A	8.5	7.4					± 2.0
LYMPHO	#	2.66	2.67	N/A	2.09	2.63					± 0.40
	%	34.5	35.6	N/A	27.5	34.1					± 5.0
MONO	#	0.59	0.62	N/A	0.77	0.58					± 0.30
	%	7.6	8.2	N/A	10.1	7.5					± 4.0
GRANULO	#	4.46	4.22	N/A	4.74	4.50					± 0.60
	%	57.9	56.2	N/A	62.4	58.4					± 6.0

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 427 H		CONTROL		H		 (Exp.)		2021-04-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE								
		MICROS		MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60														
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	21.1	20.3	20.5	20.7	21.0													± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.72	5.65	5.73	5.65	5.72													± 0.20	
HB HGB	g/dl	18.0	17.7	17.9	18.1	18.1													± 0.6	
	g/l	180	177	179	181	181													± 6	
	mmol/l	11.18	10.99	11.12	11.24	11.24													± 0.37	
HT HCT	%	50.3	49.7	50.4	49.7	50.9													± 3.0	
	l/l	0.503	0.497	0.504	0.497	0.509													± 0.030	
VGM MCV	µm ³ ; fl	88	88	88	88	89													± 4	
TGMH MCH	pg	31.5	31.3	31.2	32.0	31.6													± 2.0	
	fmol	1.95	1.95	1.94	1.99	1.97													± 0.12	
CCMH MCHC	g/dl	35.8	35.6	35.5	36.4	35.6													± 3.0	
	g/l	358	356	355	364	356													± 30	
	mmol/l	22.21	22.11	22.04	22.61	22.08													± 1.86	
IDR RDW	%	12.8	12.8	N/A	12.6	12.0													± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	550	536	545	552	543													± 45	
VPM MPV	µm ³ ; fl	7.7	7.8	N/A	8.4	7.2													± 2.0	
LYMPHO	#	3.44	3.27	N/A	2.63	3.28													± 1.00	
	%	16.3	16.1	N/A	12.7	15.6													± 5.0	
MONO	#	0.99	1.08	N/A	0.89	0.99													± 0.80	
	%	4.7	5.3	N/A	4.3	4.7													± 4.0	
GRANULO	#	16.67	15.96	N/A	17.18	16.74													± 1.40	
	%	79.0	78.6	N/A	83.0	79.7													± 7.0	