



LOT		MX 418 L		Rev 1		CONTROL		L		 (Exp.)		2019-10-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	1.9	2.1	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.32	2.29	2.32	2.30	2.32						± 0.15	
HB HGB	g/dl	6.0	6.0	6.0	6.0	6.1						± 0.4	
	g/l	60	60	60	60	61						± 4	
	mmol/l	3.73	3.73	3.73	3.73	3.79						± 0.25	
HT HCT	%	16.2	16.0	16.2	16.3	16.2						± 2.0	
	l/l	0.162	0.160	0.162	0.163	0.162						± 0.020	
VGM MCV	µm ³ ; fl	70	70	70	71	70						± 4	
TGMH MCH	pg	25.9	26.2	25.9	26.1	26.3						± 2.0	
	fmol	1.61	1.63	1.61	1.62	1.63						± 0.12	
CCMH MCHC	g/dl	36.9	37.4	36.9	36.7	37.6						± 3.0	
	g/l	369	374	369	367	376						± 30	
	mmol/l	22.94	23.24	22.94	22.82	23.33						± 1.86	
IDR RDW	%	14.0	14.0	N/A	14.0	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	76	78	70	75	78						± 20	
VPM MPV	µm ³ ; fl	8.5	8.1	N/A	8.6	7.6						± 2.0	
LYMPHO	#	1.29	1.27	N/A	1.01	1.28						± 0.40	
	%	64.5	67.0	N/A	50.5	64.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.35	0.20						± 0.20	
	%	9.0	9.5	N/A	17.5	9.5						± 6.0	
GRANULO	#	0.53	0.45	N/A	0.64	0.53						± 0.40	
	%	26.5	23.5	N/A	32.0	26.5						± 7.0	



LOT		MX 418 N		CONTROL N				 (Exp.)		2019-10-05 (YYYY - MM - DD)		
		Rev 1		ABX Lyse								
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA		
		ES60 Care ST	CRP 200	8P	18P	60 8P - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7	7.3	7.6	7.6	7.8					± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.55	4.64	4.56	4.61					± 0.18	
HB HGB	g/dl	13.3	13.2	13.3	13.3	13.4					± 0.5	
	g/l	133	132	133	133	134					± 5	
	mmol/l	8.26	8.20	8.26	8.26	8.32					± 0.31	
HT HCT	%	36.8	36.4	37.1	36.5	36.9					± 2.5	
	l/l	0.368	0.364	0.371	0.365	0.369					± 0.025	
VGM MCV	µm ³ ; fl	80	80	80	80	80					± 4	
TGMH MCH	pg	28.9	29.0	28.7	29.2	29.1					± 2.0	
	fmol	1.80	1.80	1.78	1.81	1.81					± 0.12	
CCMH MCHC	g/dl	36.1	36.3	35.8	36.5	36.3					± 3.0	
	g/l	361	363	358	365	363					± 30	
	mmol/l	22.44	22.52	22.25	22.64	22.56					± 1.86	
IDR RDW	%	14.0	14.0	N/A	14.0	13.0					± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	254	253	251	256	253					± 40	
VPM MPV	µm ³ ; fl	8.2	7.9	N/A	8.7	7.5					± 2.0	
LYMPHO	#	2.77	2.66	N/A	2.17	2.73					± 0.40	
	%	36.0	36.5	N/A	28.5	35.0					± 5.0	
MONO	#	0.50	0.55	N/A	0.72	0.55					± 0.30	
	%	6.5	7.5	N/A	9.5	7.0					± 4.0	
GRANULO	#	4.43	4.09	N/A	4.71	4.52					± 0.60	
	%	57.5	56.0	N/A	62.0	58.0					± 6.0	

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 418 H		CONTROL		H		 (Exp.)		2019-10-05 (YYYY - MM - DD)											
		Rev 1		ABX Lyse																	
PARAMETRES PARAMETERS		UNITES UNITS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE	
				ES60 Care ST		CRP 200		8P		18P		60 8P - 60 18P ADVIA 60									
GB WBC		10 ³ /mm ³ ; 10 ⁹ /l		20.8		20.0		20.2		20.2		20.7								± 1.6	
GR RBC		10 ⁶ /mm ³ ; 10 ¹² /l		5.62		5.57		5.64		5.54		5.60								± 0.20	
HB HGB		g/dl		17.9		17.5		17.8		18.0		17.9								± 0.6	
		g/l		179		175		178		180		179								± 6	
		mmol/l		11.12		10.87		11.05		11.18		11.12								± 0.37	
HT HCT		%		50.0		49.6		50.2		49.3		49.8								± 3.0	
		l/l		0.500		0.496		0.502		0.493		0.498								± 0.030	
VGM MCV		µm ³ ; fl		89		89		89		89		89								± 4	
TGMH MCH		pg		31.9		31.4		31.6		32.5		32.0								± 2.0	
		fmol		1.98		1.95		1.96		2.02		1.98								± 0.12	
CCMH MCHC		g/dl		35.8		35.3		35.5		36.5		35.9								± 3.0	
		g/l		358		353		355		365		359								± 30	
		mmol/l		22.22		21.92		22.02		22.67		22.30								± 1.86	
IDR RDW		%		13.5		13.5		N/A		13.5		12.5								± 3.0	
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		520		515		520		530		522								± 45	
VPM MPV		µm ³ ; fl		8.0		7.8		N/A		8.5		7.4								± 2.0	
LYMPHO		#		3.43		3.40		N/A		2.63		3.42								± 1.00	
		%		16.5		17.0		N/A		13.0		16.5								± 5.0	
MONO		#		1.14		1.20		N/A		0.91		1.04								± 0.80	
		%		5.5		6.0		N/A		4.5		5.0								± 4.0	
GRANULO		#		16.22		15.40		N/A		16.67		16.25								± 1.40	
		%		78.0		77.0		N/A		82.5		78.5								± 7.0	