



LOT		MX 428 L		Rev 1		CONTROL		L		 (Exp.)		2021-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1	2.0			2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.38	2.34			2.38						± 0.15	
HB HGB	g/dl	6.0	5.9			6.1						± 0.4	
	g/l	69	59			61						± 4	
	mmol/l	3.73	3.66			3.79						± 0.25	
HT HCT	%	16.4	15.9			16.4						± 2.0	
	l/l	0.164	0.159			0.164						± 0.020	
VGM MCV	µm ³ ; fl	69	68			69						± 4	
TGMH MCH	pg	25.2	25.2			25.6						± 2.0	
	fmol	1.57	1.57			1.59						± 0.12	
CCMH MCHC	g/dl	36.5	37.1			37.1						± 3.0	
	g/l	365	371			371						± 30	
	mmol/l	22.69	23.03			23.07						± 1.86	
IDR RDW	%	14.5	14.5			13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	77	81			77						± 20	
VPM MPV	µm ² ; fl	9.1	8.6			8.1						± 2.0	
LYMPHO	#	1.40	1.40			1.34						± 0.40	
	%	66.5	70.0			67.0						± 8.0	
MONO	#	0.20	0.20			0.20						± 0.20	
	%	8.5	9.0			8.5						± 6.0	
GRANULO	#	0.53	0.42			0.49						± 0.40	
	%	25.0	21.0			24.5						± 7.0	



LOT		MX 428 N		Rev 1		CONTROL		N		 (Exp.)		2021-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.7	7.5			7.7						± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.65	4.60			4.61						± 0.18	
HB HGB	g/dl	13.4	13.2			13.5						± 0.5	
	g/l	134	132			135						± 5	
	mmol/l	8.32	8.20			8.38						± 0.31	
HT HCT	%	37.2	36.8			37.3						± 2.5	
	l/l	0.372	0.368			0.373						± 0.025	
VGM MCV	µm ³ ; fl	80	80			81						± 4	
TGMH MCH	pg	28.8	28.7			29.3						± 2.0	
	fmol	1.79	1.78			1.82						± 0.12	
CCMH MCHC	g/dl	36.0	35.9			36.2						± 3.0	
	g/l	360	359			362						± 30	
	mmol/l	22.37	22.28			22.45						± 1.86	
IDR RDW	%	13.7	13.8			13.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	252	248			250						± 40	
VPM MPV	µm ³ ; fl	8.2	8.1			7.5						± 2.0	
LYMPHO	#	2.66	2.63			2.62						± 0.40	
	%	34.5	35.0			34.0						± 5.0	
MONO	#	0.54	0.60			0.58						± 0.30	
	%	7.0	8.0			7.5						± 4.0	
GRANULO	#	4.50	4.28			4.50						± 0.60	
	%	58.5	57.0			58.5						± 6.0	

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513



LOT		MX 428 H		Rev 1		CONTROL		H		 (Exp.)		2021-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.8	19.8			20.5						± 1.6	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.68	5.55			5.63						± 0.20	
HB HGB	g/dl	17.8	17.5			17.9						± 0.6	
	g/l	178	175			179						± 6	
	mmol/l	11.05	10.87			11.12						± 0.37	
HT HCT	%	50.0	48.8			50.1						± 3.0	
	l/l	0.500	0.488			0.501						± 0.030	
VGM MCV	µm ³ ; fl	88	88			89						± 4	
TGMH MCH	pg	31.3	31.5			31.8						± 2.0	
	fmol	1.95	1.96			1.97						± 0.12	
CCMH MCHC	g/dl	35.6	35.8			35.7						± 3.0	
	g/l	356	358			357						± 30	
	mmol/l	22.11	22.25			22.18						± 1.86	
IDR RDW	%	13.2	13.2			12.0						± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	498	490			493						± 45	
VPM MPV	µm ³ ; fl	8.1	8.1			7.5						± 2.0	
LYMPHO	#	3.43	3.17			3.28						± 1.00	
	%	16.5	16.0			16.0						± 5.0	
MONO	#	0.94	1.09			0.92						± 0.80	
	%	4.5	5.5			4.5						± 4.0	
GRANULO	#	16.43	15.54			16.30						± 1.40	
	%	79.0	78.5			79.5						± 7.0	