

LOT		MX 429 L		CONTROL		L		 (Exp.)		2021-08-05 (YYYY - MM - DD)			
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP 200			60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0			2.0					± 0.4		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.27	2.24			2.28					± 0.15		
HB HGB	g/dl	6.1	6.0			6.1					± 0.4		
	g/l	61	60			61					± 4		
	mmol/l	3.79	3.73			3.79					± 0.25		
HT HCT	%	16.6	16.4			16.9					± 2.0		
	l/l	0.166	0.164			0.169					± 0.020		
VGM MCV	µm ³ ; fl	73	73			74					± 4		
TGMH MCH	pg	26.9	26.8			26.8					± 2.0		
	fmol	1.67	1.66			1.66					± 0.12		
CCMH MCHC	g/dl	36.8	36.7			36.2					± 3.0		
	g/l	368	367			362					± 30		
	mmol/l	22.86	22.79			22.45					± 1.86		
IDR RDW	%	15.0	14.0			14.0					± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	72	70			73					± 20		
VPM MPV	µm ³ ; fl	8.2	8.1			7.6					± 2.0		
LYMPHO	#	1.34	1.38			1.34					± 0.40		
	%	67.0	69.0			67.0					± 8.0		
MONO	#	0.20	0.20			0.20					± 0.20		
	%	8.0	8.0			7.5					± 6.0		
GRANULO	#	0.50	0.46			0.51					± 0.40		
	%	25.0	23.0			25.5					± 7.0		



LOT		MX 429 N		CONTROL		N		 (Exp.)		2021-08-05 (YYYY - MM - DD)									
		Rev 1		ABX Lyse															
PARAMETRES PARAMETERS		UNITES UNITS												TOLERANCES TOLERANCE					
				MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS				PENTRA		PENTRA	
		ES60 Care ST		CRP 200						60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.8	7.6					7.8									± 0.8	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.55	4.50					4.55									± 0.18	
HB	HGB	g/dl	13.4	13.3					13.5									± 0.5	
		g/l	134	133					135								± 5		
		mmol/l	8.32	8.26					8.38								± 0.31		
HT	HCT	%	37.8	37.4					37.8									± 2.5	
		l/l	0.378	0.374					0.378								± 0.025		
VGM	MCV	µm ³ ; fl	83	83					83									± 4	
TGMH	MCH	pg	29.5	29.6					29.7									± 2.0	
		fmol	1.83	1.84					1.84									± 0.12	
CCMH	MCHC	g/dl	35.5	35.6					35.7									± 3.0	
		g/l	355	356					357								± 30		
		mmol/l	22.03	22.11					22.20								± 1.86		
IDR	RDW	%	14.0	13.5					13.0								± 3.0		
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	251	243					246									± 40	
VPM	MPV	µm ³ ; fl	8.1	7.9					7.5									± 2.0	
LYMPHO	#		2.73	2.70					2.73									± 0.40	
	%		35.0	35.5					35.0									± 5.0	
MONO	#		0.62	0.65					0.62									± 0.30	
	%		8.0	8.5					8.0									± 4.0	
GRANULO	#		4.45	4.26					4.45									± 0.60	
	%		57.0	56.0					57.0									± 6.0	

Ref: TEMP-0830 Rev.39 FRONT / RECTO 130060513

LOT		MX 429 H		CONTROL		H		 (Exp.)		2021-08-05 (YYYY - MM - DD)									
		Rev 1		ABX Lyse															
PARAMETRES PARAMETERS		UNITES UNITS										TOLERANCES TOLERANCE							
				MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA	
		ES60 Care ST		CRP 200				60 BP - 60 18P ADVIA 60											
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	20.6	20.2				20.7										± 1.6	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.64	5.55				5.63										± 0.20	
HB	HGB	g/dl	18.0	17.9				18.2										± 0.6	
		g/l	180	179				182										± 6	
		mmol/l	11.18	11.12				11.30										± 0.37	
HT	HCT	%	50.8	50.5				51.2										± 3.0	
		l/l	0.508	0.505				0.512										± 0.030	
VGM	MCV	µm ³ ; fl	90	91				91										± 4	
TGMH	MCH	pg	31.9	32.3				32.3										± 2.0	
		fmol	1.98	2.00				2.01										± 0.12	
CCMH	MCHC	g/dl	35.5	35.4				35.5										± 3.0	
		g/l	355	354				355										± 30	
		mmol/l	22.02	22.01				22.06										± 1.86	
IDR	RDW	%	13.7	13.5				12.5										± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	499	480				494										± 45	
VPM	MPV	µm ³ ; fl	7.8	7.8				7.4										± 2.0	
LYMPHO	#	3.19	3.33				3.21											± 1.00	
	%	15.5	16.5				15.5											± 5.0	
MONO	#	1.03	1.11				0.93											± 0.80	
	%	5.0	5.5				4.5											± 4.0	
GRANULO	#	16.38	15.76				16.56											± 1.40	
	%	79.5	78.0				80.0											± 7.0	