



LOT		MX 410 L		Rev 1		CONTROL		L		 (Exp.)		2018-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.0	2.0	2.0						± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.28	2.25	2.30	2.27	2.28						± 0.15	
HB HGB	g/dl	6.1	6.0	6.0	5.9	6.1						± 0.4	
	g/l	61	60	60	59	61						± 4	
	mmol/l	3.79	3.73	3.73	3.66	3.79						± 0.25	
HT HCT	%	16.2	16.0	16.3	16.1	16.2						± 2.0	
	l/l	0.162	0.160	0.163	0.161	0.162						± 0.020	
VGM MCV	µm ³ ; fl	71	71	71	71	71						± 4	
TGMH MCH	pg	26.8	26.7	26.1	26.0	26.8						± 2.0	
	fmol	1.66	1.66	1.62	1.61	1.66						± 0.12	
CCMH MCHC	g/dl	37.7	37.6	36.7	36.6	37.7						± 3.0	
	g/l	377	376	367	366	377						± 30	
	mmol/l	23.40	23.32	22.82	22.73	23.40						± 1.86	
IDR RDW	%	14.5	13.5	N/A	14.0	13.5						± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	83	78	73	80	80						± 20	
VPM MPV	µm ³ ; fl	8.5	8.5	N/A	9.0	8.3						± 2.0	
LYMPHO	#	1.31	1.32	N/A	1.08	1.28						± 0.40	
	%	65.5	66.0	N/A	54.0	64.0						± 8.0	
MONO	#	0.20	0.20	N/A	0.28	0.20						± 0.20	
	%	8.5	9.0	N/A	14.0	10.0						± 6.0	
GRANULO	#	0.52	0.50	N/A	0.64	0.52						± 0.40	
	%	26.0	25.0	N/A	32.0	26.0						± 7.0	

LOT		MX 410 N		CONTROL N				 (Exp.)		2018-06-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7	7.7	7.5	7.6	7.8					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.58	4.55	4.62	4.56	4.61					± 0.18
HB HGB	g/dl	13.4	13.4	13.5	13.4	13.5					± 0.5
	g/l	134	134	135	134	135					± 5
	mmol/l	8.32	8.32	8.38	8.32	8.38					± 0.31
HT HCT	%	37.1	36.9	37.4	36.9	37.3					± 2.5
	l/l	0.371	0.369	0.374	0.369	0.373					± 0.025
VGM MCV	µm ³ ; fl	81	81	81	81	81					± 4
TGMH MCH	pg	29.3	29.5	29.2	29.4	29.3					± 2.0
	fmol	1.82	1.83	1.81	1.82	1.82					± 0.12
CCMH MCHC	g/dl	36.1	36.4	36.1	36.3	36.2					± 3.0
	g/l	361	364	361	363	362					± 30
	mmol/l	22.43	22.58	22.40	22.53	22.45					± 1.86
IDR RDW	%	14.5	13.5	N/A	14.5	13.0					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	263	256	252	260	258					± 40
VPM MPV	µm ³ ; fl	8.2	8.2	N/A	8.7	7.8					± 2.0
LYMPHO	#	2.70	2.66	N/A	2.20	2.65					± 0.40
	%	35.0	34.5	N/A	29.0	34.0					± 5.0
MONO	#	0.58	0.62	N/A	0.65	0.62					± 0.30
	%	7.5	8.0	N/A	8.5	8.0					± 4.0
GRANULO	#	4.43	4.43	N/A	4.75	4.52					± 0.60
	%	57.5	57.5	N/A	62.5	58.0					± 6.0

Ref: TEMP-0830 Rev.38 FRONT / RECTO 1219930084-B



LOT		MX 410 H		CONTROL H		 (Exp.)		2018-06-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60											
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.4		20.3		20.1		20.3		20.5								± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.69		5.63		5.75		5.65		5.70								± 0.20			
HB	HGB	g/dl		18.0		17.9		18.0		18.0		18.0								± 0.6			
		g/l		180		179		180		180		180								± 6			
		mmol/l		11.18		11.12		11.18		11.18		11.18								± 0.37			
HT	HCT	%		50.1		49.5		50.6		49.7		50.2								± 3.0			
		l/l		0.501		0.495		0.506		0.497		0.502								± 0.030			
VGM	MCV	µm ³ ; fl		88		88		88		88		88								± 4			
TGMH	MCH	pg		31.6		31.8		31.3		31.9		31.6								± 2.0			
		fmol		1.96		1.97		1.94		1.98		1.96								± 0.12			
CCMH	MCHC	g/dl		35.9		36.1		35.6		36.2		35.9								± 3.0			
		g/l		359		361		356		362		359								± 30			
		mmol/l		22.32		22.44		22.09		22.48		22.28								± 1.86			
IDR	RDW	%		14.0		13.0		N/A		13.5		12.5								± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		505		495		505		500		498								± 45			
VPM	MPV	µm ³ ; fl		8.2		8.3		N/A		9.0		8.0								± 2.0			
LYMPHO	#			3.26		3.15		N/A		2.64		3.18								± 1.00			
	%			16.0		15.5		N/A		13.0		15.5								± 5.0			
MONO	#			0.92		0.91		N/A		0.81		0.92								± 0.80			
	%			4.5		4.5		N/A		4.0		4.5								± 4.0			
GRANULO	#			16.22		16.24		N/A		16.85		16.40								± 1.40			
	%			79.5		80.0		N/A		83.0		80.0								± 7.0			