



LOT

MX 113 L

Rev 1

CONTROL

L



(Exp.)

2014-02-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA		
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60		60 60 C+ ES60	80 XL80	120 - 120 RETIC DF120 - DX120		
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.0	2.0	2.1		2.0	2.1		± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.40	2.39	2.42	2.41	2.38		2.42	2.36		± 0.15	
HB HGB	g/dl	6.1	6.1	5.9	6.0	6.0		5.9	5.9		± 0.4	
	g/l	61	61	59	60	60		59	59		± 4	
	mmol/l	3.79	3.79	3.66	3.73	3.73		3.66	3.66		± 0.25	
HT HCT	%	16.3	16.3	16.7	16.9	16.4		16.5	16.8		± 2.0	
	l/l	0.163	0.163	0.167	0.169	0.164		0.165	0.168		± 0.020	
VGM MCV	µm ³ ; fl	68	68	69	70	69		68	71		± 4	
TGMH MCH	pg	25.4	25.5	24.4	24.9	25.2		24.4	25.0		± 2.0	
	fmol	1.58	1.58	1.51	1.55	1.57		1.51	1.55		± 0.12	
CCMH MCHC	g/dl	37.4	37.5	35.3	35.6	36.5		35.9	35.2		± 3.0	
	g/l	374	375	353	356	365		359	352		± 30	
	mmol/l	23.21	23.31	21.94	22.09	22.69		22.26	21.87		± 1.86	
IDR RDW	%	14.2	13.5	N/A	13.5	13.8		15.0	15.9		± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	80	78	73	75	80		* 77	* 75		± 20	*± 22
VPM MPV	µm ³ ; fl	8.4	8.5	N/A	8.8	7.8		9.9	10.0		± 2.0	
LYMPHO	#	1.26	1.28	N/A	1.00	1.34		N/A	N/A		± 0.40	
	%	63.0	64.0	N/A	50.0	64.0		N/A	N/A		± 7.0	
MONO	#	0.20	0.20	N/A	0.34	0.21		N/A	N/A		± 0.20	
	%	10.0	10.0	N/A	17.0	10.0		N/A	N/A		± 6.0	
GRANULO	#	0.54	0.52	N/A	0.66	0.55		N/A	N/A		± 0.40	
	%	27.0	26.0	N/A	33.0	26.0		N/A	N/A		± 7.0	

* For ABX Pentra 60 Pentra 80 instruments only



LOT		MX 113 N		Rev 1		CONTROL N		 (Exp.)		2014-02-05 (YYYY - MM - DD)		
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA		
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60		60 60 C+ ES60	80 XL80	120 - 120 RETIC DF120 - DX120		
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.7	7.6	7.5	7.5	7.7		7.7	7.8		± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.58	4.56	4.60	4.53	4.58		4.58	4.56		± 0.18	
HB HGB	g/dl	13.4	13.3	13.2	13.4	13.4		13.4	13.3		± 0.5	
	g/l	134	133	132	134	134		134	133		± 5	
	mmol/l	8.32	8.26	8.20	8.32	8.32		8.32	8.26		± 0.31	
HT HCT	%	37.1	36.9	37.3	37.1	37.1		36.6	36.5		± 2.5	
	l/l	0.371	0.369	0.373	0.371	0.371		0.366	0.365		± 0.025	
VGM MCV	µm ³ ; fl	81	81	81	82	81		80	80		± 4	
TGMH MCH	pg	29.3	29.2	28.7	29.6	29.3		29.3	29.2		± 2.0	
	fmol	1.82	1.81	1.78	1.84	1.82		1.82	1.81		± 0.12	
CCMH MCHC	g/dl	36.1	36.0	35.4	36.1	36.1		36.6	36.5		± 3.0	
	g/l	361	360	354	361	361		366	365		± 30	
	mmol/l	22.43	22.36	22.00	22.40	22.43		22.71	22.64		± 1.86	
IDR RDW	%	13.8	13.0	N/A	13.0	13.5		14.0	15.0		± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260	252	252	257	255		* 258	* 260		± 40 *± 40	
VPM MPV	µm ³ ; fl	7.7	7.5	N/A	7.9	7.2		8.7	8.8		± 2.0	
LYMPHO	#	2.54	2.43	N/A	1.95	2.50		N/A	N/A		± 0.40	
	%	33.0	32.0	N/A	26.0	32.5		N/A	N/A		± 5.0	
MONO	#	0.58	0.61	N/A	0.64	0.58		N/A	N/A		± 0.30	
	%	7.5	8.0	N/A	8.5	7.5		N/A	N/A		± 4.0	
GRANULO	#	4.58	4.56	N/A	4.91	4.62		N/A	N/A		± 0.60	
	%	59.5	60.0	N/A	65.5	60.0		N/A	N/A		± 6.0	

* For ABX Pentra 60 Pentra 80 instruments only

9930079-B
FRONT / RECTO
Ref TEMP-0830 Rev.32



LOT		MX 113 H		CONTROL		H		 (Exp.)		2014-02-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60				60 60 C+ ES60		80 XL80		120 - 120 RETIC DF120 - DX120		
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.2		20.0		19.7		19.5		20.3				20.8		20.8				± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.68		5.62		5.65		5.60		5.65				5.63		5.64				± 0.20
HB HGB	g/dl	17.9		17.7		17.5		17.9		17.9				17.9		17.7				± 0.6
	g/l	179		177		175		179		179				179		177				± 6
	mmol/l	11.12		10.99		10.87		11.12		11.12				11.12		10.99				± 0.37
HT HCT	%	50.6		50.0		50.3		49.8		50.3				49.5		49.1				± 3.0
	l/l	0.506		0.500		0.503		0.498		0.503				0.495		0.491				± 0.030
VGM MCV	µm ³ ; fl	89		89		89		89		89				88		87				± 4
TGMH MCH	pg	31.5		31.5		31.0		32.0		31.7				31.8		31.4				± 2.0
	fmol	1.96		1.96		1.92		1.98		1.97				1.97		1.95				± 0.12
CCMH MCHC	g/dl	35.4		35.4		34.8		35.9		35.6				36.1		36.1				± 3.0
	g/l	354		354		348		359		356				361		361				± 30
	mmol/l	21.99		21.98		21.61		22.30		22.11				22.44		22.40				± 1.86
IDR RDW	%	13.3		12.5		N/A		12.6		13.2				13.0		13.9				± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	492		477		500		505		485				* 495		* 515				± 45 *± 55
VPM MPV	µm ³ ; fl	7.7		7.5		N/A		7.8		7.2				8.9		8.8				± 2.0
LYMPHO	#	3.13		3.00		N/A		2.33		3.05				N/A		N/A				± 1.00
	%	15.5		15.0		N/A		12.0		15.0				N/A		N/A				± 5.0
MONO	#	1.11		1.10		N/A		0.80		1.02				N/A		N/A				± 0.80
	%	5.5		5.5		N/A		4.0		5.0				N/A		N/A				± 4.0
GRANULO	#	15.96		15.90		N/A		16.37		16.24				N/A		N/A				± 1.40
	%	79.0		79.5		N/A		84.0		80.0				N/A		N/A				± 7.0

* For ABX Pentra 60 Pentra 80 instruments only