



ABX Minotrol 16

LOT		MX 409 L		CONTROL		L		 (Exp.)		2018-04-05 (YYYY-MM-DD)													
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		2.1		2.1		2.1		2.1		2.1								± 0.4			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.27		2.24		2.29		2.26		2.26								± 0.15			
HB	HGB	g/dl		6.0		5.9		5.9		5.9		6.0								± 0.4			
		g/l		60		59		59		59		60								± 4			
		mmol/l		3.73		3.66		3.66		3.66		3.73								± 0.25			
HT	HCT	%		15.9		15.7		16.3		16.0		16.0								± 2.0			
		l/l		0.159		0.157		0.163		0.160		0.160								± 0.020			
VGM	MCV	µm ³ ; fl		70		70		71		71		71								± 4			
TGMH	MCH	pg		26.4		26.3		25.8		26.1		26.5								± 2.0			
		fmol		1.64		1.64		1.60		1.62		1.65								± 0.12			
CCMH	MCHC	g/dl		37.8		37.6		36.3		36.8		37.4								± 3.0			
		g/l		378		376		363		368		374								± 30			
		mmol/l		23.45		23.37		22.53		22.83		23.22								± 1.86			
IDR	RDW	%		14.5		13.5		N/A		14.0		13.5								± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		79		76		74		77		83								± 20			
VPM	MPV	µm ³ ; fl		8.9		8.6		N/A		8.5		7.9								± 2.0			
LYMPHO	#			1.33		1.34		N/A		1.11		1.31								± 0.40			
		%		63.5		64.0		N/A		53.0		62.5								± 8.0			
MONO	#			0.20		0.20		N/A		0.29		0.21								± 0.20			
		%		8.5		9.5		N/A		14.0		10.0								± 6.0			
GRANULO	#			0.59		0.56		N/A		0.69		0.58								± 0.40			
		%		28.0		26.5		N/A		33.0		27.5								± 7.0			

LOT		MX 409	N	CONTROL		N	 (Exp.)				2018-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS		UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
			MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
			ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.4	7.6	7.5	7.6	7.7				± 0.8	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.68	4.69	4.70	4.68	4.72				± 0.18	
HB	HGB	g/dl	13.6	13.5	13.4	13.7	13.6				± 0.5	
		g/l	136	135	134	137	136				± 5	
		mmol/l	8.45	8.38	8.32	8.51	8.45				± 0.31	
HT	HCT	%	37.4	37.5	37.6	37.4	37.8				± 2.5	
		l/l	0.374	0.375	0.376	0.374	0.378				± 0.025	
VGM	MCV	µm ³ ; fl	80	80	80	80	80				± 4	
TGMH	MCH	pg	29.1	28.8	28.5	29.3	28.8				± 2.0	
		fmol	1.80	1.79	1.77	1.82	1.79				± 0.12	
CCMH	MCHC	g/dl	36.3	36.0	35.6	36.6	36.0				± 3.0	
		g/l	363	360	356	366	360				± 30	
		mmol/l	22.56	22.34	22.13	22.72	22.37				± 1.86	
IDR	RDW	%	15.0	13.5	N/A	14.5	13.5				± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	250	255	252	252	255				± 40	
VPM	MPV	µm ³ ; fl	9.2	9.0	N/A	9.0	8.4				± 2.0	
LYMPHO	#	2.66	2.66	N/A	2.24	2.62					± 0.40	
	%	36.0	35.0	N/A	29.5	34.0					± 5.0	
MONO	#	0.48	0.53	N/A	0.61	0.54					± 0.30	
	%	6.5	7.0	N/A	8.0	7.0					± 4.0	
GRANULO	#	4.26	4.41	N/A	4.75	4.54					± 0.60	
	%	57.5	58.0	N/A	62.5	59.0					± 6.0	



ABX Minotrol 16



LOT		MX 409 H		CONTROL		H		 (Exp.)		2018-04-05 (YYYY-MM-DD)													
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		19.8		20.2		20.0		20.2		20.4								± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.65		5.69		5.74		5.70		5.71								± 0.20			
HB	HGB	g/dl		17.6		17.6		17.7		17.9		17.8								± 0.6			
		g/l		176		176		177		179		178								± 6			
		mmol/l		10.93		10.93		10.99		11.12		11.05								± 0.37			
HT	HCT	%		49.7		49.5		50.5		50.2		50.2								± 3.0			
		l/l		0.497		0.495		0.505		0.502		0.502								± 0.030			
VGM	MCV	µm ³ ; fl		88		87		88		88		88								± 4			
TGMH	MCH	pg		31.2		30.9		30.8		31.4		31.2								± 2.0			
		fmol		1.93		1.92		1.91		1.95		1.94								± 0.12			
CCMH	MCHC	g/dl		35.4		35.6		35.0		35.7		35.4								± 3.0			
		g/l		354		356		350		357		354								± 30			
		mmol/l		21.98		22.08		21.76		22.16		22.00								± 1.86			
IDR	RDW	%		13.5		13.0		N/A		13.5		13.0								± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		500		505		525		510		508								± 45			
VPM	MPV	µm ³ ; fl		8.3		8.0		N/A		8.2		7.5								± 2.0			
LYMPHO		#		3.27		3.13		N/A		2.63		3.16								± 1.00			
		%		16.5		15.5		N/A		13.0		15.5								± 5.0			
MONO		#		0.99		1.01		N/A		0.81		1.02								± 0.80			
		%		5.0		5.0		N/A		4.0		5.0								± 4.0			
GRANULO		#		15.54		16.06		N/A		16.77		16.22								± 1.40			
		%		78.5		79.5		N/A		83.0		79.5								± 7.0			