



LOT

MX 074 L

Rev 1

CONTROL

L



2014-10-05

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.1	2.0	2.0	2.0	2.0					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.37	2.36	2.40	2.37	2.35					± 0.15
HB HGB	g/dl	6.0	5.9	5.9	5.9	5.9					± 0.4
	g/l	60	59	59	59	59					± 4
	mmol/l	3.73	3.66	3.66	3.66	3.66					± 0.25
HT HCT	%	16.1	16.0	16.6	16.4	16.2					± 2.0
	l/l	0.161	0.160	0.166	0.164	0.162					± 0.020
VGM MCV	µm ³ ; fl	68	68	69	69	69					± 5
TGMH MCH	pg	25.3	25.0	24.6	24.9	25.1					± 2.0
	fmol	1.57	1.55	1.53	1.55	1.56					± 0.12
CCMH MCHC	g/dl	37.2	36.8	35.6	36.1	36.4					± 3.0
	g/l	372	368	356	361	364					± 30
	mmol/l	23.12	22.83	22.13	22.41	22.60					± 1.86
IDR RDW	%	14.0	13.8	N/A	13.0	13.6					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	84	80	77	77	81					± 20
VPM MPV	µm ³ ; fl	8.9	8.7	N/A	9.0	8.2					± 2.0
LYMPHO	#	1.30	1.25	N/A	1.00	1.23					± 0.40
	%	62.0	62.5	N/A	50.0	62.0					± 7.0
MONO	#	0.21	0.20	N/A	0.30	0.20					± 0.20
	%	10.0	10.0	N/A	15.0	9.5					± 6.0
GRANULO	#	0.59	0.55	N/A	0.70	0.57					± 0.40
	%	28.0	27.5	N/A	35.0	28.5					± 7.0



LOT		MX 074 N		CONTROL N				 (Exp.)		2014-10-05 (YYYY-MM-DD)											
		Rev 1		ABX Lyse																	
PARAMETRES PARAMETERS		UNITES UNITS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE	
				ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		7.9		7.7		7.6		7.5		7.8								± 0.8	
GR RBC		10 ⁶ /mm ³ ; 10 ¹² /l		4.63		4.62		4.68		4.60		4.66								± 0.18	
HB HGB		g/dl		13.5		13.4		13.4		13.4		13.4								± 0.5	
		g/l		135		134		134		134		134								± 5	
		mmol/l		8.38		8.32		8.32		8.32		8.32								± 0.31	
HT HCT		%		36.6		37.0		37.4		37.3		37.3								± 2.5	
		l/l		0.366		0.370		0.374		0.373		0.373								± 0.025	
VGM MCV		µm ³ ; fl		79		80		80		81		80								± 5	
TGMH MCH		pg		29.2		29.0		28.6		29.1		28.8								± 2.0	
		fmol		1.81		1.80		1.78		1.81		1.79								± 0.12	
CCMH MCHC		g/dl		36.9		36.3		35.8		36.0		35.9								± 3.0	
		g/l		369		363		358		360		359								± 30	
		mmol/l		22.92		22.51		22.23		22.33		22.32								± 1.86	
IDR RDW		%		13.4		13.0		N/A		12.5		13.4								± 3.0	
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		265		262		268		260		265								± 40	
VPM MPV		µm ³ ; fl		8.2		8.0		N/A		8.5		7.6								± 2.0	
LYMPHO		#		2.69		2.58		N/A		2.10		2.65								± 0.40	
		%		34.0		33.5		N/A		28.0		34.0								± 5.0	
MONO		#		0.59		0.62		N/A		0.60		0.59								± 0.30	
		%		7.5		8.0		N/A		8.0		7.5								± 4.0	
GRANULO		#		4.62		4.50		N/A		4.80		4.56								± 0.60	
		%		58.5		58.5		N/A		64.0		58.5								± 6.0	

9930084-B
FRONT / RECTO
Ref: TEMP-0830 Rev.35



LOT		MX 074 H		CONTROL H		 (Exp.)		2014-10-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse										TOLERANCES TOLERANCE									
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60											
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.4	20.0	20.0	19.8	20.3														± 1.6			
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.55	5.50	5.60	5.56	5.58														± 0.20			
HB HGB	g/dl	17.8	17.6	17.5	17.9	17.7														± 0.6			
	g/l	178	176	175	179	177														± 6			
	mmol/l	11.05	10.93	10.87	11.12	10.99														± 0.37			
HT HCT	%	49.4	49.0	49.8	49.5	49.7														± 3.0			
	l/l	0.494	0.490	0.498	0.495	0.497														± 0.030			
VGM MCV	µm ³ ; fl	89	89	89	89	89														± 5			
TGMH MCH	pg	32.1	32.0	31.3	32.2	31.7														± 2.0			
	fmol	1.99	1.99	1.94	2.00	1.97														± 0.12			
CCMH MCHC	g/dl	36.0	36.0	35.1	36.2	35.6														± 3.0			
	g/l	360	360	351	362	356														± 30			
	mmol/l	22.38	22.33	21.80	22.46	22.13														± 1.86			
IDR RDW	%	12.9	12.5	N/A	12.3	13.0														± 3.0			
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	515	507	535	520	515														± 45			
VPM MPV	µm ³ ; fl	8.0	8.0	N/A	8.4	7.5														± 2.0			
LYMPHO	#	3.16	3.00	N/A	2.48	3.05														± 1.00			
	%	15.5	15.0	N/A	12.5	15.0														± 5.0			
MONO	#	1.02	1.10	N/A	0.80	1.02														± 0.80			
	%	5.0	5.5	N/A	4.0	5.0														± 4.0			
GRANULO	#	16.22	15.90	N/A	16.52	16.23														± 1.40			
	%	79.5	79.5	N/A	83.5	80.0														± 7.0			