



LOT		MX 094 L		CONTROL		L		 (Exp.)		2014-12-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60														
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.0	1.9	2.0													± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.33	2.29	2.32	2.34	2.32													± 0.15	
HB HGB	g/dl	6.1	6.0	6.0	6.0	6.0													± 0.4	
	g/l	61	60	60	60	60													± 4	
	mmol/l	3.79	3.73	3.73	3.73	3.73													± 0.25	
HT HCT	%	16.1	15.8	16.2	16.4	16.2													± 2.0	
	l/l	0.161	0.158	0.162	0.164	0.162													± 0.020	
VGM MCV	µm ³ ; fl	69	69	70	70	70													± 4	
TGMH MCH	pg	26.2	26.2	25.9	25.6	25.9													± 2.0	
	fmol	1.63	1.63	1.61	1.59	1.61													± 0.12	
CCMH MCHC	g/dl	37.9	38.0	36.9	36.6	36.9													± 3.0	
	g/l	379	380	369	366	369													± 30	
	mmol/l	23.56	23.58	22.94	22.75	22.94													± 1.86	
IDR RDW	%	14.0	13.5	N/A	13.5	14.0													± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	80	75	74	78	79													± 20	
VPM MPV	µm ³ ; fl	9.0	8.8	N/A	9.0	8.0													± 2.0	
LYMPHO	#	1.27	1.28	N/A	0.94	1.26													± 0.40	
	%	63.5	64.0	N/A	49.5	63.0													± 7.0	
MONO	#	0.21	0.21	N/A	0.33	0.21													± 0.20	
	%	10.5	10.5	N/A	17.5	10.5													± 6.0	
GRANULO	#	0.52	0.51	N/A	0.63	0.53													± 0.40	
	%	26.0	25.5	N/A	33.0	26.5													± 7.0	



LOT		MX 094 N		CONTROL N		 (Exp.)		2014-12-05 (YYYY-MM-DD)				
		Rev 1										
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse										TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA		
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60						
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.8	7.7	7.7	7.5	7.7					± 0.8	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.63	4.55	4.64	4.60	4.60					± 0.18	
HB HGB	g/dl	13.4	13.3	13.4	13.4	13.4					± 0.5	
	g/l	134	133	134	134	134					± 5	
	mmol/l	8.32	8.26	8.32	8.32	8.32					± 0.31	
HT HCT	%	37.0	36.4	37.1	36.8	36.8					± 2.5	
	l/l	0.370	0.364	0.371	0.368	0.368					± 0.025	
VGM MCV	µm ³ ; fl	80	80	80	80	80					± 4	
TGMH MCH	pg	28.9	29.2	28.9	29.1	29.1					± 2.0	
	fmol	1.80	1.82	1.79	1.81	1.81					± 0.12	
CCMH MCHC	g/dl	36.2	36.5	36.1	36.4	36.4					± 3.0	
	g/l	362	365	361	364	364					± 30	
	mmol/l	22.47	22.69	22.42	22.61	22.61					± 1.86	
IDR RDW	%	13.8	13.0	N/A	13.0	13.7					± 3.0	
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	260	254	260	260	255					± 40	
VPM MPV	µm ³ ; fl	8.5	8.5	N/A	9.0	8.0					± 2.0	
LYMPHO	#	2.65	2.62	N/A	2.03	2.62					± 0.40	
	%	34.0	34.0	N/A	27.0	34.0					± 5.0	
MONO	#	0.59	0.58	N/A	0.68	0.54					± 0.30	
	%	7.5	7.5	N/A	9.0	7.0					± 4.0	
GRANULO	#	4.56	4.50	N/A	4.80	4.54					± 0.60	
	%	58.5	58.5	N/A	64.0	59.0					± 6.0	

9930084-B
FRONT / RECTO
Ref: TEMP-0830 Rev.35



LOT		MX 094 H		CONTROL H		 (Exp.)		2014-12-05 (YYYY-MM-DD)													
		Rev 1		ABX Lyse																	
PARAMETRES PARAMETERS		UNITES UNITS												TOLERANCES TOLERANCE							
				MICROS		MICROS		MICROS		MICROS		MICROS				ARGOS - HELIOS		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60									
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		21.4		21.0		20.8		20.8		21.3								± 1.6	
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.74		5.65		5.72		5.75		5.75								± 0.20	
HB	HGB	g/dl		18.1		17.8		17.9		18.3		18.0								± 0.6	
		g/l		181		178		179		183		180								± 6	
		mmol/l		11.24		11.05		11.12		11.36		11.18								± 0.37	
HT	HCT	%		49.9		49.2		49.8		50.0		50.0								± 3.0	
		l/l		0.499		0.492		0.498		0.500		0.500								± 0.030	
VGM	MCV	µm ³ ; fl		87		87		87		87		87								± 4	
TGMH	MCH	pg		31.5		31.5		31.3		31.8		31.3								± 2.0	
		fmol		1.96		1.96		1.94		1.98		1.94								± 0.12	
CCMH	MCHC	g/dl		36.2		36.2		36.0		36.6		36.0								± 3.0	
		g/l		362		362		360		366		360								± 30	
		mmol/l		22.51		22.49		22.34		22.72		22.34								± 1.86	
IDR	RDW	%		13.0		12.5		N/A		12.5		13.0								± 3.0	
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		520		510		540		535		520								± 45	
VPM	MPV	µm ³ ; fl		8.2		8.0		N/A		8.6		7.7								± 2.0	
LYMPHO		#		3.42		3.26		N/A		2.70		3.41								± 1.00	
		%		16.0		15.5		N/A		13.0		16.0								± 5.0	
MONO		#		1.07		1.26		N/A		0.83		1.07								± 0.80	
		%		5.0		6.0		N/A		4.0		5.0								± 4.0	
GRANULO		#		16.91		16.49		N/A		17.26		16.83								± 1.40	
		%		79.0		78.5		N/A		83.0		79.0								± 7.0	