



LOT

MX 015 L

Rev 1

CONTROL

L



(Exp.)

2015-04-05

(YYYY - MM - DD)

		Rev 1	ABX Lyse								(YYYY-MM-DD)
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	2.1	2.1	2.1	2.1	2.1					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.34	2.35	2.39	2.36	2.34					± 0.15
HB HGB	g/dl	6.1	6.1	6.1	6.0	6.1					± 0.4
	g/l	61	61	61	60	61					± 4
	mmol/l	3.79	3.79	3.79	3.73	3.79					± 0.25
HT HCT	%	16.4	16.2	16.7	16.5	16.4					± 2.0
	l/l	0.164	0.162	0.167	0.165	0.164					± 0.020
VGM MCV	µm ³ ; fl	70	69	70	70	70					± 4
TGMH MCH	pg	26.1	26.0	25.5	25.4	26.1					± 2.0
	fmol	1.62	1.61	1.58	1.58	1.62					± 0.12
CCMH MCHC	g/dl	37.2	37.6	36.5	36.3	37.2					± 3.0
	g/l	372	376	365	363	372					± 30
	mmol/l	23.13	23.36	22.64	22.55	23.13					± 1.86
IDR RDW	%	13.9	14.0	N/A	14.0	13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	78	76	70	75	77					± 20
VPM MPV	µm ³ ; fl	8.6	8.7	N/A	8.6	8.0					± 2.0
LYMPHO	#	1.39	1.37	N/A	1.11	1.35					± 0.40
	%	66.0	65.0	N/A	53.0	64.5					± 8.0
MONO	#	0.20	0.20	N/A	0.29	0.20					± 0.20
	%	9.0	9.0	N/A	14.0	9.5					± 6.0
GRANULO	#	0.53	0.55	N/A	0.69	0.55					± 0.40
	%	25.0	26.0	N/A	33.0	26.0					± 7.0



LOT		MX 015 N		CONTROL N		 (Exp.)		2015-04-05 (YYYY - MM - DD)					
		Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse										TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.6	7.6	7.4	7.5	7.7					± 0.8		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.60	4.58	4.65	4.57	4.60					± 0.18		
HB HGB	g/dl	13.7	13.5	13.5	13.5	13.6					± 0.5		
	g/l	137	135	135	135	136					± 5		
	mmol/l	8.51	8.38	8.38	8.38	8.45					± 0.31		
HT HCT	%	37.3	37.1	37.2	37.0	37.3					± 2.5		
	l/l	0.373	0.371	0.372	0.370	0.373					± 0.025		
VGM MCV	µm ³ ; fl	81	81	80	81	81					± 4		
TGMH MCH	pg	29.8	29.5	29.0	29.5	29.6					± 2.0		
	fmol	1.85	1.83	1.80	1.83	1.84					± 0.12		
CCMH MCHC	g/dl	36.8	36.4	36.3	36.5	36.5					± 3.0		
	g/l	368	364	363	365	365					± 30		
	mmol/l	22.83	22.60	22.54	22.65	22.67					± 1.86		
IDR RDW	%	13.5	13.0	N/A	13.0	13.3					± 3.0		
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	255	255	250	250	260					± 40		
VPM MPV	µm ³ ; fl	8.3	8.3	N/A	8.6	7.6					± 2.0		
LYMPHO	#	2.58	2.47	N/A	2.06	2.54					± 0.40		
	%	34.0	32.5	N/A	27.5	33.0					± 5.0		
MONO	#	0.53	0.57	N/A	0.56	0.54					± 0.30		
	%	7.0	7.5	N/A	7.5	7.0					± 4.0		
GRANULO	#	4.48	4.56	N/A	4.88	4.62					± 0.60		
	%	59.0	60.0	N/A	65.0	60.0					± 6.0		

9930084-B
FRONT / RECTO
Ref: TEMP-0830 Rev.36



LOT		MX 015 H		CONTROL H		 (Exp.)		2015-04-05 (YYYY-MM-DD)													
		Rev 1		ABX Lyse																	
PARAMETRES PARAMETERS		UNITES UNITS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE	
				ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		20.4		20.5		19.8		20.1		20.6								± 1.6	
GR RBC		10 ⁶ /mm ³ ; 10 ¹² /l		5.67		5.63		5.70		5.65		5.70								± 0.20	
HB HGB		g/dl		18.0		17.7		17.7		18.0		17.9								± 0.6	
		g/l		180		177		177		180		179								± 6	
		mmol/l		11.18		10.99		10.99		11.18		11.12								± 0.37	
HT HCT		%		49.9		49.5		50.2		49.7		50.7								± 3.0	
		l/l		0.499		0.495		0.502		0.497		0.507								± 0.030	
VGM MCV		µm ³ ; fl		88		88		88		88		89								± 4	
TGMH MCH		pg		31.7		31.4		31.1		31.9		31.4								± 2.0	
		fmol		1.97		1.95		1.93		1.98		1.95								± 0.12	
CCMH MCHC		g/dl		36.1		35.7		35.3		36.2		35.3								± 3.0	
		g/l		361		357		353		362		353								± 30	
		mmol/l		22.40		22.19		21.91		22.48		21.91								± 1.86	
IDR RDW		%		13.2		12.9		N/A		13.0		13.0								± 3.0	
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		505		505		515		510		520								± 45	
VPM MPV		µm ³ ; fl		7.8		7.8		N/A		8.2		7.3								± 2.0	
LYMPHO		#		3.06		2.87		N/A		2.41		2.88								± 1.00	
		%		15.0		14.0		N/A		12.0		14.0								± 5.0	
MONO		#		1.22		1.23		N/A		0.80		1.13								± 0.80	
		%		6.0		6.0		N/A		4.0		5.5								± 4.0	
GRANULO		#		16.12		16.40		N/A		16.88		16.58								± 1.40	
		%		79.0		80.0		N/A		84.0		80.5								± 7.0	