



LOT		MX 095 L		CONTROL		L		 (Exp.)		2015-12-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS	UNITES UNITS										TOLERANCES TOLERANCE		
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA			
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.1	2.0	2.0	2.0					± 0.4		
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.45	2.45	2.45	2.45	2.45					± 0.15		
HB HGB	g/dl	6.2	6.2	6.2	6.1	6.2					± 0.4		
	g/l	62	62	62	61	62					± 4		
	mmol/l	3.85	3.85	3.85	3.79	3.85					± 0.25		
HT HCT	%	16.7	16.7	16.9	16.9	16.7					± 2.0		
	l/l	0.167	0.167	0.169	0.169	0.167					± 0.020		
VGM MCV	µm ³ ; fl	68	68	69	69	68					± 4		
TGMH MCH	pg	25.3	25.3	25.3	24.9	25.3					± 2.0		
	fmol	1.57	1.57	1.57	1.55	1.57					± 0.12		
CCMH MCHC	g/dl	37.2	37.2	36.7	36.1	37.2					± 3.0		
	g/l	372	372	367	361	372					± 30		
	mmol/l	23.11	23.11	22.78	22.41	23.11					± 1.86		
IDR RDW	%	14.5	14.0	N/A	14.5	14.0					± 3.0		
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	78	71	67	72	75					± 20		
VPM MPV	µm ³ ; fl	9.1	8.8	N/A	8.9	8.2					± 2.0		
LYMPHO	#	1.22	1.30	N/A	0.97	1.22					± 0.40		
	%	61.0	62.0	N/A	48.5	61.0					± 8.0		
MONO	#	0.21	0.23	N/A	0.32	0.22					± 0.20		
	%	10.5	11.0	N/A	16.0	11.0					± 6.0		
GRANULO	#	0.57	0.57	N/A	0.71	0.56					± 0.40		
	%	28.5	27.0	N/A	35.5	28.0					± 7.0		

LOT		MX 095 N		CONTROL N				 (Exp.)		2015-12-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		7.7		7.7		4.50		4.50		7.6		7.6		7.7						± 0.8	
GR RBC		10 ⁹ /mm ³ ; 10 ¹² /l		4.50		4.50		4.50		4.50		7.6		7.6		7.7						± 0.18	
HB HGB		g/dl		13.4		13.3		13.3		13.3		13.4		13.4		13.4						± 0.5	
		g/l		134		133		133		133		134		134		134						± 5	
		mmol/l		8.32		8.26		8.26		8.26		8.32		8.32		8.32						± 0.31	
HT HCT		%		36.5		36.5		36.5		36.5		36.5		36.5		36.5						± 2.5	
		l/l		0.365		0.365		0.365		0.365		0.365		0.365		0.365						± 0.025	
VGM MCV		µm ³ ; fl		81		81		81		81		81		81		81						± 4	
TGMH MCH		pg		29.8		29.6		29.6		29.6		29.8		29.8		29.8						± 2.0	
		fmol		1.85		1.84		1.84		1.84		1.85		1.85		1.85						± 0.12	
CCMH MCHC		g/dl		36.8		36.5		36.5		36.5		36.8		36.8		36.8						± 3.0	
		g/l		368		365		365		365		368		368		368						± 30	
		mmol/l		22.83		22.66		22.66		22.66		22.83		22.83		22.83						± 1.86	
IDR RDW		%		13.8		13.5		N/A		13.5		13.5		13.5		13.5						± 3.0	
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		265		255		250		255		255		255		255						± 40	
VPM MPV		µm ³ ; fl		8.7		8.3		N/A		8.9		8.0		8.0		8.0						± 2.0	
LYMPHO		#		2.62		2.58		N/A		2.09		2.58		2.58		2.58						± 0.40	
		%		34.0		33.5		N/A		27.5		33.5		33.5		33.5						± 5.0	
MONO		#		0.54		0.58		N/A		0.65		0.58		0.58		0.58						± 0.30	
		%		7.0		7.5		N/A		8.5		7.5		7.5		7.5						± 4.0	
GRANULO		#		4.54		4.54		N/A		4.86		4.54		4.54		4.54						± 0.60	
		%		59.0		59.0		N/A		64.0		59.0		59.0		59.0						± 6.0	

Ref: TEMP-0830 Rev.37 FRONT / RECTO 1219930084-B



LOT		MX 095 H		CONTROL H		 (Exp.)		2015-12-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	19.9	20.1	19.8	19.9	20.2													± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.70	5.65	5.67	5.68	5.70													± 0.20			
HB	HGB	g/dl	18.2	18.0	18.0	18.2	18.1													± 0.6			
		g/l	182	180	180	182	181												± 6				
		mmol/l	11.30	11.18	11.18	11.30	11.24												± 0.37				
HT	HCT	%	50.2	49.7	49.9	50.0	50.2												± 3.0				
		l/l	0.502	0.497	0.499	0.500	0.502											± 0.030					
VGM	MCV	µm ³ ; fl	88	88	88	88	88													± 4			
TGMH	MCH	pg	31.9	31.9	31.7	32.0	31.8													± 2.0			
		fmol	1.98	1.98	1.97	1.99	1.97												± 0.12				
CCMH	MCHC	g/dl	36.3	36.2	36.1	36.4	36.1													± 3.0			
		g/l	363	362	361	364	361												± 30				
		mmol/l	22.53	22.48	22.40	22.61	22.41												± 1.86				
IDR	RDW	%	13.0	13.0	N/A	13.0	13.0												± 3.0				
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	515	505	520	510	505													± 45			
VPM	MPV	µm ³ ; fl	8.4	7.9	N/A	8.5	7.5													± 2.0			
LYMPHO	#	2.89	2.91	N/A	2.39	2.83														± 1.00			
	%	14.5	14.5	N/A	12.0	14.0														± 5.0			
MONO	#	1.09	1.11	N/A	0.80	1.11														± 0.80			
	%	5.5	5.5	N/A	4.0	5.5														± 4.0			
GRANULO	#	15.92	16.08	N/A	16.72	16.26														± 1.40			
	%	80.0	80.0	N/A	84.0	80.5														± 7.0			