

| LOT                      |                                                        | MX 115 L        |                | Rev 1  |        | CONTROL                    |                | L      |        |  (Exp.) |                         | 2016-02-05<br>( YYYY - MM - DD ) |  |
|--------------------------|--------------------------------------------------------|-----------------|----------------|--------|--------|----------------------------|----------------|--------|--------|--------------------------------------------------------------------------------------------|-------------------------|----------------------------------|--|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        | ABX Lyse        |                |        |        |                            |                |        |        |                                                                                            | TOLERANCES<br>TOLERANCE |                                  |  |
|                          |                                                        | MICROS          | MICROS         | MICROS | MICROS | MICROS                     | ARGOS - HELIOS | PENTRA | PENTRA | PENTRA                                                                                     |                         |                                  |  |
|                          |                                                        | ES60<br>Care ST | CRP<br>CRP 200 | 8P     | 18P    | 60 8P - 60 18P<br>ADVIA 60 |                |        |        |                                                                                            |                         |                                  |  |
| GB WBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 2.1             | 2.1            | 2.1    | 2.1    | 2.1                        |                |        |        |                                                                                            |                         | ± 0.4                            |  |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 2.35            | 2.37           | 2.36   | 2.35   | 2.33                       |                |        |        |                                                                                            |                         | ± 0.15                           |  |
| HB HGB                   | g/dl                                                   | 5.9             | 6.0            | 5.8    | 5.8    | 5.9                        |                |        |        |                                                                                            |                         | ± 0.4                            |  |
|                          | g/l                                                    | 59              | 60             | 58     | 58     | 59                         |                |        |        |                                                                                            |                         | ± 4                              |  |
|                          | mmol/l                                                 | 3.66            | 3.73           | 3.60   | 3.60   | 3.66                       |                |        |        |                                                                                            |                         | ± 0.25                           |  |
| HT HCT                   | %                                                      | 16.0            | 16.1           | 16.3   | 16.5   | 16.1                       |                |        |        |                                                                                            |                         | ± 2.0                            |  |
|                          | l/l                                                    | 0.160           | 0.161          | 0.163  | 0.165  | 0.161                      |                |        |        |                                                                                            |                         | ± 0.020                          |  |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 68              | 68             | 69     | 70     | 69                         |                |        |        |                                                                                            |                         | ± 4                              |  |
| TGMH MCH                 | pg                                                     | 25.1            | 25.3           | 24.6   | 24.7   | 25.3                       |                |        |        |                                                                                            |                         | ± 2.0                            |  |
|                          | fmol                                                   | 1.56            | 1.57           | 1.53   | 1.53   | 1.57                       |                |        |        |                                                                                            |                         | ± 0.12                           |  |
| CCMH MCHC                | g/dl                                                   | 36.9            | 37.2           | 35.6   | 35.3   | 36.7                       |                |        |        |                                                                                            |                         | ± 3.0                            |  |
|                          | g/l                                                    | 369             | 372            | 356    | 353    | 367                        |                |        |        |                                                                                            |                         | ± 30                             |  |
|                          | mmol/l                                                 | 22.93           | 23.12          | 22.12  | 21.90  | 22.79                      |                |        |        |                                                                                            |                         | ± 1.86                           |  |
| IDR RDW                  | %                                                      | 14.5            | 13.7           | N/A    | 14.6   | 14.0                       |                |        |        |                                                                                            |                         | ± 3.0                            |  |
| PLAQ. PLTS               | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 71              | 70             | 68     | 66     | 71                         |                |        |        |                                                                                            |                         | ± 20                             |  |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 8.9             | 8.9            | N/A    | 9.0    | 8.0                        |                |        |        |                                                                                            |                         | ± 2.0                            |  |
| LYMPHO                   | #                                                      | 1.33            | 1.32           | N/A    | 1.10   | 1.32                       |                |        |        |                                                                                            |                         | ± 0.40                           |  |
|                          | %                                                      | 63.5            | 63.0           | N/A    | 52.5   | 63.0                       |                |        |        |                                                                                            |                         | ± 8.0                            |  |
| MONO                     | #                                                      | 0.20            | 0.20           | N/A    | 0.27   | 0.20                       |                |        |        |                                                                                            |                         | ± 0.20                           |  |
|                          | %                                                      | 8.5             | 9.5            | N/A    | 13.0   | 9.0                        |                |        |        |                                                                                            |                         | ± 6.0                            |  |
| GRANULO                  | #                                                      | 0.59            | 0.58           | N/A    | 0.72   | 0.59                       |                |        |        |                                                                                            |                         | ± 0.40                           |  |
|                          | %                                                      | 28.0            | 27.5           | N/A    | 34.5   | 28.0                       |                |        |        |                                                                                            |                         | ± 7.0                            |  |

| LOT                      |                                                        | MX 115 N        |                | CONTROL |        | N                          |                |  (Exp.) |        | 2016-02-05<br>(YYYY - MM - DD) |                         |
|--------------------------|--------------------------------------------------------|-----------------|----------------|---------|--------|----------------------------|----------------|--------------------------------------------------------------------------------------------|--------|--------------------------------|-------------------------|
| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                        | ABX Lyse        |                |         |        |                            |                |                                                                                            |        |                                | TOLERANCES<br>TOLERANCE |
|                          |                                                        | MICROS          | MICROS         | MICROS  | MICROS | MICROS                     | ARGOS - HELIOS | PENTRA                                                                                     | PENTRA | PENTRA                         |                         |
|                          |                                                        | ES60<br>Care ST | CRP<br>CRP 200 | 8P      | 18P    | 60 8P - 60 18P<br>ADVIA 60 |                |                                                                                            |        |                                |                         |
| GB WBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 7.9             | 8.0            | 7.8     | 7.8    | 7.9                        |                |                                                                                            |        |                                | ± 0.8                   |
| GR RBC                   | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 4.54            | 4.55           | 4.51    | 4.50   | 4.52                       |                |                                                                                            |        |                                | ± 0.18                  |
| HB HGB                   | g/dl                                                   | 13.4            | 13.4           | 13.3    | 13.3   | 13.4                       |                |                                                                                            |        |                                | ± 0.5                   |
|                          | g/l                                                    | 134             | 134            | 133     | 133    | 134                        |                |                                                                                            |        |                                | ± 5                     |
|                          | mmol/l                                                 | 8.32            | 8.32           | 8.26    | 8.26   | 8.32                       |                |                                                                                            |        |                                | ± 0.31                  |
| HT HCT                   | %                                                      | 36.8            | 36.9           | 37.0    | 36.5   | 37.1                       |                |                                                                                            |        |                                | ± 2.5                   |
|                          | l/l                                                    | 0.368           | 0.369          | 0.370   | 0.365  | 0.371                      |                |                                                                                            |        |                                | ± 0.025                 |
| VGM MCV                  | µm <sup>3</sup> ; fl                                   | 81              | 81             | 82      | 81     | 82                         |                |                                                                                            |        |                                | ± 4                     |
| TGMH MCH                 | pg                                                     | 29.5            | 29.5           | 29.5    | 29.6   | 29.6                       |                |                                                                                            |        |                                | ± 2.0                   |
|                          | fmol                                                   | 1.83            | 1.83           | 1.83    | 1.84   | 1.84                       |                |                                                                                            |        |                                | ± 0.12                  |
| CCMH MCHC                | g/dl                                                   | 36.4            | 36.4           | 36.0    | 36.5   | 36.2                       |                |                                                                                            |        |                                | ± 3.0                   |
|                          | g/l                                                    | 364             | 364            | 360     | 365    | 362                        |                |                                                                                            |        |                                | ± 30                    |
|                          | mmol/l                                                 | 22.63           | 22.58          | 22.33   | 22.66  | 22.45                      |                |                                                                                            |        |                                | ± 1.86                  |
| IDR RDW                  | %                                                      | 13.7            | 12.9           | N/A     | 14.0   | 13.5                       |                |                                                                                            |        |                                | ± 3.0                   |
| PLAQ. PLTS               | 10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 255             | 255            | 250     | 245    | 250                        |                |                                                                                            |        |                                | ± 40                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                   | 8.2             | 8.0            | N/A     | 8.3    | 7.5                        |                |                                                                                            |        |                                | ± 2.0                   |
| LYMPHO                   | #                                                      | 2.53            | 2.56           | N/A     | 2.03   | 2.49                       |                |                                                                                            |        |                                | ± 0.40                  |
|                          | %                                                      | 32.0            | 32.0           | N/A     | 26.0   | 31.5                       |                |                                                                                            |        |                                | ± 5.0                   |
| MONO                     | #                                                      | 0.51            | 0.56           | N/A     | 0.59   | 0.51                       |                |                                                                                            |        |                                | ± 0.30                  |
|                          | %                                                      | 6.5             | 7.0            | N/A     | 7.5    | 6.5                        |                |                                                                                            |        |                                | ± 4.0                   |
| GRANULO                  | #                                                      | 4.86            | 4.88           | N/A     | 5.19   | 4.90                       |                |                                                                                            |        |                                | ± 0.60                  |
|                          | %                                                      | 61.5            | 61.0           | N/A     | 66.5   | 62.0                       |                |                                                                                            |        |                                | ± 6.0                   |

Ref: TEMP-0830 Rev.37 FRONT / RECTO 1219930084-B



| LOT                      |      | MX 115 H                                               |       | CONTROL H       |       |  (Exp.) |       | 2016-02-05<br>( YYYY - MM - DD ) |  |        |  |                         |  |                            |  |                |  |        |         |        |  |        |  |
|--------------------------|------|--------------------------------------------------------|-------|-----------------|-------|--------------------------------------------------------------------------------------------|-------|----------------------------------|--|--------|--|-------------------------|--|----------------------------|--|----------------|--|--------|---------|--------|--|--------|--|
|                          |      | Rev 1                                                  |       |                 |       |                                                                                            |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         |        |  |        |  |
| PARAMETRES<br>PARAMETERS |      | UNITES<br>UNITS                                        |       | ABX Lyse        |       |                                                                                            |       |                                  |  |        |  | TOLERANCES<br>TOLERANCE |  |                            |  |                |  |        |         |        |  |        |  |
|                          |      |                                                        |       | MICROS          |       | MICROS                                                                                     |       | MICROS                           |  | MICROS |  |                         |  | MICROS                     |  | ARGOS - HELIOS |  | PENTRA |         | PENTRA |  | PENTRA |  |
|                          |      |                                                        |       | ES60<br>Care ST |       | CRP<br>CRP 200                                                                             |       | 8P                               |  | 18P    |  |                         |  | 60 8P - 60 18P<br>ADVIA 60 |  |                |  |        |         |        |  |        |  |
| GB                       | WBC  | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 20.1  | 20.3            | 19.8  | 20.0                                                                                       | 20.1  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 1.6  |  |        |  |
| GR                       | RBC  | 10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l | 5.80  | 5.77            | 5.72  | 5.75                                                                                       | 5.76  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 0.20 |  |        |  |
| HB                       | HGB  | g/dl                                                   | 18.3  | 18.3            | 18.1  | 18.4                                                                                       | 18.3  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 0.6  |  |        |  |
|                          |      | g/l                                                    | 183   | 183             | 181   | 184                                                                                        | 183   |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 6     |        |  |        |  |
|                          |      | mmol/l                                                 | 11.36 | 11.36           | 11.24 | 11.43                                                                                      | 11.36 |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 0.37  |        |  |        |  |
| HT                       | HCT  | %                                                      | 51.0  | 51.4            | 50.9  | 51.2                                                                                       | 51.3  |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 3.0   |        |  |        |  |
|                          |      | l/l                                                    | 0.510 | 0.514           | 0.509 | 0.512                                                                                      | 0.513 |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 0.030 |        |  |        |  |
| VGM                      | MCV  | µm <sup>3</sup> ; fl                                   | 88    | 89              | 89    | 89                                                                                         | 89    |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 4    |  |        |  |
| TGMH                     | MCH  | pg                                                     | 31.6  | 31.7            | 31.6  | 32.0                                                                                       | 31.8  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 2.0  |  |        |  |
|                          |      | fmol                                                   | 1.96  | 1.97            | 1.97  | 1.99                                                                                       | 1.97  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 0.12 |  |        |  |
| CCMH                     | MCHC | g/dl                                                   | 35.9  | 35.6            | 35.6  | 36.0                                                                                       | 35.7  |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 3.0  |  |        |  |
|                          |      | g/l                                                    | 359   | 356             | 356   | 360                                                                                        | 357   |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 30    |        |  |        |  |
|                          |      | mmol/l                                                 | 22.27 | 22.13           | 22.08 | 22.33                                                                                      | 22.17 |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 1.86  |        |  |        |  |
| IDR                      | RDW  | %                                                      | 13.3  | 12.7            | N/A   | 13.5                                                                                       | 13.0  |                                  |  |        |  |                         |  |                            |  |                |  |        | ± 3.0   |        |  |        |  |
| PLAQ.                    | PLTS | 10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l  | 525   | 520             | 530   | 520                                                                                        | 510   |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 45   |  |        |  |
| VPM                      | MPV  | µm <sup>3</sup> ; fl                                   | 7.9   | 7.7             | N/A   | 8.0                                                                                        | 7.2   |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 2.0  |  |        |  |
| LYMPHO                   | #    | 3.22                                                   | 3.05  | N/A             | 2.60  | 3.12                                                                                       |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 1.00 |  |        |  |
|                          | %    | 16.0                                                   | 15.0  | N/A             | 13.0  | 15.5                                                                                       |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 5.0  |  |        |  |
| MONO                     | #    | 1.11                                                   | 1.32  | N/A             | 0.80  | 1.11                                                                                       |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 0.80 |  |        |  |
|                          | %    | 5.5                                                    | 6.5   | N/A             | 4.0   | 5.5                                                                                        |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 4.0  |  |        |  |
| GRANULO                  | #    | 15.78                                                  | 15.94 | N/A             | 16.60 | 15.88                                                                                      |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 1.40 |  |        |  |
|                          | %    | 78.5                                                   | 78.5  | N/A             | 83.0  | 79.0                                                                                       |       |                                  |  |        |  |                         |  |                            |  |                |  |        |         | ± 7.0  |  |        |  |