



LOT		MX 016 L		CONTROL		L		 (Exp.)		2016-04-05 (YYYY - MM - DD)			
		Rev 1		ABX Lyse									
PARAMETRES PARAMETERS		UNITES UNITS		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	TOLERANCES TOLERANCE
				ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.0	2.0	1.9	2.0						± 0.4
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.41	2.43	2.46	2.42	2.42						± 0.15
HB	HGB	g/dl	6.0	6.1	6.0	5.9	6.0						± 0.4
		g/l	60	61	60	59	60						± 4
		mmol/l	3.73	3.79	3.73	3.66	3.73						± 0.25
HT	HCT	%	16.4	16.5	17.0	16.5	16.5						± 2.0
		l/l	0.164	0.165	0.170	0.165	0.165						± 0.020
VGM	MCV	µm ³ ; fl	68	68	69	68	68						± 4
TGMH	MCH	pg	24.9	25.1	24.4	24.4	24.8						± 2.0
		fmol	1.55	1.56	1.51	1.51	1.54						± 0.12
CCMH	MCHC	g/dl	36.6	36.9	35.3	35.9	36.5						± 3.0
		g/l	366	369	353	359	365						± 30
		mmol/l	22.74	22.92	21.95	22.26	22.64						± 1.86
IDR	RDW	%	14.3	14.0	N/A	14.8	14.0						± 3.0
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	78	75	73	74	75						± 20
VPM	MPV	µm ³ ; fl	9.1	8.8	N/A	9.3	8.5						± 2.0
LYMPHO	#	1.30	1.35	N/A	1.03	1.31							± 0.40
	%	65.0	67.5	N/A	54.0	65.5							± 8.0
MONO	#	0.20	0.20	N/A	0.28	0.20							± 0.20
	%	10.0	9.5	N/A	14.5	10.0							± 6.0
GRANULO	#	0.50	0.46	N/A	0.60	0.49							± 0.40
	%	25.0	23.0	N/A	31.5	24.5							± 7.0

LOT		MX 016 N		CONTROL N				 (Exp.)		2016-04-05 (YYYY - MM - DD)	
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.3	7.2	7.3	7.2	7.4					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.63	4.62	4.66	4.60	4.64					± 0.18
HB HGB	g/dl	13.4	13.4	13.4	13.4	13.4					± 0.5
	g/l	134	134	134	134	134					± 5
	mmol/l	8.32	8.32	8.32	8.32	8.32					± 0.31
HT HCT	%	36.6	37.0	37.3	36.3	37.1					± 2.5
	l/l	0.366	0.370	0.373	0.363	0.371					± 0.025
VGM MCV	µm ³ ; fl	79	80	80	79	80					± 4
TGMH MCH	pg	28.9	29.0	28.8	29.1	28.9					± 2.0
	fmol	1.80	1.80	1.79	1.81	1.79					± 0.12
CCMH MCHC	g/dl	36.6	36.3	35.9	36.9	36.1					± 3.0
	g/l	366	363	359	369	361					± 30
	mmol/l	22.75	22.51	22.32	22.90	22.42					± 1.86
IDR RDW	%	13.7	13.5	N/A	14.2	13.5					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	257	253	250	255	256					± 40
VPM MPV	µm ³ ; fl	8.4	8.2	N/A	8.5	7.8					± 2.0
LYMPHO	#	2.37	2.34	N/A	2.02	2.44					± 0.50
	%	32.5	32.5	N/A	28.0	33.0					± 6.0
MONO	#	0.62	0.61	N/A	0.58	0.59					± 0.40
	%	8.5	8.5	N/A	8.0	8.0					± 5.0
GRANULO	#	4.31	4.25	N/A	4.61	4.37					± 0.60
	%	59.0	59.0	N/A	64.0	59.0					± 7.0

Ref: TEMP-0830 Rev.37 FRONT / RECTO 1219930084-B



LOT		MX 016 H		CONTROL H		 (Exp.)		2016-04-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l	19.7	19.6	19.6	19.6	19.8													± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.66	5.62	5.68	5.66	5.66													± 0.20			
HB	HGB	g/dl	18.1	18.2	18.1	18.3	18.1													± 0.6			
		g/l	181	182	181	183	181												± 6				
		mmol/l	11.24	11.30	11.24	11.36	11.24												± 0.37				
HT	HCT	%	50.4	50.0	50.6	50.4	50.4												± 3.0				
		l/l	0.504	0.500	0.506	0.504	0.504												± 0.030				
VGM	MCV	µm ³ ; fl	89	89	89	89	89													± 4			
TGMH	MCH	pg	32.0	32.4	31.9	32.3	32.0													± 2.0			
		fmol	1.99	2.01	1.98	2.01	1.99													± 0.12			
CCMH	MCHC	g/dl	35.9	36.4	35.8	36.3	35.9													± 3.0			
		g/l	359	364	358	363	359												± 30				
		mmol/l	22.31	22.60	22.23	22.56	22.31												± 1.86				
IDR	RDW	%	13.3	13.0	N/A	13.1	13.0												± 3.0				
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l	512	503	520	520	502													± 45			
VPM	MPV	µm ³ ; fl	8.2	8.0	N/A	8.3	7.7													± 2.0			
LYMPHO	#	3.15	2.94	N/A	2.55	3.07														± 1.00			
	%	16.0	15.0	N/A	13.0	15.5														± 5.0			
MONO	#	1.08	1.18	N/A	0.80	1.09														± 0.80			
	%	5.5	6.0	N/A	4.0	5.5														± 4.0			
GRANULO	#	15.46	15.48	N/A	16.27	15.64														± 1.40			
	%	78.5	79.0	N/A	83.0	79.0														± 7.0			