



LOT		MX 036 L		CONTROL		L		 (Exp.)		2016-06-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS											TOLERANCES TOLERANCE								
		MICROS		MICROS		MICROS		MICROS		MICROS			ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60														
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	2.0	2.1		2.0	2.1													± 0.4	
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.44	2.44	2.46	2.43	2.42													± 0.15	
HB HGB	g/dl	6.3	6.2	6.2	6.1	6.2													± 0.4	
	g/l	63	62	62	61	62													± 4	
	mmol/l	3.91	3.85	3.85	3.79	3.85													± 0.25	
HT HCT	%	17.1	17.1	17.2	17.0	16.9													± 2.0	
	l/l	0.171	0.171	0.172	0.170	0.169													± 0.020	
VGM MCV	µm ³ ; fl	70	70	70	70	70													± 4	
TGMH MCH	pg	25.8	25.4	25.2	25.1	25.6													± 2.0	
	fmol	1.60	1.58	1.57	1.56	1.59													± 0.12	
CCMH MCHC	g/dl	36.9	36.3	36.0	35.9	36.6													± 3.0	
	g/l	369	363	360	359	366													± 30	
	mmol/l	22.91	22.54	22.36	22.27	22.73													± 1.86	
IDR RDW	%	14.5	14.0	N/A	14.4	14.0													± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	75	70	77	77													± 20	
VPM MPV	µm ³ ; fl	9.2	9.0	N/A	9.0	8.2													± 2.0	
LYMPHO	#	1.26	1.34	N/A	1.04	1.33													± 0.40	
	%	63.0	64.0	N/A	52.0	63.5													± 8.0	
MONO	#	0.21	0.23	N/A	0.31	0.21													± 0.20	
	%	10.5	11.0	N/A	15.5	10.0													± 6.0	
GRANULO	#	0.53	0.53	N/A	0.65	0.56													± 0.40	
	%	26.5	25.0	N/A	32.5	26.5													± 7.0	

LOT		MX 036 N		CONTROL N		 (Exp.)		2016-06-05 (YYYY - MM - DD)			
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	7.5	7.5	7.3	7.2	7.5					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.66	4.66	4.67	4.64	4.64					± 0.18
HB HGB	g/dl	13.6	13.4	13.5	13.5	13.5					± 0.5
	g/l	136	134	135	135	135					± 5
	mmol/l	8.45	8.32	8.38	8.38	8.38					± 0.31
HT HCT	%	37.3	37.3	37.8	37.6	37.6					± 2.5
	l/l	0.373	0.373	0.378	0.376	0.376					± 0.025
VGM MCV	µm ³ ; fl	80	80	81	81	81					± 4
TGMH MCH	pg	29.2	28.8	28.9	29.1	29.1					± 2.0
	fmol	1.81	1.79	1.80	1.81	1.81					± 0.12
CCMH MCHC	g/dl	36.5	35.9	35.7	35.9	35.9					± 3.0
	g/l	365	359	357	359	359					± 30
	mmol/l	22.65	22.32	22.16	22.31	22.31					± 1.86
IDR RDW	%	14.0	13.5	N/A	14.0	13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	261	250	250	250	255					± 40
VPM MPV	µm ³ ; fl	8.4	8.0	N/A	8.3	7.7					± 2.0
LYMPHO	#	2.44	2.33	N/A	1.93	2.40					± 0.40
	%	32.5	31.0	N/A	26.8	32.0					± 5.0
MONO	#	0.56	0.60	N/A	0.56	0.53					± 0.30
	%	7.5	8.0	N/A	7.8	7.0					± 4.0
GRANULO	#	4.50	4.58	N/A	4.71	4.58					± 0.60
	%	60.0	61.0	N/A	65.4	61.0					± 6.0

Ref: TEMP-0830 Rev.37 FRONT / RECTO 1219930084-B



LOT		MX 036 H		CONTROL H		 (Exp.)		2016-06-05 (YYYY - MM - DD)															
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB	WBC	10 ³ /mm ³ ; 10 ⁹ /l		19.6		19.7		19.3		19.3		19.8								± 1.6			
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.66		5.66		5.66		5.64		5.63								± 0.20			
HB	HGB	g/dl		18.0		17.8		17.9		18.0		17.9								± 0.6			
		g/l		180		178		179		180		179								± 6			
		mmol/l		11.18		11.05		11.12		11.18		11.12								± 0.37			
HT	HCT	%		49.8		49.8		50.4		50.2		50.1								± 3.0			
		l/l		0.498		0.498		0.504		0.502		0.501								± 0.030			
VGM	MCV	µm ³ ; fl		88		88		89		89		89								± 4			
TGMH	MCH	pg		31.8		31.4		31.6		31.9		31.8								± 2.0			
		fmol		1.97		1.95		1.96		1.98		1.97								± 0.12			
CCMH	MCHC	g/dl		36.1		35.7		35.5		35.9		35.7								± 3.0			
		g/l		361		357		355		359		357								± 30			
		mmol/l		22.44		22.19		22.07		22.27		22.18								± 1.86			
IDR	RDW	%		13.5		13.0		N/A		13.3		13.0								± 3.0			
PLAQ.	PLTS	10 ³ /mm ³ ; 10 ⁹ /l		524		515		525		515		510								± 45			
VPM	MPV	µm ³ ; fl		8.0		7.8		N/A		7.8		7.5								± 2.0			
LYMPHO	#			3.04		2.96		N/A		2.41		3.07								± 1.00			
	%			15.5		15.0		N/A		12.5		15.5								± 5.0			
MONO	#			0.98		0.99		N/A		0.77		0.89								± 0.80			
	%			5.0		5.0		N/A		4.0		4.5								± 4.0			
GRANULO	#			15.58		15.76		N/A		16.12		15.84								± 1.40			
	%			79.5		80.0		N/A		83.5		80.0								± 7.0			