



LOT		MX 056 L		CONTROL		L		 (Exp.)		2016-08-05 (YYYY - MM - DD)										
		Rev 1		ABX Lyse																
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS		MICROS		MICROS		MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA		TOLERANCES TOLERANCE
		ES60 Care ST		CRP CRP 200		8P		18P		60 8P - 60 18P ADVIA 60										
GB	WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		2.0		2.0		2.0		2.0										± 0.4
GR	RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.39		2.39		2.41		2.39										± 0.15
HB	HGB	g/dl		6.2		6.1		6.0		6.1										± 0.4
		g/l		62		61		60		61										± 4
		mmol/l		3.85		3.79		3.73		3.79										± 0.25
HT	HCT	%		16.5		16.5		16.9		16.7										± 2.0
		l/l		0.165		0.165		0.169		0.167										± 0.020
VGM	MCV	µm ³ ; fl		69		69		70		70										± 4
TGMH	MCH	pg		25.9		25.5		24.9		25.5										± 2.0
		fmol		1.61		1.58		1.55		1.58										± 0.12
CCMH	MCHC	g/dl		37.6		37.0		35.6		36.5										± 3.0
		g/l		376		370		356		365										± 30
		mmol/l		23.35		22.97		22.09		22.64		22.97								± 1.86
IDR	RDW	%		14.6		14.5		N/A		15.0										± 3.0
PLAQ.	PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l		75		71		68		73										± 20
VPM	MPV	µm ³ ; fl		9.4		9.0		N/A		8.8										± 2.0
LYMPHO		#		1.26		1.27		N/A		1.01		1.25								± 0.40
		%		63.0		63.5		N/A		50.5		62.5								± 8.0
MONO		#		0.21		0.21		N/A		0.33		0.22								± 0.20
		%		10.5		10.5		N/A		16.5		11.0								± 6.0
GRANULO		#		0.53		0.52		N/A		0.66		0.53								± 0.40
		%		26.5		26.0		N/A		33.0		26.5								± 7.0

LOT		MX 056 N		CONTROL N				 (Exp.)		2016-08-05 (YYYY - MM - DD)													
		Rev 1																					
PARAMETRES PARAMETERS		UNITES UNITS		ABX Lyse								TOLERANCES TOLERANCE											
				MICROS		MICROS		MICROS		MICROS				MICROS		ARGOS - HELIOS		PENTRA		PENTRA		PENTRA	
				ES60 Care ST		CRP CRP 200		8P		18P				60 8P - 60 18P ADVIA 60									
GB WBC		10 ⁹ /mm ³ ; 10 ⁹ /l		8.0		7.8		8.0		7.7		8.0									± 0.8		
GR RBC		10 ⁹ /mm ³ ; 10 ¹² /l		4.63		4.59		4.63		4.60		4.60									± 0.18		
HB HGB		g/dl		13.4		13.1		13.3		13.3		13.2									± 0.5		
		g/l		134		131		133		133		132									± 5		
		mmol/l		8.32		8.14		8.26		8.26		8.20									± 0.31		
HT HCT		%		37.0		36.7		37.0		36.8		36.8									± 2.5		
		l/l		0.370		0.367		0.370		0.368		0.368									± 0.025		
VGM MCV		µm ³ ; fl		80		80		80		80		80									± 4		
TGMH MCH		pg		28.9		28.5		28.7		28.9		28.7									± 2.0		
		fmol		1.80		1.77		1.78		1.80		1.78									± 0.12		
CCMH MCHC		g/dl		36.2		35.7		35.9		36.1		35.9									± 3.0		
		g/l		362		357		359		361		359									± 30		
		mmol/l		22.47		22.15		22.30		22.44		22.28									± 1.86		
IDR RDW		%		14.0		13.8		N/A		14.2		13.5									± 3.0		
PLAQ. PLTS		10 ³ /mm ³ ; 10 ⁹ /l		261		258		258		255		255									± 40		
VPM MPV		µm ³ ; fl		8.9		8.5		N/A		8.7		7.8									± 2.0		
LYMPHO		#		2.56		2.46		N/A		1.96		2.56									± 0.40		
		%		32.0		31.5		N/A		25.5		32.0									± 5.0		
MONO		#		0.64		0.66		N/A		0.65		0.68									± 0.30		
		%		8.0		8.5		N/A		8.5		8.5									± 4.0		
GRANULO		#		4.80		4.68		N/A		5.08		4.76									± 0.60		
		%		60.0		60.0		N/A		66.0		59.5									± 6.0		

Ref: TEMP-0830 Rev.37 FRONT / RECTO 1219930084-B

LOT		MX 056 H		CONTROL H		 (Exp.)		2016-08-05 (YYYY - MM - DD)			
		Rev 1									
PARAMETRES PARAMETERS	UNITES UNITS	ABX Lyse									TOLERANCES TOLERANCE
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA	PENTRA	
		ES60 Care ST	CRP CRP 200	8P	18P	60 8P - 60 18P ADVIA 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	20.2	19.8	20.0	19.9	20.2					± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.84	5.80	5.85	5.83	5.80					± 0.20
HB HGB	g/dl	18.4	18.1	18.2	18.4	18.2					± 0.6
	g/l	184	181	182	184	182					± 6
	mmol/l	11.43	11.24	11.30	11.43	11.30					± 0.37
HT HCT	%	52.0	51.6	51.5	51.9	51.6					± 3.0
	l/l	0.520	0.516	0.515	0.519	0.516					± 0.030
VGM MCV	µm ³ ; fl	89	89	88	89	89					± 4
TGMH MCH	pg	31.5	31.2	31.1	31.6	31.4					± 2.0
	fmol	1.96	1.94	1.93	1.96	1.95					± 0.12
CCMH MCHC	g/dl	35.4	35.1	35.4	35.5	35.3					± 3.0
	g/l	354	351	354	355	353					± 30
	mmol/l	21.98	21.77	21.95	22.02	21.90					± 1.86
IDR RDW	%	13.4	13.2	N/A	13.6	13.0					± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	523	523	530	530	510					± 45
VPM MPV	µm ³ ; fl	8.4	8.0	N/A	8.2	7.4					± 2.0
LYMPHO	#	3.23	3.07	N/A	2.59	3.23					± 1.00
	%	16.0	15.5	N/A	13.0	16.0					± 5.0
MONO	#	1.11	1.19	N/A	0.80	1.11					± 0.80
	%	5.5	6.0	N/A	4.0	5.5					± 4.0
GRANULO	#	15.86	15.54	N/A	16.52	15.86					± 1.40
	%	78.5	78.5	N/A	83.0	78.5					± 7.0